



Lo-Back Contour Wrap-Around in Stingray Gray
Available in 3 sizes

FEATURES

- Lo-Back™ and Hi-Back™ options available (Hi-Back model includes four-point harness for added trunk support)
- Bacteria and fungal resistant ComfortMesh™ seat
- Latex and silicone free



Hi-Back Contour Wrap-Around in Stingray Gray
Available in 3 sizes

BENEFITS

- Orthopedic back height adjusts vertically
- Wide-stance stable base fits into any tub
- Easily adjustable back heights
- Compact design for easy storage or travel
- Rustproof frame and tear resistant mesh fabric for extra durability
- Mounting Base allows you to upgrade or replace the support back to the next size for a fraction of the cost of a new support

ACCESSORIES



Pelvic Support Belt

ORDER FORM FAX TO 310-305-1718 OR SUBMIT ONLINE

Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO#: _____ Ordered By: _____

Phone: _____ Fax: _____ E-mail: _____

Ship to address: _____

BASE UNIT	LO-BACK WRAP-AROUND			HI-BACK WRAP-AROUND	
	BCWL-7210S \$660	BCWL-7310M \$720.50	BCWL-7410L \$737	BCWH-7610S \$770	BCWH-7710M \$775.50
Product Size	Small	Medium	Large	Small	Medium
Weight Capacity	N/A	N/A	N/A	N/A	N/A
Overall Depth	29"	29"	29"	29"	29"
Seat Depth	14.5"	14.5"	14.5"	14.5"	14.5"
Overall Seat Width	18"	18"	18"	18"	18"
Inside Seat Width	15"	15"	15"	15"	15"
Inside Back Width at Top	7"	11"	13"	7"	8"
Inside Back Width at Bottom	7"	11"	13"	12"	12"
Back Height Range	9" - 14"	10" - 14"	11" - 16"	17" - 22"	19" - 24"
Product Weight	5.5 lb	6.5 lb	7.75 lb	9.5 lb	9.75 lb

ACCESSORIES (fit all sizes of bath chair unless noted)

Pelvic Positioning Belt \$77 BCWZ-A7830