

## Immedia E-Board

Immedia E-Board is a thin and flexible transfer board that allows easier application and removal.



Immedia E-Board with wings.

### E-Board with wings

Immedia E-Board with wings has wings along the sides, which can be bent up or down. This provides more sturdiness and stability. The upward wing is used to protect against the wheel of the wheelchair. Available in two lengths.



Immedia E-Board without wings.

### E-Board Compact

A smaller and lighter board that is easily carried and suitable for the more active and independent user. Immedia E-Board Compact is available in two lengths.



ONLINE TUTORIAL!

### Learn more ...

... about how to use the transferboards. Our new animated film *Transferboards – E-Board* gives you useful instructions for your daily work.



## Immedia Dolphin

Dolphin is a sturdy and stable transfer board uniquely designed for transfer to the toilet. The board has cut-outs for positioning against the wheel of the wheelchair and non-slip pads, which allows good positioning and avoids sliding out of place. The cut-out allows the user to remain on the board while on the toilet.



Immedia Dolphin is uniquely designed for transfer to toilet.

### INTENDED USE

- Transfer from bed to wheelchair / shower chair
- Transfer in and out of the car
- Transfer between wheelchair, chair, toilet, shower chair, commode

DESCRIPTION	WIDTH X LENGTH X HEIGHT		ITEM NO.	SPECIFICATIONS
IMMEDIA	CM	INCH		
~ 3B-Board, curved, white	23 x 77	9.1 x 30.3	IM403	<b>Immedia E-board and 3B-board</b> MATERIAL: Polyethylene
~ E-Board Small, black	20 x 45 x 0,5	7.9 x 17.7 x 0.2	IM410	
~ E-Board Medium, black	25 x 60 x 0,5	9.8 x 23.6 x 0.2	IM408	<b>Immedia Butterfly and Dolphin</b> MATERIAL: Fibreglass
~ E-Board Wide, with wings black	33 x 60 x 0,5	13 x 23 x 0.2	IM404	
~ E-Board Long, with wings, black	33 x 75 x 0,6	13 x 29.5 x 0.2	IM405	
~ Butterfly	32 x 66	12.6 x 26	IM401	<b>Immedia Butterfly GlideWing</b> MATERIAL OUTSIDE & INSIDE: Nylon PADDING: Polyester
~ Butterfly GlideWing L	35 x 40	13.8 x 15.7	IM4001	
~ Butterfly GlideWing XL	35 x 47.5	13.8 x 18.7	IM4002	
~ Dolphin Toilet Board, white	Lenght: 69	27.2	IM409	

### WASHING INSTRUCTIONS

TRANSFER BOARDS: Clean with soapy water or disinfectant.



BUTTERFLY GLIDEWING: Wipe off with soapy water or disinfectant. Do not use fabric softener or bleach.

