


# Immedia MiniBoard


## SOFT TRANSFER STRETCHER

MiniBoard is suitable for transfer between bed and operating table, turning and positioning in bed and when lifting from the floor. The board in combination with the low friction material makes it easier to insert under the user.

There are two options with head support which allow safe and comfortable support for the head during transfer.



*The low friction material facilitates applying Immedia MiniBoard under the user.*



### INTENDED USE

- Horizontal transfer
- Lift from the floor

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
<b>Immedia MiniBoard</b>	<b>CM</b>	<b>INCH</b>		
~ 4 handles, rectangular	44 X 100	17.3 x 39.4	IM8010044	TOP AND UNDERSIDE: Nylon INSIDE: Polyethylene HANDLES: Polyester Zipper/velcro
~ 4 handles, rectangular	54 X 110	21.3 x 43.3	IM8011054	
~ 6 handles, rectangular	54 X 150	21.3 x 59.1	IM8015054	
~ 7 handles, with head support	54 x 124	21.3 x 48.8	IM8012454	
~ 11 handles, with head support	54 x 200	21.3 x 78.7	IM8020054	
<b>WASHING INSTRUCTIONS</b>				
Can be wiped with disinfectant				
Remember to remove the stuffing and the plastic plate before washing.				