

Immedia MultiGlide SG

NYLON SLIDE SHEET WITH LOW-LOW FRICTION

Immedia MultiGlide SG is suitable for all types of situations when there is a necessity to minimize friction under pressure points. It is ultra low friction both inside and outside with the double layer minimising friction.

MultiGlide SG is easy to apply and remove, and should always be removed after use. MultiGlide SG tubular is available in several sizes and also in a 10 meter roll that can be cut to any desired length.



Immedia MultiGlide SG roll that can be cut down to any desired length. Outside: ultra low friction, inside: ultra low friction.



Immedia MultiGlide SG low-low friction is easy to apply and remove.


Immedia MultiGlide SG tubular has a black loop strap for easy storage and differentiation. Immedia MultiGlide SG should always be removed after use.

ONLINE TUTORIAL!



Learn more about ...

... how to use the Immedia MultiGlide. Our new animated film Immedia MultiGlide gives you useful instructions for your daily work.



INTENDED USE

- Minimize friction under the pressure points
- Turning
- Positioning
- Higher up in bed
- Bed to bed transfer
- Dressing

DESCRIPTION

WIDTH X LENGTH

ITEM NO.

SPECIFICATIONS

Immedia MultiGlide SG

	CM	INCH	
~ Tubular nylon sheet, black strap	40 x 40	15.7 x 15.7	IM0040
~ Tubular nylon sheet, black strap	60 x 60	23.6 x 23.6	IM0066
~ Tubular nylon sheet, black strap	72 x 70	28.3 x 27.6	IM0077
~ Tubular nylon sheet, black strap	72 x 100	28.3 x 39.4	IM0080
~ Tubular nylon sheet, black strap	70 x 200	27.6 x 78.7	IM0082
~ Tubular nylon sheet, black strap	72 x 130	28.3 x 51.2	IM0083
~ Tubular nylon sheet, black strap	100 x 70	39.4 x 27.6	IM00100
~ Tubular nylon sheet, black strap	108 x 80	42.5 x 31.5	IM00108
~ Tubular nylon sheet, black strap	120 x 90	47.2 x 35.4	IM00120/1
~ Tubular nylon sheet, black strap	130 x 80	51.2 x 31.5	IM00130
~ Tubular nylon sheet, black strap	200 x 70	78.7 x 27.6	IM00200
~ Tubular nylon sheet, open	150 x 70	59.1 x 27.6	IM00150/70
~ Nylon sheet, on roll	150 x 1000	59.1 x 293.7	IM2245

MultiGlide SG

MATERIAL: Nylon



WASHING INSTRUCTIONS

Can be wiped off.

Do not use fabric softener or bleach.



Width = open side