



Immedia MultiGlide


NYLON TUBULAR SLIDE SHEET WITH HIGH-LOW FRICTION

Immedia MultiGlide is a thin tubular slide sheet with ultra low friction on the inside and high friction on the outside. The outside surface also provides a waterproof barrier. MultiGlide is suitable when there is a need to minimize friction under the pressure points, for instance when transferring and turning in bed, between bed and wheelchair, on the operating table or x-ray table.

MultiGlide is to be applied when necessary and removed after use. Immedia MultiGlide is available in several sizes. MultiGlide has a loop strap for easy storage and colour codes for differentiation. MultiGlide should always be removed after use.












Immedia MultiGlide with high friction outside ensures that it remains in place under the user during transfer and turning. Inside is ultra low friction.



Immedia MultiGlide is a thin tubular slide sheet with ultra low friction on the inside and high friction on the outside.

INTENDED USE

- Minimize friction under the pressure points
- Turning
- Positioning
- Higher up in bed
- Bed to bed transfer

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
Immedia MultiGlide Waterproof	CM	INCH		MultiGlide Waterproof
~ Tubular, blue strap	40 x 40	15.7 x 15.7	IM40	MATERIAL: Nylon with polyurethane coating
~ Tubular, silver strap	60 x 25	23.6 x 9.8	IM62	  
~ Tubular, yellow strap	60 x 60	23.6 x 23.6	IM66	WASHING INSTRUCTIONS
~ Tubular, brown strap	72 x 70	28.3 x 27.6	IM77	Can be wiped with rubbing alcohol.
~ Tubular, red strap	72 x 100	27.6 x 39.4	IM80	Do not use fabric softener or bleach.
~ Tubular, lightbrown strap	70 x 200	27.6 x 78.7	IM82	  
~ Tubular, green strap	72 x 130	27.6 x 51.2	IM83	  
~ Tubular, silver strap	80 x 80	31.5 x 31.5	IM88	
~ Tubular, white strap	100 x 70	39.4 x 27.6	IM100	
~ Tubular, silver strap	100 x 200	39.4 x 78.7	IM102	
~ Tubular, green strap	100 x 135	39.4 x 53.1	IM103/1	
~ Tubular, grey strap	108 x 80	42.5 x 31.5	IM108	
~ Tubular, white strap	120 x 70	47.2 x 27.6	IM120	
~ Tubular, pink strap	130 x 80	51.2 x 31.5	IM130	
~ Tubular, lightblue strap	160 x 80	63 x 31.5	IM168	
~ Tubular, purple strap	200 x 70	78.7 x 27.6	IM200	
~ Open	100 x 140	39.4 x 55.1	IM100/140	Width = open side