

# Immedia Positioning Wedge

## SUPPORT AND POSITIONING CUSHIONS WITH NON-SLIP MATERIAL

Immedia Positioning Wedge supports the user on the side after turning, when resting or when performing hygiene or wound care, saving the caregiver from maintaining a static holding position. The Positioning Wedge can be combined with 4WayGlide System, SatinSheet System, MultiGlide or other slide sheets. The cover is a non-slip Polyurethane material that helps to keep the Positioning Wedge in place independent of surface material.



### INTENDED USE

- Support static positioning
- Resting
- Care and personal hygiene
- Can be used alone or in combination with other positioning aids



*The Immedia Positioning Wedge have unequally angled sides, turn the wedge around and use the angle that fits best for the situation.*