



Support

Immedia SupportBelt


SUPPORT BELT WITH HANDLES

Immedia SupportBelt is used to offer support and assistance during mobilisation and transfers. The diagonally placed handles allows the caregiver to maintain a more ergonomic grip.

SupportBelt is a useful aid when helping a person to stand up or sit down, to support walking or to assist during transfer between bed, wheelchair, chair, etc. The SupportBelt has non-slip material on the inside and an adjustable belt buckle.



Immedia SupportBelt has diagonally placed handles that allows the caregiver to maintain a more ergonomic grip.



INTENDED USE

- Support stand up / sit down
- Support walking training
- Support transfer between bed, wheelchair, chair, etc
- Transfer from wheelchair to toilet

DESCRIPTION

IMMEDIA SupportBelt

- ~ XXS, 4 handles, size 45 x 90 cm
- ~ XS, 4 handles, size 65-105 cm
- ~ S, 4 handles, size 65-105 cm
- ~ M, 6 handles, size 90-130 cm
- ~ L, 8 handles, size 125-150 cm
- ~ XL, 10 handles, size 145-170 cm

WIDTH X LENGTH

CM	INCH
9 x 90	3.5 x 35.4
9 x 105	3.5 x 41.3
13 x 105	5.1 x 41.3
13 x 130	5.1 x 51.2
13 x 150	5.1 x 59.1
15 x 170	5.9 x 66.9

ITEM NO.

- IM8009009
- IM8010509
- IM8010513
- IM8013013
- IM8015013
- IM8017015

SPECIFICATIONS

OUTSIDE: Nylon
 INSIDE: Nylon with PU-coating
 PADDING: Polyethylene
 HANDLES: Polyester strap with polypropylene
 BELT BUCKLE: Acetal, nylon

COLOUR CODE



WASHING INSTRUCTIONS

Can be disinfected with LT-CIN isopropylalcohol 72 %.

Do not use fabric softener or bleach.

