

Immedia Swivel Cushions

FOR TURNING IN SITTING POSITION

Immedia AutoTurn

The Immedia AutoTurn facilitates getting in and out of a car. Its low friction allows the user to slide and turn into the desired position in the car seat. This minimises strain and twisting of feet, knees, hips and back. AutoTurn is fixed around the car seat with a belt and buckle. Handles at the sides can be used to assist the user into place.



Immedia AutoTurn


Immedia EasyTurn

Immedia EasyTurn is a soft swivel cushion that makes it easier for the user to get in and out of a car or bed. EasyTurn is easy to apply and use.



Immedia EasyTurn

Immedia HandySwing Seat



Immedia HandySwing Seat

Immedia HandySwing Seat is a soft swivel cushion specifically suitable to support turning in and out of bed or chair. HandySwing Seat is softer and foldable which makes it easier to place under the user.

INTENDED USE

- In and out of a car
- In and out of bed

DESCRIPTION

WIDTH X LENGTH

ITEM NO.

SPECIFICATIONS

Immedia

- ~ AutoTurn, swivel cushion for car
- ~ EasyTurn, swivel cushion
- ~ EasyTurn, swivel cushion
- ~ EasyTurn, swivel cushion
- ~ HandySwing Seat

CM

Ø45

INCH

Ø17.7

IM99999

IM999405

IM990454

IM999505

IM106010

Immedia AutoTurn

TOP: Velour (polyester)

INSIDE: Nylon

BACKSIDE: Nylon with PU-coating

Immedia EasyTurn

TOP: Velour (polyester)

INSIDE: Nylon

BACKSIDE: Nylon with polyurethane coating

MIDDLE: Polyethylene

Immedia HandySwing Seat

TOP: Polyester

INSIDE: Nylon

BACKSIDE: Nylon

with polyurethane

coating.

WASHING INSTRUCTIONS

Can be wiped off.



Immedia PediTurn

SWIVEL CUSHION FOR THE FLOOR

Immedia PediTurn is a swivel cushion that facilitates transfer between sitting positions, such as bed to wheelchair. It is used when the user cannot move their feet during the transfer.

Available in two options:


- Immedia PediTurn SOFT: a soft cushion suitable when the user is not wearing shoes. It has a non-slip surface on both sides to avoid sliding out of position.
- Immedia PediTurn HARD: a hard plate suitable when the user is wearing shoes.



Immedia PediTurn HARD



Immedia PediTurn SOFT



INTENDED USE

- Support moving feet during seated transfer
- Transfer between bed and wheelchair
- Transfer between wheelchair, chair, shower chair, toilet, commode
- Transfer in and out of a car

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
Immedia PediTurn	CM	INCH		
~ SOFT, swivel cushion for floor, black	Ø40	Ø15.7	IM99836GT	Immedia PediTurn SOFT TOP: Toughtek INSIDE: Nylon UNDERSIDE: Toughtek
~ HARD, swivel cushion for floor, grey	Ø38	Ø15	IM99700	Immedia PediTurn HARD MATERIAL: ABS/TPU INSIDE: HDPE FIXING: Brass/steel
WASHING INSTRUCTIONS Can be wiped off.				