



Immedia TransferMattress

SOFT PADDED MATTRESS FOR HORIZONTAL TRANSFER


Immedia TransferMattress facilitates horizontal transfer between two supine positions, such as bed to bed, bed to shower trolley, bed to x-ray table, etc. It can also be used for turning in bed, for instance from “face up” to “face down”. Immedia TransferMattress is made of a soft padded cotton/polyester material and available in two widths, 60 cm and 70 cm.

Combined with a cover

Immedia TransferMattress can be combined with a waterproof nylon cover or a disposable cover. Immedia TransferMattress should not be used if there is a gap between the two transfer surfaces.



The quilted TransferMattress provides a soft and comfortable horizontal transfer.



INTENDED USE

- Transfer between beds
- Transfer between bed and shower trolley
- Transfer between bed and stretcher, x-ray table, operating table
- Turning in bed

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
Immedia TransferMattress	CM	INCH		OUTSIDE: 65 % polyester / 35 % cotton INSIDE: Nylon PADDING: Polyester
~ Narrow	60 x 195	23.6 x 77.6	IM60/195	
~ Wide	70 x 195	27.6 x 77.6	IM70/195	
Accessories – Immedia				
~ Nylon cover for IM60/195	60 x 220	23.6 x 86.6	IM60220N	
~ Nylon cover for IM70/195	70 x 200	27.6 x 78.7	IM7000N	
~ Disposable cover, roll with 75 pcs	60 x 200	23.6 x 78.7	IM4030	
Set – Immedia TransferMattress				
~ Narrow Set – includes IM60/195 + IM60220N			IM4304	
~ Wide Set – includes IM70/195 + IM7000N			IM4305	
WASHING INSTRUCTIONS Do not use fabric softener or bleach.				