




Immedia TwinSheet

THE IN-BED SOLUTION FOR INCONTINENT USERS

TwinSheet4Glide system

Immedia TwinSheet4Glide is a 2-part glidemattress suitable for turning and positioning of incontinent users with severely reduced mobility. Suitable also for users with dementia, or for those who feel unsafe, have brain damage, pain, pressure ulcers or are heavy. Disconnecting the lower part in this position provides better access for care of pressure ulcers, hygiene, dressing and even to change the bedding, both the lower part and the top part.



*Immedia
Complete 2-part
TwinSheet4Glide system
including nylon sheet.*

DESCRIPTION	WIDTH X LENGTH		ITEM NO.	SPECIFICATIONS
IMMEDIA TwinSheet4Glide	CM	INCH		Glide Mattress
~ Top part	150 x 75	59.1 x 29.5	IM150/75N	TOP MICROFIBER: Polyester
~ Down part	150 x 140	59.1 x 55.1	IM150/140N	MIDDLE: 80% polyester / 20% viscose
~ LPL Nylon Sheet with NS-lock, NS-sides	85 x 200	33.5 x 78.7	IM85/200LPLNS	PROTECTION BARRIER: Polyurethane
~ LPL Nylon Sheet – NS-lock	85 x 200	33.5 x 78.7	IM85/200LPL	COVER: Polyurethane
~ LPL Nylon Sheet – 3xNS-lock, NS-sides	85 x 200	33.5 x 78.7	IM85/2003LPL	BACKSIDE: Nylon
				HANDLES: Polyester
				ABSORPTION: 2000ml/m ²
				NYLON SHEET: Nylon
				LOCKING SYSTEM: Polyurethane covered nylon

