

# SlingOn

## Short Instruction

**immedia**<sup>®</sup>  
a part of Etac



### Material

Nylon, polyester

Class 1 - Medical device directive  
MDD 93/42/EEC

Do not add fabric conditioner



**Item no**

**EAN no**



**Size in mm**

IM437  
IM437/S  
IM438

5708012446872  
5708012481132  
5708012001828

W1200 x L1500  
W800 x L1000  
W800 x L1700

# SlingOn



**Check always:**

The icons represent: 1. Visual inspection (eye). 2. Proper strap placement and tension (hand holding strap). 3. Weight limit (weight scale labeled 'KG'). 4. Proper strap attachment (hand holding strap).