

INVACARE®  
**9000 XDT**  
WHEELCHAIR



The **Invacare® 9000 XDT Wheelchair** is designed for individuals who demand an extra wide, heavy-duty or tall chair with seat-to-floor heights ranging from hemi through tall. The high-performance 9000 XDT wheelchair provides a 350 lb. weight capacity and the components are interchangeable with the other 9000 series wheelchairs.

## Features

Durable, low-maintenance carbon steel frame is long-lasting

Tall frames are designed to accommodate taller individuals and provide a more comfortable ride

Multiple frame color variations for customizing

Triple-axle positions allow wide variation of seat-to-floor heights

Urethane rear tires mounted on "no flex" wheels offer superior performance

Urethane casters for better ride and lighter weight

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

HCPCS Code: K0007



1. Triple-axle positions
2. Urethane casters
3. Bariatric sizes available
4. Seating available
5. Builder options

Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height (2" increments)	Back Height
9XDT	23"-29"	36"	29"	13"	16"-20", 22" x 16"-20"	17.5"-21.5"	14"-19"
9XDT	8"	24"	37 lb.	350 lb.	55 lb.	34"x13"x 39"	5 Years

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.

