



INVACARE® MARINER™
REHAB SHOWER CHAIR



Yes, you can:

Features

Aluminum frame and stainless steel hardware
 Fits over most standard and elongated toilet bowls
 Padded four-position seat
 Flip-back padded arms that lock down for safety
 Swing-away front riggings

Tool-less adjustable height footrests
 Folds into a compact unit for storage and transport
 Removable pan and pan holder
 Safety strap included
 23" treaded urethane tire and 5" caster models



1. Padded four-position seat
2. Swing-away front riggings with tool-less adjustable height footrests
3. Flip-back padded arms that lock down for safety



| Model # | Dimensions HxWxD | Width Between Arms | Width Folded | Seat Dimensions | Seat-to-Floor Height | Seat-to-Arm Height | Back Dimensions |
|-------------------------|---------------------|-----------------------|-----------------|--------------------|-------------------------|-----------------------|--------------------|
| 6895 w/ 24" Rear Wheels | 39"x26.5"x32" | 18.5" | 12" | 18.25"x18.25" | 22.75" | 8" | 15.5"Hx18"W |
| 6891w/ 5" Casters | 39"x22.5"x29" | 18.5" | 8" | 18.25"x18.25" | 22.75" | 8" | 15.5"Hx18"W |



| | Footrest Extensions | Clearance Over Toilet HxWxD | Weight | Weight Capacity | Units Per Carton | Limited Warranty |
|-------------------------|------------------------|--------------------------------|--------|--------------------|---------------------|---------------------|
| 6895 w/ 24" Rear Wheels | 13"-18" | 20.25"x16"x26.5" | 41 lb. | 300 lb. | 1 | 3 Years |
| 6891w/ 5" Casters | 13"-18" | 20.75"x16"x23.5" | 32 lb. | 300 lb. | 1 | 3 Years |

