


INVACARE®
SOLACE® PREVENTION
SERIES THERAPEUTIC SUPPORT SURFACES



The **Invacare® Solace® Prevention Series Mattresses** are unique foam mattress replacements that are designed to help minimize the risk of pressure injuries by distributing the weight as uniformly as possible to allow for pressure redistribution. The unique and patented designs in the Prevention Series combine comfort and technology for resident and caregiver confidence.

Invacare® Solace® Prevention 1080 Mattress

The convoluted "fortress" cut design offers more surface area for resident support and relief. This high quality foam includes a soft head and heel section. Available in 80" standard length or custom 84" (SPS1084) and optional raised side rail design (SPS1080RSR). The three zones help relieve pressure in the important susceptible body zones.

Features


- 70 Denier nylon top
- 11 oz. tri-laminate vinyl bottom
- Fluid resistant
- Raised side rail design available
- Compliant with federal fire regulations*

Invacare® Solace® Prevention 2080 Mattress

Added to the design are built-in side supports. The head and heel soft sections add comfort and help with pressure redistribution, while the side supports provide stability and durability. Also available in bariatric sizes; SPS2080B42, SPS2080B48 (450 lb. weight capacity on bariatric models).

Features

- Built-in side supports
- 70 Denier nylon top cover
- 11 oz. tri-laminate vinyl bottom
- Soft head and heel section
- Fortress cut convoluted foam



1. Test shows subject in supine position demonstrating contact measurement between key pressure points using BodiTrak mapping technology on the SPS1084 support surface. Subject: Male 5'11", 185 lb.

Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
SPS1080	80" L 36" W 6" H	350 lb.	16 CFR 1632 16 CFR 1633	5 Years (Mattress) 1 Year (Cover)
SPS1084	84" L 36" W 6" H	350 lb.	16 CFR 1632 16 CFR 1633	5 Years (Mattress) 1 Year (Cover)
SPS1080RSR	80" L 36" W 6" H 3" Raised Side Rails	350 lb.	16 CFR 1632 16 CFR 1633	5 Years (Mattress) 1 Year (Cover)

Model #	Dimensions	Weight Capacity	Regulation Compliance	Limited Warranty
SPS2080B42/ SPS2080B48	80" L 42"/48" W 6" H	400 lb.	16 CFR 1633 16 CFR 1633	5 Years (Mattress) 1 Year (Cover)

INVACARE® MATTRESSES



Now you can offer your customers the comfort, protection and durability of superior **Homecare Mattresses** by Invacare, along with **Invacare® Mattress Covers**. Choose from innerspring or foam mattresses. Mattress covers will keep your Invacare® mattress clean and fresh, and protect it from stains, odors, dust and incontinence. Covers are made of non-allergenic waterproof vinyl and are available in zippered and contoured styles.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment.

Invacare® Innerspring Mattress

Model # 5185 , 5185XL

High quality .312 innerspring design

Durable 11-ounce vinyl cover is anti-bacterial, anti-static, acid-resistant and waterproof for easy use and care

5185XL (not shown) is available for use with extender kit

Invacare® Economy Foam Mattress

Model # 5180 , 5184

High quality 1.5 lb. density foam provides comfort and durability in a one-piece construction

Invacare® Bariatric Foam Mattress

Model # BARMATT(39/42/48/EXP)

Constructed of latex-free materials

Durable 11-ounce vinyl cover is anti-bacterial and waterproof for easy use and care

Solid inner core

Mattress is reversible and flippable

Constructed of high resilience foam for extended life and service for the homecare provider

Easy to clean

HCPCS Code: E0271

HCPCS Code: E01272



5185



BARMATT39

1. Invacare® MC0195-1 zippered mattress cover comes individually wrapped; 12 per case

2. Invacare® MC0065-1 contoured mattress cover comes individually wrapped; 12 per case



Model #	Dimensions LxWxD	Gauge	Regulation Compliance	Limited Warranty
MC0165-1, MC0195-1	80"x36"x6"	4"	N/A	30 Days
5180	80"x36"x5.5"	N/A	16CFR1632/33	1 Year
5184	84"x36"x5.5"	N/A	16CFR1632/33	1 Year
5185	80"x16"x6"	N/A	16CFR1632/33	1 Year
5185XL	84"x36"x6"	N/A	16CFR1632/33	1 Year
BARMATT39	80"x39"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT42	80"x42"x6.5"	N/A	16CFR1632/33	1 Year
BARMATT48	80"x48"x6.5"	N/A	16CFR1632/33	1 Year
BARMATTEXP	80"x39"-48"x6.5"	N/A	16CFR1632/33	1 Year



Yes, you can.®

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.