

INVACARE®
TRACER® IV
WHEELCHAIR



The **Invacare® Tracer® IV Wheelchair** features superior durability, rollability and streamlined looks. It is designed for the active individual who needs the strength of a reinforced heavy-duty frame. Its many options include components interchangeable within the Tracer series. The Tracer IV wheelchair is long-lasting and low-maintenance.

Features

Dual axle position allows variation of seat-to-floor heights

Durable, low-maintenance, triple chrome-plated carbon steel frame is long-lasting

Reinforced frame includes gusseted sides, caster journals and dual crossbraces

Corded urethane rear tires (350 lb.) and pneumatic flat free tires (450 lb.) mounted on "no flex" wheels offer superior performance

Low-maintenance, dual-sealed precision bearings

Heavy-duty inner liners keep seat and back from stretching

Upholstery is dual-embossed

HCPCS Code: K0007



1. Dual axle positions
2. Corded urethane tires
3. Dual-embossed upholstery
4. Bariatric sizes available
5. Seating available
6. Builder options



Model #	Overall Width	Overall Height	Overall Depth w/o riggings	Folded Width	Seat Dimensions	Seat-to-Floor Height	Back Height
T420, T422, T424, T4X22, T4X24	29"	36"	30"	13"	18"-24"x18", 20"	17.5", 19.5"	16"



	Caster Height	Wheel Height	Weight w/o Rigging	Weight Capacity	Shipping Weight	Carton Dimensions	Limited Warranty
T420, T422, T424, T4X22, T4X24	8"	24"	42 lb.	350 lb., 450 lb.	53 lb.	32"x13"x 37"	5 Years



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.