

INVACARE®
BAR600IVC
BARIATRIC BED



The **Invacare® BAR600IVC Bariatric Bed** is a heavy-duty full-electric bed frame designed for bariatric individuals. It is capable of supporting patients who weigh up to 600 lb. The bed frame extends to the edges of the bed deck for better support than most other split-spring bariatric beds. Using the hand pendant, the patient or caregiver is able to change the positioning of the head and foot sections, as well as the bed height. Despite being substantial in size, bed setup within the home is easy – no special tools are needed.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new and unused after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment.

Features

Heavy-duty frame design ensures added strength for client support

Heavy-duty bed sleep surface is 17% larger than conventional beds to ensure comfortable positioning

Split-spring design with universal bed ends for ease in setup

No special tools necessary for in-home setup

Easy-release motor mounts allow quick motor removal

Robotically welded, steel frame construction

Pendant allows for dynamic client positioning

Quiet, smooth operation

Use with Invacare bariatric mattress (BARMATT42) and bed rails (BAR6640IVC):

- BARMATT42 meets fire safety regulations 16CFR1633 and 16CFR1632

600 lb. Client weight capacity

HCPCS Code: E0301



1. Hand pendant



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
BAR600IVC	88" L 14"-22" H 42" W	80" L 42" W	Safe Working Load: 700 lb.* Max. Patient Weight: 600 lb.	N/A	N/A	UL962 CSA Certified	5 Years Welds/ Frame/Structure 1 Years Mechanical/Electrical
BAR5000IVC Headspring				114 lb.	42" L 10" W 43" D		5 Years Welds/ Frame/Structure 2 Years Mechanical/Electrical
BAR5490IVC Footspring				97 lb.	48" L 9" W 43" D		5 Years Welds/ Frame/Structure 2 Years Mechanical/Electrical
BAR5301IVC Bed End				55 lb.	46" L 7" W 22" D		

*Includes patient, mattress and all accessories.

©2015
Invacare Corporation.
All rights reserved.
Trademarks are
identified by the symbols
™ and ®. All trademarks
are owned by or licensed
to Invacare Corporation
unless otherwise
noted. Specifications
are subject to change
without notification.
Form No. 00-258
Rev. 08/15
1506107



All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.