

INVACARE® MATRIX®

PSP (Posture Seat Polymer) CUSHION



The **Invacare Matrix PSP** (Posture Seat Polymer) cushion combines the stability and support of the PS foam base with a polymer overlay for extra protection and shear reduction. An ultra-lightweight cushion providing zero maintenance, worry free skin protection and positioning.



matrix
SEATING SERIES



Features and Functional Benefits

- Startex™ outer cover – Reversible can be used smooth or fabric side out
- Moisture resistant and breathable to keep cool and dry
- Includes hook and loop strips
- Startex™ inner cover for extra moisture protection
- Molded high resilient polyurethane foam made with Ultra-fresh™ anti-microbial
- Rear cushion radius and waterfall front edge
- Gently contoured to provide neutral leg position
- Molded polymer insert for extra protection and shear reduction
- Unique waffled relief zone allows immersion in ischial area
- Designed to provide off-loading to hips and thighs for support and stability

HCPCS Code: E2607


Optional Accessories



POK – Pelvic Obliquity Kit (2 pieces) accommodates up to one inch of obliquity



Cushion Rigidizer – lightweight and low profile



Model #	Width	Depth	Height Abductor/Adductor	Height at Leg Trough	Product Weight	Weight Capacity	Limited Warranty
PSP (Standard)	14", 16", 17", 18", 20"	16", 18", 20", 22"	3.5"/3.5"	2.75"	2.2 lb. and up	300 lb.	3 Years

