



Yes, you can.®

TRANSACTION TYPE:

- Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

The Invacare IVC 9000 Jymni pediatric wheelchair offers seat widths and depths to accommodate the pediatric client for short-term rentals.

FRAMES

FRAME TYPE

- 9JYLT 9000 Jymni Light Weight Chair \$885.00

FRAME WIDTH AND DEPTH

- WD00 10" x 10".....No Charge
 Not available with 17" seat-to-floor height, unless chair is a recliner.
- WD20 12" x 10".....No Charge
 Not available with 17" seat-to-floor height, unless chair is a recliner.
- WD22 12" x 12".....No Charge
 Not available with 17" seat-to-floor height, unless chair is a recliner.
- WD42 14" x 12".....No Charge
 Not available with 17" seat-to-floor height, unless chair is a recliner.
- WD44 14" x 14".....No Charge
- WD64 16" x 14".....No Charge

SEAT TO FLOOR HEIGHT

- 15STF 15" Seat to Floor Height No Charge
 Only available with 20" rear wheels.
- 16STF 16" Seat to Floor Height Standard
 Only available with 20" rear wheels.
- 17STF 17" Seat to Floor Height No Charge
 Only available with 22" rear wheels.

ARM TYPE

All above selections come standard with fixed height arms 7 1/2", adjustable height arms 8" - 10" in 1" increments. Removable flip-back arms.

- 18JY Full Length Fixed Height Conventional No Charge
 Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 19JY Full Length Adjustable Ht Conventional \$130.00
 Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 28JY Desk Length Fixed Height Conventional No Charge
 Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 29JY Desk Length Adjustable Ht Conventional \$130.00
 Recliner and angle back must use conventional arms. Arms do not flip back on the recliner.
- 78JY Full Length Fixed Height Space-Saver No Charge
- 79JY Full Length Adjustable Height Space-Saver \$130.00
- 88JY Desk Length Fixed Height Space-Saver No Charge
- 89JY Desk Length Adjustable Height Space-Saver \$130.00

BACK TYPE

All below selections come standard with fixed height arms 7.5", adjustable height arms 8" - 10" in 1" increments. Removable flip-back arms.

- ADJBH Adjustable Ht Back with 10° Bend Standard
 Adjustable back height 14"-15". Factory set at 14".
- R Recliner Frame \$390.00
 Recliner is only available in U240 black nylon upholstery. Recliner back height 16" plus 11" with headrest extension. This selection comes standard with rear anti-tippers. This selection must be made when WD42 and 17STF are selected.
- ANGLE Adjustable Angle Back with 10° Bend \$100.00
 Recliner and angle back must use conventional arms. Arms do not flip back on recliner. This selection comes standard with rear anti-tippers.

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1311BK Airline Buckle Seat Positioning Strap \$50.00
- 1515 Push Button Style Seat Positioning Strap \$50.00 ^{A9900}_{E0978}
- 1321BK Hook and Loop Seat Positioning Strap \$40.00

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap \$100.00 ^{E0960}

Invacare® IVC™ 9000 Jymni™ Wheelchair

CASTER AND WHEEL OPTIONS

CASTER OPTIONS

- 1687 6" x 1" Dark Grey Comp Urethane Casters Standard
- 1686 6" x 1" Composite Urethane Casters No Charge
- AT902 6" Composite Urethane 3 Spoke No Charge

WHEEL OPTIONS

- 1763BK 20" Composite Pneumatic Tires \$117.00 A9900
- 4763BK 20" Composite Pneumatic/Flat Free Inserts \$232.00 E2213
- U1763BK 20" Composite Urethane Tires Standard A9900
- 1764BK 22" Composite Pneumatic Tires \$117.00 A9900
- 4764BK 22" Composite Pneumatic/Flat Free Inserts \$232.00 E2213
- U1764BK 22" Composite Urethane Tires \$113.00 A9900

HANDRIM OPTIONS

- BHDRM Black Handrims Standard

ACCESSORIES

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$10.00
- 9120 Solid Seat Insert \$98.00 A9900

BACK ACCESSORIES

- 1483JY Stroller Handle \$91.00
Not available on Recliner or Angle back type.
- 9130 Solid Back Insert \$107.00
Only available in 14" back height.
- L8601 Omit Back Upholstery Less \$10.00
- 1483R Headrest w/Push Handle Extensions \$95.00

MISCELLANEOUS ACCESSORIES

- 9357JY 1" Forward Positioner \$120.00
This comes standard with rear anti-tippers. Not available with Recliner or Angle back type or 10" and 12" deep chairs.
- A1576 Wheel Lock Extensions \$39.00
- 9691 Push to Lock Wheel Locks No Charge
Pull-to-Lock wheel locks are standard.
- 1577QR Quick Release Axles \$86.00 K0108
- 1487 IV/Oxygen Holder Package Assembly \$158.00
Not available on Recliner or Angle back type.
- 1488 Oxygen Holder Package Assembly \$151.00
Not available on Recliner or Angle back type.
- 1360 Adjustable Rear Anti-Tipper \$82.00
Anti-Tippers for 22" rear wheels.
- 1360B Adjustable Rear Anti-Tippers- 20" Wheels \$82.00
Anti-Tippers for 20" rear wheels.

COLOR OPTIONS

FRAME FINISH

- 104P Electric Teal No Charge
- 105P Midnight Blue No Charge
- 08P Deep Purple No Charge
- 24P Wet Black No Charge
- 34P Satin (Flat) Black No Charge
- 61P Electric Red No Charge
- 62P Electric Blue No Charge
- 66P Electric Purple No Charge
- 71P Silver Vein No Charge

UPHOLSTERY COLORS

- U240 Black Nylon Upholstery Standard
Recliner is only available with U240 black nylon upholstery. Recliner back height 16" plus 11" with headrest extension.
- U550 Midnight Blue Vinyl Upholstery No Charge
- U67 Black Vinyl Upholstery No Charge

Invacare® IVC™ 9000 Jymni™ Wheelchair

FRONT RIGGING OPTIONS

FRONT RIGGINGS

- J932A 70° Footrests with Aluminum Footplates \$190.00
- J932A-14 80° Footrests with Alum Ftplts - 12" and 14" deep ... \$190.00
Not Available on 10" deep frames.
- J933A 90° Footrests with Aluminum Footplates \$190.00
- J934A Elevating Legrests with Aluminum Footplates \$270.00
Knee-to-Heel range is required with J934A
- J934A-14 Elevating Legrests with Alum Ftplts - 14" deep \$270.00
Not Available on 10" and 12" deep frames.

FOOTREST/LEGREST ACCESSORIES

- 1820 Non-Removable Device \$12.00 A9900
- 1346 Heel Loops \$39.00 E0951
Standard on Swingaway Footrests.

KNEE-TO-HEEL RANGE

Knee-to-Heel Range must be selected if J934A is chosen.

- SHEXT 7" - 10" (Short Extended) No Charge
- MDEXT 9" - 12" (Medium Extended) No Charge
- LGEXT 11" - 14" (Long Extended) No Charge

Seat to Floor Matrix 6" Casters			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
15"	Top	20"	Top
16"	Middle	20"	Bottom
17"	Bottom	22"	Bottom

Weight Limitation 125 Pounds
 Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Comments

