

Invacare® IVC™ 9000 XT Recliner Wheelchair

Price List and Order Form



TRANSACTION TYPE:

Order Quote

Date _____ Account # _____ P.O. # _____
 Company Name _____ Phone # _____
 Address _____ City _____ State _____ Zip Code _____ Country _____
 Contact Name _____ Comments _____

Special Note

The IVC 9000XT Recliner is standard with attendant operated wheel locks

PRE-CONFIGURED STOCK WHEELCHAIRS

STOCK CHAIRS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	9RC	9RC/WD87/R/28/2R/34P/AT903/U2222/COM/U240 9RC 9000XT Recliner (9153638254)	\$1,759.49

STOCK OPTIONS

Quantity	Model Number	Description	List Price
<input type="checkbox"/>	1100 Tool Kit	1100 Tool Kit (9153629686)	No Charge

FRAME AND SEAT SELECTIONS

FRAME TYPE

9RC 9000 XT Recliner \$1,759.49

FRAME WIDTH AND DEPTH

- WD46 14" x 16" No Charge A9900
- WD48 14" x 18" \$126.38 A9900
- WD66 16" x 16" No Charge
- WD67 16" x 17" No Charge A9900
- WD68 16" x 18" \$126.38 A9900
- WD86 18" x 16" No Charge
- WD87 18" x 17" No Charge A9900
- WD88 18" x 18" \$126.38 A9900
- WD06 20" x 16" \$152.08 E2201
- WD07 20" x 17" \$272.03
- WD08 20" x 18" \$272.03
- WD26 22" x 16" \$152.08
- WD27 22" x 17" \$272.03
- WD28 22" x 18" \$272.03 E2201

ARM AND BACK OPTIONS

ARM TYPE

Fixed Height Arms 9", Adjustable Height Arms 9"-13", removable style.

- 18 Full Length Fxd Ht Conv No Charge A9900
- 19 Full Length Adj Ht Conv \$134.94 A9900
- 28 Desk Length Fxd Ht Conv No Charge A9900
- 29 Desk Length Adj Ht Conv \$134.94 E0973
- 19JY Full Length Adj Ht Conv (8" - 10") \$139.23 E0973
- 29JY Desk Length Adj Ht Conv (8" - 10") \$139.23 E0973

BACK TYPE

Add an additional 11" for the headrest extension.

- 2R Reclining Back 23" High No Charge
- L Reclining Back 23" High with Lift Pkg \$447.68 E1226

SEAT-TO-FLOOR HEIGHT

R Recliner No Charge

**All prices subject to change without notice.
 NOTE: All specifications and dimensions are approximate.**

Invacare® IVC™ 9000 XT Recliner Wheelchair

COLOR OPTIONS

FRAME FINISH

- 17P Cobalt BlueNo Charge
- 61P Electric RedNo Charge
- 24P Wet BlackNo Charge
- 34P Satin (flat) BlackStandard
- 71P Silver VeinNo Charge

UPHOLSTERY COLORS

All back upholstery is screwed on.

- U240 Black Nylon Upholstery Standard
- U245 Black Sail Cloth Adj Tension Seat \$10.71
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon upholstery only.
- U67 Black Vinyl No Charge
- U550 Midnight Blue Vinyl No Charge
- U86 Ruby Red Vinyl No Charge

CASTERS

CASTER OPTIONS

- 1228 8" x 1.75" Semi Pneu Casters Light Grey \$59.40 A9900
- AT903 8" x 1.25" Comp Urethane Standard
- 1290 8" x 1" Solid Rubber Casters Charcoal No Charge
- 1292 8" x 1.25" Pneu Casters \$93.18 A9900
- 1492 8" x 1.25" Pneu with Flat Free Inserts \$104.00 E2217
- 1299 8" x 1" Solid Rubber Casters Dark GreyStandard

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite Pneu Tires \$125.31 A9900
- 4763BK 20" Composite Pneu. Flat Free Inserts \$248.47 E2213
- U1763BK 20" Composite Urethane Tires \$113.00 A9900
- 1764BK 22" Composite Pneu Tires \$125.31 A9900
- 4764BK 22" Composite Pneu. Flat Free Inserts \$248.47 E2213
- U1764BK 22" Composite Urethane Tires \$113.00 A9900
- 2244BK 24" Composite PneumaticTires \$140.00 A9900
- 4244BK 24" Composite Pneu. Flat Free Inserts \$225.00 E2213
- U2244BK 24" Comp Urethane Treaded Tires \$120.00 A9900
- U2222C 24" Comp Urethane Corded Tires Standard

HANDRIM OPTIONS

- AHDRM Aluminum HandrimsNo Charge
Not available with U2222C rear wheels.
- COM Composite Handrims Standard
Available with 24" rear wheels only.
- 1512 Plastic Coated Aluminum Handrims \$91.04 A9900
- 1224-8 8 Alum Oblique Proj Handrims, Plst Coated ... \$269.89 A9900
Available with 24" rear wheels only.
- 1225-8 8 Alum Vert Proj Hndrms, Plst Coated \$269.89 A9900
Available with 24" rear wheels only.

WHEEL LOCKS

- Attendant Operated Standard
Attendant operated wheel locks are standard and not available with patient operated wheel locks.
- 9696R Patient Operated Push To Lock (Recliner) \$96.39 A9900
- 1329BK Toggle Wheel Lock Extension (8") \$41.77
- 1245 Patient Operated Pull to Lock (Recliner) No Charge

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap \$53.55 A9900 E0978
- 1321BK Hook and Loop Seat Positioning strap \$40.00 E0978

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap \$53.55 E0960

ACCESSORIES AND STOCK OPTIONS

BACK ACCESSORIES

- 1160 Carrying Pocket \$57.83 A9270
- 3046 Head Support For Recliner \$117.00 E0955
- STC Stencil Back Upholstery \$29.00 E9270
3/4" block letters. Only available with vinyl upholstery.
- RB Dual Back Liner \$52.48 E9270

SEAT ACCESSORIES

- 9120 Solid Seat Insert \$98.00 A9900
- RS Dual Seat Liner \$52.48 E9270

ARM ACCESSORIES

- NRAS Non - Removable Arm With Screws \$10.00 A9900
- 1300 Utility Tray \$168.15 E0950

MISCELLANEOUS ACCESSORIES

- 2219 IV Rod, Telescopic Recliner \$117.81
- 2226 Front Anti - Tippers \$83.00
- 1100 Tool Kit No Charge

All prices subject to change without notice.

NOTE: All specifications and dimensions are approximate.

Invacare® IVC™ 9000 XT Recliner Wheelchair

FRONT RIGGING OPTIONS

FRONT RIGGING

- T93HA Hemi Sw. Frst Alum Ftplt \$93.00
- T93HC Hemi Sw. Frst Comp Ftplt \$93.00
- T94HA Hemi Sw. Elevating Legrests Alum Footplate .. \$180.00
- T94HC Hemi Sw. Elevating Legrests Comp Ftplt \$180.00
- T93HS Hemi Sw. Frst Aluminum Ftplts 14" w \$90.00 A9900
- T94HS Hemi Sw. Elevating Legrests Aluminum 14" .. \$175.00 K0195
Only available on 14" wide chair.
- AHL4A Hemi Smartleg® Art. Legrests 16" and above .. \$341.09
- AHL4C Hemi Smartleg® Art Legrests 16" and Above .. \$341.09
- AHS4A Hemi Smartleg Art. Legrest 14" - 15" wide \$341.09
Only available on 14" wide chair.
- H93M Hemi Minimum Swingaway Footrests \$185.00 A9900
- H93E Hemi Ext Swingaway Footrests \$185.00 A9900
- T93HAB Hemi Bolt - Thru Footrests, Alum Footplates \$90.00 A9900
- T94HAB Hemi Bolt - Thru Legrests, Alum Footplates .. \$175.00 A9900

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates \$79.00 K0041
- AT5543 Adjustable Angle Flip - Up Footplates \$67.00

FOOTREST / LEGREST ACCESSORIES

- 1820 Non - Removable Device \$12.00 A9900
- 1337 Calf Strap \$27.00 K0038
Not available with AHL4A or AHL4C elevating legrests.
- 1348 Legrest Panel \$66.00 K0108
- 1346 Heel Loops \$39.00 E0951
Standard on footrests.
- L8689 Omit Heel Loops Less \$4.00
- 1319 Black Non - Padded Calf Pads \$10.00 A9900
Not available with AHL4A or AHL4C elevating legrests.

Seat - to - Floor Matrix 8" Casters			
Seat to Floor Height	Front Caster Mounting Position	Rear Wheel Size	Rear Wheel Mounting Position
17 1/2"	Top	20"	1 Position Only
18 1/2"	Middle	22"	1 Position Only
19 1/2"	Bottom	24"	1 Position Only

Weight Limitation 250 Pounds
Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.

Comments

All prices subject to change without notice.
NOTE: All specifications and dimensions are approximate.

Invacare Corporation

©2020 Invacare Corporation
All rights reserved. Trademarks are identified by the symbols™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications and prices are subject to change without notification.

