



Invacare® IVC™ 9000 XT Recliner Wheelchair

Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.invacare.com/eForms



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682



Yes, you can.®

Account Information

Request Type: Quote Order

Date: _____

Purchase Order #: _____

Account #: _____

CONTACT

Company: _____

Name: _____

SHIP TO

Back Up Contact: _____

Name: _____

Phone: _____

Address: _____

Email: _____

City: _____ State: _____ Zip: _____

Comments: _____

Country: _____

Special Note

The IVC 9000XT Recliner is standard with attendant operated wheel locks

FRAME AND SEAT SELECTIONS

FRAME TYPE

9RC 9000 XT Recliner\$1,759.49 K0004 + E1226

FRAME WIDTH AND DEPTH

WD46 14" x 16" No Charge E1236
14" Seat width changes base code to E1236

WD88 18" x 18" \$126.38

WD48 14" x 18" \$126.38 E1236
14" Seat width changes base code to E1236

WD06 20" x 16" \$152.08 E2201

WD66 16" x 16" No Charge

WD07 20" x 17" \$272.03

WD67 16" x 17" No Charge

WD08 20" x 18" \$272.03 E2201

WD68 16" x 18" \$126.38

WD26 22" x 16" \$152.08 E2201

WD86 18" x 16" No Charge

WD27 22" x 17" \$272.03 E2201

WD87 18" x 17" No Charge

WD28 22" x 18" \$272.03 E2201



Yes, you can.®

SEAT-TO-FLOOR HEIGHT

R ReclinerNo Charge

ARM AND BACK OPTIONS

ARM TYPE

Fixed Height Arms 9", Adjustable Height Arms 9"-13", removable style.

18 Full Length Fxd Ht Conv No Charge
19 Full Length Adj Ht Conv\$134.95 E0973
28 Desk Length Fxd Ht ConvNo Charge

29 Desk Length Adj Ht Conv\$134.95 E0973
19JY Full Length Adj Ht Conv (8" - 10") \$139.23 E0973
29JY Desk Length Adj Ht Conv (8" - 10") ... \$139.23 E0973

BACK TYPE

Add an additional 11" for the headrest extension.

2R Reclining Back 23" High No Charge

Lift package for Reclining Back 23" High ... \$447.68 K0108

COLOR OPTIONS

FRAME FINISH

17P Cobalt Blue No Charge
61P Electric Red No Charge
24P Wet BlackNo Charge

34P Satin (flat) Black Standard
71P Silver Vein No Charge

UPHOLSTERY COLORS

All back upholstery is screwed on.

U240 Black Nylon Upholstery Standard
U245 Black Sail Cloth Adj Tension Seat \$10.71
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon upholstery only.
U67 Black Vinyl No Charge

U550 Midnight Blue Vinyl No Charge
U86 Ruby Red Vinyl No Charge

CASTERS

CASTER OPTIONS

1228 8" x 1.75" Semi Pneu Casters Light Grey \$59.40
AT903 8" x 1.25" Comp Urethane Standard
1290 8" x 1" Solid Rubber Casters CharcoalNo Charge
1292 8" x 1.25" Pneu Casters\$93.18 E2214+E2215

1492 8" x 1.25" Pneu with Flat Free Inserts
.....\$104.00 E2217
1299 8" x 1" Solid Rubber Casters Dark Grey Standard



Yes, you can.®

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite Pneu Tires\$125.31 E2211 + E2212
- 4763BK 20" Composite Pneu. Flat Free Inserts
.....\$248.47 E2211 + E2213
- U1763BK 20" Composite Urethane Tires\$113.00
- 1764BK 22" Composite Pneu Tires
.....\$125.31 E2211 + E2212
- 4764BK 22" Composite Pneu. Flat Free Inserts
.....\$248.47 E2211 + E2213

- U1764BK 22" Composite Urethane Tires
.....\$113.00
- 2244BK 24" Composite Pneumatic Tires
.....\$140.00 E2211 + E2212
- 4244BK 24" Composite Pneu. Flat Free Inserts
.....\$225.00 E2211 + E2213
- U2244BK 24" Comp Urethane Treaded Tires\$120.00
- U2222C 24" Comp Urethane Corded Tires Standard

HANDRIM OPTIONS

- AHDRM Aluminum HandrimsNo Charge
Not available with U2222C rear wheels.
- COM Composite Handrims Standard
Available with 24" rear wheels only.
- 1512 Plastic Coated Aluminum Handrims \$91.04

- 1224-8 8 Alum Oblique Proj Handrims, Plst Coated
.....\$269.89
Available with 24" rear wheels only.
- 1225-8 8 Alum Vert Proj Hndrms, Plst Coated\$269.89
Available with 24" rear wheels only.

WHEEL LOCKS

- Attendant Operated Standard
Attendant operated wheel locks are standard and not available with patient operated wheel locks.
- 9696R Patient Operated Push To Lock (Recliner)
.....\$96.39

- 1329BK Toggle Wheel Lock Extension (8")\$41.77 E0961
- 1245 Patient Operated Pull to Lock (Recliner)
..... No Charge

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap
.....\$53.55 E0978

- 1321BK Hook and Loop Seat Positioning Strap
.....\$40.00 E0978

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap
.....\$53.55 E0960



Yes, you can.®

ACCESSORIES AND STOCK OPTIONS

BACK ACCESSORIES

- 1160 Carrying Pocket \$57.83
- 3046 Head Support For Recliner \$117.00 E0966

- STC Stencil Back Upholstery \$29.00
3/4" block letters. Only available with vinyl upholstery.
- RB Dual Back Liner \$52.48

SEAT ACCESSORIES

- 9120 Solid Seat Insert \$98.00

- RS Dual Seat Liner \$52.48

ARM ACCESSORIES

- NRAS Non - Removable Arm With Screws \$10.00

- 1300 Utility Tray \$168.15 E0950

MISCELLANEOUS ACCESSORIES

- 2219 IV Rod, Telescopic Recliner \$117.81 K0105
- 2226 Front Anti - Tippers \$83.00

- 1100 Tool Kit No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

- T93HA Hemi Sw. Footrests Aluminum Footplates w/ Heel Loops \$93.00
- T93HC Hemi Sw. Footrests Composite Footplates w/ Heel Loops \$93.00
- T94HA Hemi Sw. Elev Legrests Aluminum Footplates & Padded Calf Pad \$180.00 K0195
- T94HC Hemi Sw. Elev Legrests Composite Footplates & Padded Calf Pad \$180.00 K0195
- T93HS Hemi Sw. Footrests Aluminum Footplates w/ Heel Loops- 14"W Chairs \$90.00
- T94HS Hemi Sw. Elevating Legrests Aluminum Footplates & Padded Calf Pad- 14" W Chairs \$175.00 K0195
- AHL4A Hemi Smart Leg Articulating Legrests Aluminum Footplates \$341.09 K0053

- AHL4C Hemi Smart Leg Articulating Legrests Composite Footplates \$341.09 K0053
- AHS4A Hemi Smart Leg Articulating Legrestst- 14-15" W Chairs \$341.09 K0053
- H93M Hemi Minimum Swingaway Footrests \$185.00
- H93E Hemi Ext Swingaway Footrests \$185.00
- T93HAB Hemi Bolt-Thru Footrests, Aluminum Footplates \$90.00
- T93HAB Hemi Bolt-Thru Elevating Legrests, Aluminum Footplates \$175.00

FOOTPLATE OPTIONS

- 1350 Extra Large Aluminum Footplates \$79.00 K0041

- AT5543 Adjustable Angle Flip - Up Footplates \$67.00 K0040



Yes, you can.®

FOOTREST / LEGREST ACCESSORIES

- 1820 Non - Removable Device\$12.00
- 1337 Calf Strap\$27.00 K0038
Not available with AHL4A or AHL4C elevating legrests.
- 1348 Legrest Panel\$66.00

- 1346 Heel Loops\$39.00 E0951
Standard on footrests.
- L8689 Omit Heel LoopsLess \$4.00
- 1319 Black Non - Padded Calf Pads\$10.00
Not available with AHL4A or AHL4C elevating legrests.

Seat to Floor Matrix 8" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2"	Top	20"	1 Position Only
18 1/2"	Middle	22"	1 Position Only
19 1/2"	Bottom	24"	1 Position Only

Weight Limitation 250 Pounds

Invacare recommends all wheelchairs to be fitted to the consumer under the guidance of a healthcare professional knowledgeable of the particular limitations of the product.



Yes, you can.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682