



Invacare® IVC™ 9000 XT Wheelchair

Thank you for using Invacare eForms.

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.invacare.com/eForms



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682



Yes, you can.®

Account Information

Request Type: Quote Order

Date: _____

Account #: _____

Company: _____

SHIP TO
Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Purchase Order #: _____

CONTACT
Name: _____

Back Up Contact: _____

Phone: _____

Email: _____

Comments: _____

Special Note

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.

Existing information (serial number) will need to be supplied to Customer Service to ensure correct configuration.

PRE-CONFIGURED STOCK WHEELCHAIRS

STOCK CHAIRS

9XT_PTO_29154 9000XT 16" x 16" Adult, Desk Length
Fixed Ht Space Saver Arm\$1,355.19 K0004

9XT_PTO_29153 9000XT 18" x 16" Adult, Desk Length
Fixed Ht Space Saver Arm\$1,355.19 K0004

9XT_PTO_29155 9000XT 18" x 16" Adult, Full Length
Fixed Ht Space Saver Arm\$1,355.19 K0004

9XT_PTO_33131 9000XT 20" x 16" Adult, Desk Length
Fixed Ht Space Saver Arm\$1,355.19 K0004+E2201

STOCK FRONT RIGGINGS

T93HAP Hemi Sw. Footrests Aluminum Footplates
w/ Heel Loops \$99.13 E0951

T93HCP Hemi Sw. Footrests Composite Footplates
w/ Heel Loops \$99.13 E0951

T94HAP Hemi Sw. Elev Legrests Aluminum Footplates
& Padded Calf Pad \$191.86 K0195

T94HCP Hemi Sw. Elev Legrests Composite Footplates
& Padded Calf Pad \$191.86 K0195

AHL4A/U67 Hemi Smart Leg Articulating Legrests
Aluminum Footplates \$341.09 K0053

STOCK OPTIONS

1497 Oxygen Holder Pkg Assy \$161.72 E2208

1496 IV / Oxygen Hldr Pkg Assy \$169.22 K0105+E2208

1360 Adjustable Rear Anti-Tippers- 8" Caster
..... \$87.82 E0971

1360A Adjustable Rear Anti-Tippers- 6" Caster.....
..... \$87.82 E0971

A1576 Wheel Lock Extensions \$41.77 E0961



Yes, you can.®

FRAME AND SEAT SELECTIONS

FRAME TYPE

9XT 9000XT\$1,355.19 K0004

FRAME WIDTH AND DEPTH

WD46 14" x 16" No Charge E1236
14" Seat width changes base code to E1236

WD48 14" x 18" \$126.38 E1236
14" Seat width changes base code to E1236

WD56 15" x 16"\$126.38

WD58 15" x 18"\$126.38

WD66 16" x 16"No Charge

WD67 16" x 17"\$126.38

WD68 16" x 18"\$126.38

WD76 17" x 16"\$126.38

WD77 17" x 17"\$126.38

WD78 17" x 18"\$126.38

WD86 18" x 16"No Charge

WD87 18" x 17"\$126.38

WD88 18" x 18"\$126.38

WD96 19" x 16"\$152.08
Not available with Super Low seat-to-floor.

WD97 19" x 17" \$272.03
Not available with Super Low seat-to-floor.

WD98 19" x 18" \$272.03
Not available with Super Low seat-to-floor.

WD06 20" x 16" \$152.08 E2201
Not available with Super Low seat-to-floor.

WD07 20" x 17" \$272.03 E2201
Not available with Super Low seat-to-floor.

WD08 20" x 18" \$272.03 E2201
Not available with Super Low seat-to-floor.

WD26 22" x 16" \$152.08 E2201
Not available with Super Low seat-to-floor.

WD27 22" x 17" \$272.03 E2201
Not available with Super Low seat-to-floor.

WD28 22" x 18" \$272.03 E2201
Not available with Super Low seat-to-floor.

SEAT-TO-FLOOR HEIGHT

A Adult (19 1/2" STF)No Charge
Based on 24" wheels and 8" casters. See seat-to-floor matrix.

H Hemi (17 1/2" STF)No Charge
Based on 24" wheels and 8" casters. See seat-to-floor matrix.

SH Super Low (14 1/2" STF)\$91.04 K0056

SHE Super Hemi (15 1/2" STF)No Charge K0056
Based on 20" wheels and 6" casters. See seat-to-floor matrix.
Requires AT902 (6" casters), U1763BK (20" composite wheels, urethane tires) AHDRM (aluminum handrims) and A1576 (wheel lock extensions).
Seat-to-floor height is approximately 15 1/2".

CASTERS

CASTER OPTIONS

AT902 6" Composite Urethane 3 SpokeNo Charge

1228 8" x 1 3/4" Semi. Pneu Casters Lt Grey \$59.98
Not available with Super Low seat-to-floor.

AT903 8" x 1 1/4" Composite UrethaneStandard
Not available with Super Low seat-to-floor.

1299 8" x 1" Solid Rubber CastersNo Charge
Not available with Super Low seat-to-floor.

1292 8" x 1 1/4" Pneu Casters\$93.18 E2214+E2215
Not available with Super Low seat-to-floor.

1492 8" x 1 1/4" Pneu. with Flat Free Inserts ...\$111.38 E2217
Not available with Super Low seat-to-floor.

CASTER JOURNAL OPTION

Not available with Adjustable Angle or Adjustable Angle-S Back Type.
Not available with Super Low seat-to-floor.
Fixed caster journals are standard.

ADJC Adjustable Caster Journal \$31.06

ADJCPB Adjustable Precision Bearing Caster Journal
..... \$43.91



Yes, you can.®

ARM AND BACK OPTIONS

ARM TYPE

Arm types are removable, flip-back.

- 18 Full Length Fixed Ht Conv (10")No Charge
- 19 Full Length Adj. Ht. Conv (10" - 14") \$134.95 E0973
- 28 Desk Length Fixed Ht. Conv (10")No Charge
- 29 Desk Length Adj. Ht. Conv (10" - 14") \$134.95 E0973
- 78 Full Length Fixed Ht. SS (10")No Charge
Not available with adjustable angle back type.
- 79 Full Length Adj. Ht. SS (10" - 14") \$134.95 E0973
Not available with adjustable angle back type.

- 88 Desk Length Fixed Ht. SS (10")No Charge
Not available with adjustable angle back type.
- 89 Desk Length Adj. Ht. SS (10" - 14") \$134.95 E0973
Not available with adjustable angle back type.
- 79JY Full Length Adj. Ht. SS (8" - 10")\$139.23 E0973
Not available with adjustable angle back type.
- 89JY Desk Length Adj. Ht. SS (8" - 10")\$139.23 E0973
Not available with adjustable angle back type.

BACK TYPE

- ADJ Adj. Back with 10° BendNo Charge
- S1418 Straight Back w/ Push Handles (15" - 18")
.....No Charge

- ANGLE Adj. Angle Back with 10° Bend\$107.10 K0108
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.
- ANGLE-S Adj. Angle Back with Straight\$107.10 K0108
Not available with Super Low seat-to-floor. Standard with rear anti-tippers.

COLOR OPTIONS

FRAME FINISH

- 08P Deep PurpleNo Charge
- 17P Cobalt BlueNo Charge
- 24P Wet BlackNo Charge
- 34P Satin (Flat) BlackNo Charge

- 71P Silver VeinNo Charge
- 88P Wet Black with Twilight SparkleNo Charge
- 105P Midnight BlueNo Charge

CANDY COLORS

- 61P Electric RedNo Charge
- 62P Electric BlueNo Charge

- 66P Electric PurpleNo Charge
- 104P Electric TealNo Charge

UPHOLSTERY COLORS

- U67 Black VinylNo Charge
Vinyl upholstery must choose WRAP or SOUPH.
- U240 Black Nylon UpholsteryStandard

- U245 Black Sail Cloth Adj. Tension Seat \$10.71
U245 sail cloth adjustable tension seat upholstery available with U240 black nylon back upholstery only.

NAUGAHYDE

- WRAP Wrap Around Upholstery\$82.47

- SOUPH Screwed On Back Upholstery\$95.32
Not available with Adjustable Angle or Adjustable Angle-S Back Type.



Yes, you can.®

WHEELS

WHEEL OPTIONS

- 1763BK 20" Composite Pneu. Tires\$125.31 E2211 + E2212
- 4763BK 20" Composite Pneu. Flat Free Inserts
.....\$248.47 E2213 + E2211
- U1763BK 20" Composite Urethane Tires\$121.02
- 1764BK 22" Composite Pneu. Tires\$125.31 E2211 + E2212
- 4764BK 22" Composite Pneu Flat Free Inserts
.....\$248.47 E2211+E2213
- 5244BK 24" HD Comp. Pneu. with FF Inserts
.....\$319.16 E2211+E2213
Not available with Super Low seat-to-floor.
- U1764BK 22" Composite Urethane Tires\$121.02

- 2244BK 24" Composite Pneu. Tires ...\$149.94 E2211 + E2212
Not available with Super Low seat-to-floor.
- 4244BK 24" Composite Pneu. with FF Inserts
.....\$240.98 E2211+E2213
Not available with Super Low seat-to-floor.
- U2244BK 24" Composite Urethane Treaded Tires
.....\$128.52
Not available with Super Low seat-to-floor.
- U2222C 24" Composite Urethane Corded Tires
.....Standard
Not available with Super Low seat-to-floor.

WHEEL OPTIONS (ONE ARM DRIVE - OAD)

The Cyclical Lever drive requires treaded rear wheels and is only available on 16", 17", and 18" deep. Hemi seat-to-floor must use space saver arms. Weight limit is 250 lbs. The above options are not available with amputee attachment or offset axle.

- 9456L LH One Arm Drive with U2222C\$649.03 E0958
- 9456R RH One Arm Drive with U2222C\$649.03 E0958
- LDR96R1 6" Right Hand Cyclical Lever One Arm Drive ...
.....\$853.59 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.
- LDR96L1 6" Left Hand Cyclical Lever One Arm Drive
.....\$853.59 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.

- LDR98R1 8" Right Hand Cyclical Lever One Arm Drive
.....\$853.59 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.
- LDR98L1 8" Left Hand Cyclical Lever One Arm Drive
.....\$853.59 E0958
Not available with Super Low seat-to-floor.
Not available with adjustable caster journal.

HANDRIM OPTIONS

- AHDRM Aluminum HandrimsNo Charge
Not available with U2222C rear wheels.
- COM Composite HandrimsStandard
Not available on 20" and 22" composite rear wheels.
- 1512 Plstc Coated Aluminum Handrims\$91.04

- 1224-8 8 Alum Oblique Proj. Hdrm, Plstc Coat ...\$269.89
Not available on 20" and 22" composite rear wheels.
- 1225-8 8 Alum Vert. Proj. Hdrm, Plstc Coat\$269.89
Not available on 20" and 22" composite rear wheels.

REAR AXLE

- OFFAXL Offset Axle\$42.84 K0108

Not available with Adjustable Angle or Adjustable Angle-S Back Type.
Not available with Super Low seat-to-floor.



Yes, you can.®

WHEEL LOCKS

The above selections are not available with Adjustable Angle or Adjustable Angle-S Back Type.

- 9690B Pull To Lock Whl Locks (Std on Hemi) No Charge
- 9691B Push To Lock Whl Locks (Std On Adult) No Charge

- A1567 Pull To lock Clamp On Whl Locks \$81.40
- 1721 High Push To Lock Clamp On Whl Locks \$76.04

POSITIONING STRAP OPTIONS

SEAT POSITIONING STRAPS

- 1515 Push Button Style Seat Positioning Strap \$53.55 E0978

- 1321BK Hook and Loop Seat Positioning Strap \$42.84 E0978

CHEST POSITIONING STRAPS

- 9515 Push Button Style Chest Positioning Strap \$53.55 E0960

ACCESSORIES AND STOCK OPTIONS

BACK ACCESSORIES

- L8601 Omit Back Upholstery Less \$0.00
- ADJTENS Adjustable Tension Back Upholstery \$119.95
- 1160 Carrying Pocket \$57.83

- STC Stencil Back Upholstery \$31.06
Only available with vinyl upholstery. Block letters are .75" in height.
- RB Dual Back Liner \$49.00
- BPO2 Back Pack \$53.55

SEAT ACCESSORIES

- L8651 Omit Seat Upholstery Less \$0.00
- 9120 Solid Seat Insert \$104.96
Not available on 22" wide chairs.

- RS Dual Seat Liner \$52.48

ARM ACCESSORIES

- NRAS Non-Removable Arm With Screws \$10.71

- 1300 Utility Tray \$168.15 E0950
Not available with 24" rear wheel set in Hemi position.

MISCELLANEOUS ACCESSORIES

- 1305L Crutch and Cane Carrier \$46.05 E2207
Not available with traditional one arm drive, but is available with CLD one arm drive.
- 1577QR Quick Release Axles \$92.11
Not available with offset axle.
Not available with any one arm drive.
- 2218 IV Rod, Telescopic \$160.65 K0105
Not available with traditional one arm drive, but is available with CLD one arm drive.

- 2226 Front Anti - Tippers \$88.89 E0971
- 9357SH Super Low Amputee Attachment . \$100.78 E0959
Only available on Super Low Seat-to-Floor.
- 1360SH Super Low Adj Rear Anti - Tipper \$87.82 E0971
Only available on Super Low Seat-to-Floor.
- 9357L Amputee Attachment \$100.78 E0959
Not available with offset axle.
Not available with any one arm drive.



STOCK OPTIONS

1497 Oxygen Holder Package Assembly \$161.72 E2208
Not available with traditional one arm drive, but is available with CLD one arm drive.
Not available with offset axle.
For 0° and 6° seat angles. See service manual for details.

1496 IV / Oxygen Holder Package Assembly
..... \$169.22 K0105+E2208
Not available with offset axle.
Not available with any one arm drive.

1360 Adjustable Rear Anti-Tippers- 8" Caster
..... \$87.82 E0971
Only available with 8" casters.
Standard with Adjustable Angle Back.

1360A Adjustable Rear Anti - Tipper (6" Casters)
..... \$87.82 E0971
Only available with 6" casters.

A1576 Wheel Lock Extensions \$41.77 E0961
Elastic cord wheel lock extensions. Hemi seat to floor come standard with non-cord wheel lock extensions.

1100 Tool Kit No Charge

FRONT RIGGING OPTIONS

FRONT RIGGING

T93HA Hemi Sw. Footrests Aluminum Footplates
w/ Heel Loops \$99.13 E0951

T93HC Hemi Sw. Footrests Composite Footplates
w/ Heel Loops \$99.13 E0951

T94HA Hemi Sw. Elev Legrests Aluminum Footplates
& Padded Calf Pad \$191.86 K0195

T94HC Hemi Sw. Elev Legrests Composite Footplates
& Padded Calf Pad \$191.86 K0195

SH904M Super Low Sw. Elevating Legrests 11"-17"
..... \$266.00 K0195

SH904MS Super Low Sw. Elevating Leg rests 14"W
..... \$266.00 K0195

SH93E Super Low Sw. Ext. Footrests 11" - 17" \$185.00

SH93ES Super Low Sw. Ext. Footrests 14"W \$185.00

AHL4A Hemi Smart Leg Articulating Legrests
Aluminum Footplates \$320.00 K0053

AHL4C Hemi Smart Leg Articulating Legrests
Composite Footplates \$320.00 K0053

AHS4A Hemi Smart Leg Articulating Legrestst-14-15"
W Chairs \$320.00 K0053

T93HS Hemi Sw. Footrests Aluminum Footplates
w/ Heel Loops- 14"W Chairs \$95.93

T94HS Hemi Sw. Elevating Legrests Aluminum
Footplates & Padded Calf Pad- 14"W Chairs
..... \$186.53 K0195

T963A 60° Sw. Footrests Aluminum Footplates ... \$95.93

T963C 60° Sw. Footrests Composite Footplates.. \$95.93

T963S 60° Sw. Footrests Composite Footplate 14"W
..... \$95.93 K0041

FOOTPLATE OPTIONS

1350 Extra Large Aluminum Footplates \$84.21 K0041

AT5543 Adjustable Angle Flip - up Footplates
..... \$71.42 K0040

FOOTREST / LEGREST ACCESSORIES

1820 Non - Removable Device \$12.79

1337 Calf Strap \$38.37 K0038

1346 Heel Loops \$41.57 E0951

IMPACTG Impact Guards \$54.36




Yes, you can.®

Seat-to-Floor 6" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
15 1/2"	Top	20"	Top
16 1/2"	Middle	22"	Top
17 1/2"	Bottom	24"	Top
17 1/2"	Bottom	20"	Bottom

Seat-to-Floor 8" Casters

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION
17 1/2"	Top	20"	Bottom
17 1/2"	Top	24"	Top
18 1/2"	Middle	22"	Bottom
19 1/2"	Bottom	24"	Bottom



Weight Limit 250 lbs

Invacare recommends all wheelchairs be fitted to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Seat-to-Floor 6" Casters and Offset Axle Bracket

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION	AXLE PLATE MOUNTING ON FRAME	
				Top Lug	Bottom Lug
15 1/2"	Top Hole	22"	#5 Hole	#4 Hole	#1 Hole
15 1/2"	Top Hole	20"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	24"	#5 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	22"	#3 Hole	#4 Hole	#1 Hole
16 1/2"	Middle Hole	20"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	24"	#3 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Bottom Hole	20"	#3 Hole	#4 Hole	#1 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.

Seat-to-Floor 8" Casters and Offset Axle Bracket

SEAT-TO-FLOOR HEIGHT	FRONT CASTER MOUNTING POSITION	REAR WHEEL SIZE	REAR WHEEL MOUNTING POSITION	AXLE PLATE MOUNTING ON FRAME	
				Top Lug	Bottom Lug
17 1/2"	Top Hole	20"	#3 Hole	#1 Hole	#4 Hole
17 1/2"	Top Hole	22"	#2 Hole	#4 Hole	#1 Hole
17 1/2"	Top Hole	24"	#3 Hole	#4 Hole	#1 Hole
18 1/2"	Middle Hole	22"	#3 Hole	#1 Hole	#4 Hole
18 1/2"	Middle Hole	24"	#2 Hole	#4 Hole	#1 Hole
19 1/2"	Bottom Hole	20"	#3 Hole	#1 Hole	#4 Hole

For a complete list of seat-to-floor heights and angles that are configurable by the dealer refer to the service manual.



Yes, you can.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-678-4682