

INVACARE® MATRX® SEATING SERIES

INVACARE® MATRX® VI CUSHION

Skin Protection and Positioning Cushion

Actual cushion width is 1/2" narrower than stated size
Actual cushion depth is 1 1/2" longer than stated size due to rear cushion radius

<input type="radio"/>	MA-VI	Vi Cushion	\$425.00	E2607 E2608
		Use this model number to configure the Matrx Vi cushion		
<input type="checkbox"/>	MA1415-Vi	Vi Cushion 14x15	\$425.00	E2607
<input type="checkbox"/>	MA1417-Vi	Vi Cushion 14x17	\$425.00	E2607
<input type="checkbox"/>	MA1616-Vi	Vi Cushion 16x16	\$425.00	E2607
<input type="checkbox"/>	MA1618-Vi	Vi Cushion 16x18	\$425.00	E2607
<input type="checkbox"/>	MA1620-Vi	Vi Cushion 16x20	\$425.00	E2607
<input type="checkbox"/>	MA1716-Vi	Vi Cushion 17x16	\$425.00	E2607
<input type="checkbox"/>	MA1718-Vi	Vi Cushion 17x18	\$425.00	E2607
<input type="checkbox"/>	MA1720-Vi	Vi Cushion 17x20	\$425.00	E2607
<input type="checkbox"/>	MA1816-Vi	Vi Cushion 18x16	\$425.00	E2607
<input type="checkbox"/>	MA1818-Vi	Vi Cushion 18x18	\$425.00	E2607
<input type="checkbox"/>	MA1820-Vi	Vi Cushion 18x20	\$425.00	E2607
<input type="checkbox"/>	MA2016-Vi	Vi Cushion 20x16	\$425.00	E2607
<input type="checkbox"/>	MA2018-Vi	Vi Cushion 20x18	\$425.00	E2607
<input type="checkbox"/>	MA2020-Vi	Vi Cushion 20x20	\$425.00	E2607
<input type="checkbox"/>	MA2218-Vi	Vi Cushion 22x18	\$475.00	E2608
<input type="checkbox"/>	MA2220-Vi	Vi Cushion 22x20	\$475.00	E2608
<input type="checkbox"/>	CV(size)	Vi Cushion Outer Cover	\$100.00	E2619

INVACARE® MATRX® VI BASE CUSHION

<input type="radio"/>	MDB-VIBASE	Matrx-Vi Base	\$200.00	E2231
-----------------------	------------	---------------------	----------	-------

NOTE

To order Invacare Matrx Vi Base, please select MDB-VI BASE within the MA-VI configurator. Matrx Vi Bases include plywood drop base and mounting hardware. Must pick drop hook sizes 7/8" or 1" diameter; 1" or 2" drop, or flushmount.

Specify Drop Base Mount:

Diameter	7/8" <input type="checkbox"/>	1" <input type="checkbox"/>	
Front Set	2" <input type="checkbox"/>	1" <input type="checkbox"/>	Flush(0") <input type="checkbox"/>
Rear Set	2" <input type="checkbox"/>	1" <input type="checkbox"/>	Flush(0") <input type="checkbox"/>

MATRX VI CUSHION BUILDER SIZES

SEAT WIDTH

<input type="radio"/>	W14	14" Wide	No Charge
<input type="radio"/>	W15	15" Wide	No Charge
		Invacare Matrx Vi Cushion Base Only	
<input type="radio"/>	W16	16" Wide	No Charge
<input type="radio"/>	W17	17" Wide	No Charge
		Invacare Matrx Vi Cushion Base Only	
<input type="radio"/>	W18	18" Wide	No Charge
<input type="radio"/>	W19	19" Wide	No Charge
		Invacare Matrx Vi Cushion Base Only	
<input type="radio"/>	W20	20" Wide	No Charge
<input type="radio"/>	W22	22" Wide	No Charge

SEAT DEPTH

<input type="radio"/>	D15	15" Deep	No Charge
<input type="radio"/>	D16	16" Deep	No Charge
<input type="radio"/>	D17	17" Deep	No Charge
<input type="radio"/>	D18	18" Deep	No Charge
<input type="radio"/>	D20	20" Deep	No Charge

INVACARE® MATRX® KIDABRA VI CUSHION

Skin Protection and Positioning Cushion

Actual cushion width is 1/2" narrower than stated size
Actual cushion depth is 1 1/2" longer than stated size due to rear cushion radius

<input type="radio"/>	KA-VI	Kidabra Vi Cushion	\$425.00	E2607
		Use this model to configure the Kidabra Vi cushion		
<input type="checkbox"/>	KA1212Vi	Kidabra Cushion 12x12	\$425.00	E2607
<input type="checkbox"/>	KA1214Vi	Kidabra Cushion 12x14	\$425.00	E2607
<input type="checkbox"/>	KA1412Vi	Kidabra Cushion 14x12	\$425.00	E2607
<input type="checkbox"/>	KA1414Vi	Kidabra Cushion 14x14	\$425.00	E2607
<input type="checkbox"/>	KC(size)	Kidabra Cushion Outer Cover	\$95.00	E2619

INVACARE® MATRX® KIDABRA VI BASE CUSHION

<input type="radio"/>	KDB-VIBASE	Kidabra Vi Cushion Base	\$200.00	E2231
-----------------------	------------	-------------------------------	----------	-------

NOTE

To order Invacare KidabraVi Base, please select KDB-VIBASE within the KA-VI configurator. Kidabra Vi Bases include plywood drop base and mounting hardware. Must pick drop hook sizes 7/8" or 1" diameter; 1" or 2" drop, or flushmount.

Specify Drop Base Mount:

Diameter	7/8" <input type="checkbox"/>	1" <input type="checkbox"/>	
Front Set	2" <input type="checkbox"/>	1" <input type="checkbox"/>	Flush(0") <input type="checkbox"/>
Rear Set	2" <input type="checkbox"/>	1" <input type="checkbox"/>	Flush(0") <input type="checkbox"/>

MATRX KIDABRA VI CUSHION BUILDER SIZES

SEAT WIDTH

<input type="radio"/>	W12	12" Wide	No Charge
<input type="radio"/>	W13	13" Wide	No Charge
		Invacare Kidabra Vi Cushion Base Only	
<input type="radio"/>	W14	14" Wide	No Charge

SEAT DEPTH

<input type="radio"/>	D12	12" Deep	No Charge
<input type="radio"/>	D14	14" Deep	No Charge

NOTE: All specifications and dimensions are approximate.