

Invacare® Solara® 3G SmartSelect Wheelchair Configurations

CHART 1: CONFIGURATION OPTIONS					
MODEL NUMBER	SOLARA3GLS1	SOLARA3GLS2	SOLARA3GLS3	SOLARA3GLS4	SOLARA3GLS5
Frame Style	Tilt in Space	Tilt in Space	Tilt in Space	Tilt in Space	Tilt in Space
Transport Selection	TRBKTS	TRBKTS	TRBKTS	TRBKTS	TRBKTS
Tilting Option	Cable Free Foot Release Right Mount	Cable Free Foot Release Right Mount	Cable Free Foot Release Right Mount	Traditional Hand Release	Traditional Hand Release
Seat Width Range	14-18"	16-20"	18-22"	16-20"	18-22"
Seat Width	16"	18"	20"	18"	20"
Seat Depth Range	15-19"	15-19"	18-22"	15-19"	18-22"
Seat Depth	16"	16"	19"	18"	19"
Seat to Floor Height	17"	17"	17"	18"	18"
Frame Color	Wet Black	Wet Black	Wet Black	Electric Red	Electric Blue
Caster	6 x 2" Semi Pneumatic	6 x 2" Semi Pneumatic	6 x 2" Semi Pneumatic	8 x 1" Urethane	8 x 1" Urethane
Rear Wheel	12" Composite	12" Composite	12" Composite	24" Composite Aluminum Handrims	24" Composite Aluminum Handrims
Axle Type	Quick Release	Quick Release	Quick Release	Quick Release	Quick Release
Rear Tire	Flat Free Inserts	Flat Free Inserts	Flat Free Inserts	Flat Free Inserts	Flat Free Inserts
Wheel Locks	Foot Operated Push to Lock	Push to Lock	Push to Lock	Push to Lock	Push to Lock
Seat Pan	Flush	Flush	Flush	Flush	Flush
Arm Height	High	High	Low	High	High
Arm Style	Conventional Full Length Two Point Adjustable Height Arms	T-Arms Adjustable Height Full Length	Flip Back Full Length Two Point Adjustable Arms	Two Point Fixed Height Desk Length Arms	Conventional Desk Length Two Point Adjustable Height Arms
Back Type	Standard Aluminum Back Cane Fold Down	Standard Aluminum Back Cane Fold Down	Standard Aluminum Back Cane Fold Down	Standard Aluminum Back Cane Fold Down	Standard Aluminum Back Cane Fold Down
Back Height	20" (Straight Cane)	24" (Straight Cane)	24" (Straight Cane)	24" (Straight Cane)	24" (Straight Cane)
Back Accessories	Adj. Angle Stroller Handle Inside Mount	None	None	Adjustable Angle Stroller Handle Inside Mount	Adjustable Angle Stroller Handle Outside Mount
Front Riggings	70° Non Tapered Swingaway	70° Non Tapered Swingaway	70° Non Tapered Swingaway	70° Swingaway MFX	90° Non Tapered Swingaway
Positioning Accessories	Push Button Buckle Seat Strap Heel Loops	Push Button Buckle Seat Strap Heel Loops	Push Button Buckle Seat Strap Heel Loops	Push Button Buckle Seat Strap Heel Loops	Push Button Buckle Seat Strap Heel Loops
Misc. Accessories	Rear Anti-tippers	Rear Anti-tippers	Rear Anti-tippers	Rear Anti-tippers	Rear Anti-tippers