



JAY® Extremity Supports

January 2024

Mark For: _____
 Date: _____
 Dealer Acct #: _____
 Dealer: _____
 Dealer Contact: _____
 Dealer Address: _____
 Dealer City: _____ ST: _____ ZIP: _____
 Dealer Phone: () _____ Fax: () _____
 Confirmation Email: _____
 Confirm Via: Fax Email

Submitting for: Quote Order
 PO#: _____
 Ship To: _____
 Attention: _____
 Address: _____
 Address: _____
 Ship To City: _____ ST: _____ ZIP: _____
 Ship To Phone: () _____ Fax: () _____

ADDITIONAL SHIPPING INFORMATION


The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

Hand and Arm Supports

Available on Dual Post and Dual Post Height Adjustable. Requires 3/4" round tubing for mounting the hardware.

Arm Trough with Fixed Hardware

✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JARM14F	E2209	Arm Trough 14" w/Fx Mnt Hdw	\$234
<input type="checkbox"/>	JARM18F	E2209	Arm Trough 18" w/ Fx Mnt Hdw	\$234
<input type="checkbox"/>	JARM22F	E2209	Arm Trough 22" w/ Fx Mnt Hdw	\$234



Fixed Hardware



Swing Away

Arm Trough with Swing Away Hardware

<input type="checkbox"/>	JARM14SA		Arm Trough 14" L w/ SA Hdw	\$447
	JARM14	E2209	Arm Trough 14" Long	\$161
	SAARMHDW	E1028	Swing Away Hardware	\$286
<input type="checkbox"/>	JARM18SA		Arm Trough 18" L w/ SA Hdw	\$447
	JARM18	E2209	Arm Trough 18" Long	\$161
	SAARMHDW	E1028	Swing Away Hardware	\$286
<input type="checkbox"/>	JARM22SA		Arm Trough 22" L w/ SA Hdw	\$447
	JARM22	E2209	Arm Trough 22" Long	\$161
	SAARMHDW	E1028	Swing Away Hardware	\$286

Hardware		
A	Fixed	Swing-Away
14"	JARM14F	JARM14SA
18"	JARM18F	JARM18SA
22"	JARM22F	JARM22SA

Wrist Straps

<input type="checkbox"/>	ATS7510	K0108	Wrist Straps - Small	\$77
<input type="checkbox"/>	ATS7520	K0108	Wrist Straps - Medium	\$77
<input type="checkbox"/>	ATS7530	K0108	Wrist Straps - Large	\$77
<input type="checkbox"/>	ATS7590	K0108	Wrist Straps - Made To Order	\$172



Wrist Circumference		
S	ATS7510	5.5" - 7"
M	ATS7520	7.5" - 9"
L	ATS7530	9.5" - 11"
MTO	ATS7590	Specify

Gel Pad Covers

<input type="checkbox"/>	131066	K0108	Gel Pad Cover - Desk Length	\$71
<input type="checkbox"/>	131067	K0108	Gel Pad Cover - Full Length	\$71
<input type="checkbox"/>	131071	K0108	Gel Pad Cover - Recline Right	\$77
<input type="checkbox"/>	131072	K0108	Gel Pad Cover - Recline Left	\$77



Armrest Gel Pad Cover

JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts. Please refer to the "JAY YOUR WAY" product catalog for complete listing of modifications

Arm Trough Modifications

<input type="checkbox"/>	Add Pommel to Hand Pad	\$28	<input type="checkbox"/>	Reduce height of side walls	NC
<input type="checkbox"/>	Foam overlay - add visco foam	\$28	<input type="checkbox"/>	Reduce Height by 1/2"	<input type="checkbox"/>
<input type="checkbox"/>	Add additional padded strap	\$18	<input type="checkbox"/>	Remove side walls	NC

Arm Trough Colored Side Band Options

<input type="checkbox"/>	Cover Options						NC						
<input type="checkbox"/>	Electric Lime	<input type="checkbox"/>	Grey	<input type="checkbox"/>	Hunter Orange	<input type="checkbox"/>	Navy Blue	<input type="checkbox"/>	Pink Sorbet	<input type="checkbox"/>	Shocking Pink	<input type="checkbox"/>	Turquoise
<input type="checkbox"/>	University Red	<input type="checkbox"/>	Very Violet	<input type="checkbox"/>	Volt Yellow	<input type="checkbox"/>	Digital Camo						

Arm Trough Contact Surface Options

<input type="checkbox"/>	Micro Climatic Contact Surface	NC	<input type="checkbox"/>	Naugahyde	NC
<input type="checkbox"/>	Sheep Skin	NC			

Lower Extremity Support

JAY® Lower Extremity Support

<input type="checkbox"/>	161565	E1020	Large Pad 10"L x 8"W - Right	\$646
<input type="checkbox"/>	161566	E1020	Large Pad 10"L x 8"W - Left	\$646
<input type="checkbox"/>	161567	E1020	Std Pad 8"L x 6"W - Right	\$646
<input type="checkbox"/>	161568	E1020	Std Pad 8"L x 6"W - Left	\$646



Knee Adductors

JAY® Knee Adductors with Pads

<input type="checkbox"/>	KAK21XXL	Knee Adductor Short w/ Pad - L		\$406
		<input type="checkbox"/> 3" Pad	<input type="checkbox"/> 4" Pad	
<input type="checkbox"/>	KAK21XXR	Knee Adductor Short w/ Pad - R		\$406
		<input type="checkbox"/> 3" Pad	<input type="checkbox"/> 4" Pad	
<input type="checkbox"/>	KAK31XXL	Knee Adductor Long w/ Pad - L		\$406
		<input type="checkbox"/> 3" Pad	<input type="checkbox"/> 4" Pad	
<input type="checkbox"/>	KAK31XXR	Knee Adductor Long w/ Pad - R		\$406
		<input type="checkbox"/> 3" Pad	<input type="checkbox"/> 4" Pad	



JAY® Knee Adductors without Pads

<input type="checkbox"/>	HDW4853	Knee Adductor Short - Left	\$441
<input type="checkbox"/>	HDW4854	Knee Adductor Short - Right	\$441
<input type="checkbox"/>	HDW4851	Knee Adductor Long - Left	\$441
<input type="checkbox"/>	HDW4852	Knee Adductor Long - Right	\$441

JAY® Knee Adductors Pads Only

<input type="checkbox"/>	KAD10XX	3" Round Pad	\$62
<input type="checkbox"/>	KAD11XX	4" Round Pad	\$62
<input type="checkbox"/>	KAD19XX	Made To Order Pad	\$134

Foot Positioning

	Calculating Foot Box Width
Flared Frame	Foot Box Width = Seat Width + 2"
Offset Hanger	Foot Box Width = Seat Width + 2"
0' Offset	Foot Box Width = Seat Width
1' Offset	Foot Box Width = Seat Width - 2"
2' Offset	Foot Box Width = Seat Width - 4"

Single Foot Box

<input type="checkbox"/>	JFOOTSS	E0954	Footbox Single Std Ht	\$258
<input type="checkbox"/>	JFOOTST	E0954	Footbox Single Tall Ht	\$258



Single Foot Box			
	Height	A	B
JFOOTSS	Standard	8"	7.5"
JFOOTST	Tall	12"	7.5"

Full Foot Box

<input type="checkbox"/>	JFOOT12S	E0954	Full Foot Box 12"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT14S	E0954	Full Foot Box 14"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT16S	E0954	Full Foot Box 16"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT18S	E0954	Full Foot Box 18"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT20S	E0954	Full Foot Box 20"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT22S	E0954	Full Foot Box 22"W Std Ht.	\$459
<input type="checkbox"/>	JFOOT16T	E0954	Full Foot Box 16"W Tall Ht.	\$459
<input type="checkbox"/>	JFOOT18T	E0954	Full Foot Box 18"W Tall Ht.	\$459
<input type="checkbox"/>	JFOOT20T	E0954	Full Foot Box 20"W Tall Ht.	\$459
<input type="checkbox"/>	JFOOT22T	E0954	Full Foot Box 22"W Tall Ht.	\$459

	Full Foot Box			
	Height	Width	A	B
JFOOT12S	Standard	12"	8"	9.5"
JFOOT14S	Standard	14"	8"	11.5"
JFOOT16S	Standard	16"	8"	13.5"
JFOOT18S	Standard	18"	8"	15.5"
JFOOT20S	Standard	20"	8"	17.5"
JFOOT22S	Standard	22"	8"	19.5"
JFOOT16T	Tall	16"	12"	13.5"
JFOOT18T	Tall	18"	12"	15.5"
JFOOT20T	Tall	20"	12"	17.5"
JFOOT22T	Tall	22"	12"	19.5"




Foot Positioning (Continued)

Full Foot Box with Divider

<input type="checkbox"/>	JFOOT16SD	E0954	Full Footbox 16"W Std Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT18SD	E0954	Full Footbox 18"W Std Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT20SD	E0954	Full Footbox 20"W Std Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT22SD	E0954	Full Footbox 22"W Std Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT16TD	E0954	Full Footbox 16"W Tall Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT18TD	E0954	Full Footbox 18"W Tall Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT20TD	E0954	Full Footbox 20"W Tall Ht. w/ Div.	\$515
<input type="checkbox"/>	JFOOT22TD	E0954	Full Footbox 22"W Tall Ht. w/ Div.	\$515

Full Foot Box with Divider				
	Height	Width	A	B
JFOOT16SD	Standard	16"	8"	13.5"
JFOOT18SD	Standard	18"	8"	15.5"
JFOOT20SD	Standard	20"	8"	17.5"
JFOOT22SD	Standard	22"	8"	19.5"
JFOOT16TD	Standard	16"	12"	13.5"
JFOOT18TD	Standard	18"	12"	15.5"
JFOOT20TD	Tall	20"	12"	17.5"
JFOOT22TD	Tall	22"	12"	19.5"



JAY® Ankle Positioning System

<input type="checkbox"/>	TST6501L	Ankle Positioning X-Small - L	\$81
<input type="checkbox"/>	TST6501R	Ankle Positioning X-Small - R	\$81
<input type="checkbox"/>	TST6511L	Ankle Positioning Small - Left	\$81
<input type="checkbox"/>	TST6511R	Ankle Positioning Small - Right	\$81
<input type="checkbox"/>	TST6521L	Ankle Positioning Medium - L	\$81
<input type="checkbox"/>	TST6521R	Ankle Positioning Medium - R	\$81
<input type="checkbox"/>	TST6541L	Ankle Positioning Large - Left	\$81
<input type="checkbox"/>	TST6541R	Ankle Positioning Large - Right	\$81
<input type="checkbox"/>	TST6561L	Ankle Positioning X-Large - L	\$81
<input type="checkbox"/>	TST6561R	Ankle Positioning X-Large - R	\$81
<input type="checkbox"/>	TST6591L	Ankle Positioning Made To Order - L	\$188
<input type="checkbox"/>	TST6591R	Ankle Positioning Made To Order - R	\$188



Ankle Positioning System

JAY® Foot & Ankle Positioning System

<input type="checkbox"/>	TST6510L	Foot & Ankle Pos Small - Left	\$104
<input type="checkbox"/>	TST6510R	Foot & Ankle Pos Small - Right	\$104
<input type="checkbox"/>	TST6520L	Foot & Ankle Pos Medium - L	\$104
<input type="checkbox"/>	TST6520R	Foot & Ankle Pos Medium - R	\$104
<input type="checkbox"/>	TST6540L	Foot & Ankle Pos Large - Left	\$104
<input type="checkbox"/>	TST6540R	Foot & Ankle Pos Large - Right	\$104
<input type="checkbox"/>	TST6560L	Foot & Ankle Pos X-Large - L	\$104
<input type="checkbox"/>	TST6560R	Foot & Ankle Pos X-Large - R	\$104



Foot & Ankle Positioning

JAY® Padded Ankle / Toe Strap Only

<input type="checkbox"/>	TST4000	E0952	Padded Ankle/Toe Strap - S	\$33
<input type="checkbox"/>	TST4200	E0952	Padded Ankle/Toe Strap - M	\$33
<input type="checkbox"/>	TST4500	E0952	Padded Ankle/Toe Strap - L	\$33



Padded Ankle Toe Strap Only

JAY® Flexi-Foot System

<input type="checkbox"/>	32600XS-RT	E0951/E0952	Flexi-Foot System - X-Small	\$261
<input type="checkbox"/>	32600S-RT	E0951/E0952	Flexi-Foot System - Small	\$261
<input type="checkbox"/>	32600M-RT	E0951/E0952	Flexi-Foot System - Medium	\$261
<input type="checkbox"/>	32600L-RT	E0951/E0952	Flexi-Foot System - Large	\$261

Flexi-Foot			
	Item#	A	B
XS	32600XS-RT	4.5"	2.75"
S	32600S-RT	6.25"	3"
M	32600M-RT	8.25"	3.5"
L	32600LS-RT	9"	4.5"



JAY® Shoe Holders

<input type="checkbox"/>	TST6000	E0951/E0952	Shoe Holder - Small	\$73
<input type="checkbox"/>	TST6100	E0951/E0952	Shoe Holder - Medium	\$73
<input type="checkbox"/>	TST6200	E0951/E0952	Shoe Holder - Large	\$73
<input type="checkbox"/>	TST6300	E0951/E0952	Shoe Holder - X-Large	\$73



Shoe Holder



Shoe Holder with Padded Straps

JAY® Shoe Holder with Padded Strap

<input type="checkbox"/>	TST6001	E0951/E0952	Shoe Holder w/ Pad Strap - S	\$81
<input type="checkbox"/>	TST6101	E0951/E0952	Shoe Holder w/ Pad Strap - M	\$81
<input type="checkbox"/>	TST6201	E0951/E0952	Shoe Holder w/ Pad Strap - L	\$81
<input type="checkbox"/>	TST6301	E0951/E0952	Shoe Holder w/ Pad Strap - XL	\$81

Foot Positioning (Continued)

JAY® Shoe Holder with Tendon Relief

<input type="checkbox"/>	TST6061	E0951/E0952	Shoe Holder w/ Tendon Relief - S	\$374
<input type="checkbox"/>	TST6161	E0951/E0952	Shoe Holder w/ Tendon Relief - M	\$379
<input type="checkbox"/>	TST6261	E0951/E0952	Shoe Holder w/ Tendon Relief - L	\$384
<input type="checkbox"/>	TST6361	E0951/E0952	Shoe Holder w/ Tendon Relief - XL	\$399




Shoe Holder with Tendon Relief

JA® Padded Calf Panel

<input type="checkbox"/>	PLS3010	E0951/E0952	Padded Calf Panel - Small	\$148
<input type="checkbox"/>	PLS3020	E0951/E0952	Padded Calf Panel - Medium	\$148
<input type="checkbox"/>	PLS3030	E0951/E0952	Padded Calf Panel - Large	\$148
<input type="checkbox"/>	PLS3099	E0951/E0952	Padded Calf Panel - Made To Order	\$314

Padded Calf Panel



JAY YOUR WAY MODIFICATIONS

Note: To order sizes offered via JAY Your Way add a "M" to the front of the part number and follow the prompts.
Please refer to the "JAY YOUR WAY" product catalog for complete listing of modifications

JAY Foot Boxes Modifications

<input type="checkbox"/>	Add Leg Strap	\$18	<input type="checkbox"/>	Reduce Leg Through Size	NC
<input type="checkbox"/>	Foam overlay - add visco foam	\$28	<input type="checkbox"/>	Reduce Opening by 1/2"	<input type="checkbox"/>
			<input type="checkbox"/>	Reduce Opening by 1"	
			<input type="checkbox"/>	Increase Leg Through Size	NC
			<input type="checkbox"/>	Increase Opening by 1/2"	<input type="checkbox"/>
			<input type="checkbox"/>	Increase Opening by 1"	

JAY Foot Boxes Colored Side Band Options

<input type="checkbox"/>	Cover Options						NC						
<input type="checkbox"/>	Electric Lime	<input type="checkbox"/>	Grey	<input type="checkbox"/>	Hunter Orange	<input type="checkbox"/>	Navy Blue	<input type="checkbox"/>	Pink Sorbet	<input type="checkbox"/>	Shocking Pink	<input type="checkbox"/>	Turquoise
<input type="checkbox"/>	University Red	<input type="checkbox"/>	Very Violet	<input type="checkbox"/>	Volt Yellow	<input type="checkbox"/>	Digital Camo						

JAY Foot Boxes Contact Surface Options

<input type="checkbox"/>	Micro Climatic Contact Surface	NC	<input type="checkbox"/>	Naughahyde	NC
<input type="checkbox"/>	Sheep Skin	NC			

