

Seating

Standing

Sleeping

Therapy

JENX[®]
 Positioning for life


Multistander Size 1&2 by Jenx

NEW Multistander Size 1 & 2 features

- New entry level Head Support
- New white frame
- Angle-adjustable Footplate as standard
- Antimicrobial agent in all PU surfaces
- Improved knee support bracket
- New fastenings to allow use of a Jenx Waistcoat Harness



Multistander Size 1

Shown at the smallest height.

NEW Multistander Size 2 features

- Multigrip™ Body Support
- Growth range from 3-13 years!
- Max user weight up to 132 pounds.
- Unlike other standers, very large overlap, no more confusion!



Now Available!
Multistander Size 2

Shown at the largest height.



distributors of



INNOVATIVE FEATURE **Smooth & Simple Angle Adjustment**
 Angle can be adjusted with the child in the standing frame, and is easily set using the angle gauge on the frame. Frame can be tilted to be completely horizontal, aiding transfers when used in supine. Concentrate on safe transfers first then positioning second.

This highly versatile standing frame is suitable for children from around **9 Months** right up to **13 Years** across 2 sizes and provides highly supportive standing in prone, supine or upright.

An unrivalled amount of adjustment for height and width and a maximum user weight of 132LBS on the Size 2, make the Multistander the perfect choice for growing children in multi-use environments. Smooth and easy angle adjustment takes the Multistander from horizontal right through to upright, facilitating easy transfers in supine and unlimited scope for setting the precise standing angle required for every child.

The Multistander has a wide variety of head, trunk, leg and foot support options to enable you to exactly customize the support to any individual user. The deliberately open design ensures you can always see precisely how the user is positioned each and every time. The capacity for outstanding thoracic and pelvic support, coupled with highly adjustable leg positioning options, make the Multistander a great option for almost every child. If you need flexibility the Jenx Multistander has got it!

Jenx Multistander Technical Specifications

	Units	MUS01 Size 1 Prone	MUS01 Size 1 Supine	MUS02 Size 2 Prone	MUS02 Size 2 Supine
Angle Range	degrees	10°-90°	90°-10°	10°-90° (Powered base) 0°-80° (Manual base)	10°-90° (Powered base) 0°-80° (Manual base)
Age Range (approx)	mths/hrs	9mths-6yrs	9mths-6yrs	3yrs-13yrs	3yrs-13yrs
Max User Height	inch	46"	46"	69"	69"
Frame Width	inch	22¼"	22¼"	25½"	25½"
Frame Length	inch	28"	28"	39"	39"
Minimum Height (footplate to top of upper support board)	inch	19½"	22½"	27½"	30½"
Maximum Height (footplate to top of upper support board)	inch	34"	34"	46"	46"
Minimum Height (footplate to top of head support)	inch	-	26"	-	33"
Maximum Height (footplate to top of head support)	inch	-	45" h	-	53" (No shoulder support) 59" (w/ shoulder support)
Chest Pad Width (min/max)	inch	6" - 10"	6" - 10"	6¾" - 12½"	6¾" - 12½"
Hip Pad Width (min/max)	inch	4¾" - 10¼"	4¾" - 10¼"	6¾" - 12½"	6¾" - 12½"
Kneeblock Width (centre to centre) (min/max)	inch	4¾" - 10¼"	4¾" - 10¼"	6¾" - 12½"	6¾" - 12½"
Height Of Stander (floor to top of upper pad) (min/max)	inch	23¾" - 35½"	26" - 35½"	23¾" - 35½"	26" - 35½"
Footplate Angle Range	degrees	-15°-+15°	-15°-+15°	-15°-+15°	-15°-+15°

Multistander by Jenx

A highly versatile standing frame suitable for children from around 9 months right up to 6 years, standing in prone, supine or upright.

An unrivaled amount of adjustment for height and width and a maximum user weight of 66 lbs make the Multistander the perfect choice for growing children. Smooth and easy to operate angle adjustment takes the Multistander from horizontal right through to a natural standing position of just off upright. A wide variety of head, trunk, leg and foot support options allow you to exactly customize the Multistander to any individual user. The capacity for outstanding thoracic and pelvic support, coupled with highly adjustable leg positioning options, make the Multistander a great option for almost every child. If you need flexibility, the Jenx Multistander has got it!

INNOVATIVE FEATURE

Independently Adjustable Thoracic and Pelvic Boards

Allowing the individual proportions of each specific child to be catered for, while also giving scope to avoid pressure on delicate areas such as feeding buttons.

INNOVATIVE FEATURE

Smooth & Simple Angle Adjustment

Angle can be adjusted with the child in the standing frame, and is easily set using the angle gauge on the frame.

INNOVATIVE FEATURE

Knee blocks

Available in two sizes and provided with specially designed straps, the Multistander knee blocks are independently adjustable for height, width, depth and angle allowing you to precisely position the knees of every child.

INNOVATIVE FEATURE

Incorporated de-rotational Thoracic and Pelvic Straps with Lateral Support Pads

Available in different sizes and styles these support pads offer lateral stability and come with a support strap which enables corrective positioning and security.

INNOVATIVE FEATURE

PU Support Pads

Flexible, easy to clean, impervious to fluids and with an antimicrobial agent, the PU boards and support pads make the Multistander simple to clean and keep hygienic. They also come in a range of vibrant colors to suit all tastes.

INNOVATIVE FEATURE

Innovative Prone/Supine Adjustable Tray

Thoughtfully designed to work effectively in both prone and supine, this universal tray addresses many of the issues reported when using 3:1 standing frames. Clever inbuilt adjustments prevent gaps for arm traps while the included overlay ensures a perfect fit when used in prone.





Supine with tray



Supine angle adjustment



Prone with tray

Multistander Color Range



Multistander

The appropriate option must be selected for both the trunk and pelvis



Multistander Frame

Comes complete with:

- Base with 4 lockable swivel casters
- PU chest & hip boards
- Knee block mounting brackets*
- Angle adjustable footplate

*NB. Knee blocks must be specified separately. Available in Blue, Pink, Green, Red or Black boards

MUS01



Description	Code	Dimensions	Pelvis	Trunk	Prone	Supine
Small support pad w/ de-rotational strap	MB06 XX	3.7"(d) 3.5"(w) 3.9"(h)	✓	✓	✓	✓
Medium support pad w/de-rotational strap	MB07 XX	5.3"(d) 3.5"(w) 3.9"(h)	✓	✓	✓	✓
Thin thoracic pad w/ de-rotational strap	MB21 XX	6.5"(d) 0.8"(w) 3"(h)	✗	✓	✓	✓
High rise thoracics with de-rotational strap (pair)	MB02 XX	3.7"(d) 3.5"(w) 7.9"(h)	✗	✓	✓	✗

For use in Prone only

Always Configurable Accessories



Molded Head Support

Anatomically molded polyurethane head support for children who require minimal - moderate head support. Comes mounted on a flexible fixing bracket giving adjustment for height depth and angle.

MH07 XX

7.5"(w) x 1.4"(d) x 4"(h)



Multigrip Head Support (supine standing only)

The original and innovative Jenx Multigrip™ head support is now available in 2 sizes for all Jenx seating systems.

1 Link MH06 11.5"(w) x 2"(d) x 3.5"(h)

2 Link MH0101 18.7"(w) x 2"(d) x 4.7"(h)



Chest Prompt

For use in Prone*.

MB03

*Must be used with a below cover option:

Blue MU04

Red MU05

Pink MU06

Green MU07



Knee Support Pads

Innovative and flexible PU knee cups available in two sizes. Complete with padded knee straps, and are adjustable for height, width, angle and depth.

Size 1 MB09 XX

2.2"(d) x 3.1"(w) x 3.7"(h)

Size 2 MB04 XX

3.9"(d) x 3.7"(w) x 4.5"(h)



Sandal Raising Blocks (Set of 4)

The raising block kits give the sandals that extra bit of height where needed! Use 1 set to accommodate for leg length discrepancy or use 2 sets to raise the height of the footplate on any product to maximize the available growth.

NB: One kit needed per individual sandal.

Sandal Risers Size 1 SRB1 1 Riser = 0.7" 2 Risers = 1.3"
Sandal Risers Size 2 SRB2 3 Risers = 2" 4 Risers = 2.7"



Sandals (Pair)

Jenx sandals are interchangeable across the entire range of Jenx seating and standing products.

Size 1 S16

6.8"(d) x 3.9"(w)

Size 2 S2

8.6"(d) x 4.3"(w)



Multi-use Tray

Tray for use in prone or supine. Easy to attach and adjust for height, angle and depth, complete with a removable tray infill for use in prone standing.

MT02

25.6"(w) x 17.7"(d) x 1.4"(h)



Tray Elbow Blocks (Pair)

Used to encourage arms into a midline position and/or to block excessive retraction, made from gray polyurethane foam and come as a pair.

UBEB

6" (h) x 4"(d) x 3.5"(w)



Code: USJM/Jenx/2014

distributors of

