

BILLING/SHIPPING INFORMATION			
Date:			Name:
PO Number:			Address:
Branch:			
Mark for:			City:
Dealer Contact:			State:
Email Address:			Zip Code:
Phone:			Phone:
QUOTE		<input type="checkbox"/>	ORDER
		<input type="checkbox"/>	



Modern, child-centered styling, combined with a comprehensive range of support options make ATOM the ideal seating system to develop correct seating posture and function



KEY FEATURES:

- Adjustable flip-away hip support pads
- Tool-free adjustments
- Memory foam ischial well
- Hi Lo base offering tilt-in-space
- Wipe clean covers

Age range (approx.)	6 Months - 5 years
User weight (max.)	55 Lbs.
Back angle (posterior - anterior)	15° - 15°
Back pad height from seat to top (min - max)	11" - 15"
Chest width (min - max)	8" - 11"
Hip width (min - max)	7" - 11"
Seat depth (min - max)	6" - 13"
Lower leg length (min - max)	5" - 14"
Tilt in space angle range (posterior - anterior)	30° - 10°
Y Base footprint	31" (L) x 26" (W)
Seat height from floor (min - max)	11" - 30"
Suggested HCPCS Code E1399*	

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes are for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

CHOOSE YOUR COLOR: Please note all part numbers containing ** will color match



REQUIRED COMPONENTS

SEAT		<p style="text-align: center;">Atom seat unit</p> <p>Includes: Padded seat with memory foam Ischial well, shaped back pad, adjustable, flip away hip support pads, upper back pad and push handle</p>	
	AT01-05-**	\$1,849	<input checked="" type="checkbox"/>

BASE		<p style="text-align: center;">Y Base</p> <p>The most flexible base in the Jenx range, the Y-Base offers a height range of 11" - 30" alongside tilt-in-space functionality and lockable swivel casters. Gray</p>	
	JM01-05	\$3,955	<input checked="" type="checkbox"/>

FOOTPLATE		<p style="text-align: center;">Footplate</p> <p>Effortlessly flips to aid in transfers. Ideal for where plantar or dorsi flexion accommodation is required</p>	
	AT01 003	\$713	<input checked="" type="checkbox"/>

PELVIC STRAP <small>*One selection required</small>	<p style="text-align: center;">3 point Lap Strap</p>	<p>With additional safety of a groin harness. Ideal for smaller users</p>	
	HB63	\$155	<input type="checkbox"/>
	<p style="text-align: center;">4 point Lap Strap</p>	<p>Additional pelvic control with a more tailored position</p>	
	HB17	\$155	<input type="checkbox"/>

HEAD SUPPORT <small>*One selection required</small>			<p style="text-align: center;">Multigrip™ Headrest</p> <p>Provides height, depth and angle adjustment. Includes removable black fabric cover. Available with 1 (A+B) or 2 (A+B+C) link sets. Gray hardware</p>
	MH10-05 (1 Link, A-B)	\$329	<input type="checkbox"/>
	MH11-05 (2 Link, A-C)	\$477	<input type="checkbox"/>

OPTIONAL COMPONENTS

UPPER BODY SUPPORT		<p style="text-align: center;">High-rise Flexible Thoracic Support Pads Slim, flexible support pads offer straps, flip-away action and independent adjustability for height and width. PAIR</p>	
	A: 6" B: 3"	THC019	\$567
	A: 7" B: 3"	THC08	\$567
	A: 2" B: 5" C: 6"	SS00	\$368

	<p style="text-align: center;">Zip Waistcoat Harness</p> <p>Provides additional trunk support to children with low tone. Zipper allows for easy access to G tube site</p>	
A: 5" B: 5"	HB61	\$181
A: 6" B: 6"	HB62	\$181

	<p style="text-align: center;">Back Pad Infill</p> <p>Reduces the seat depth for smaller children. Attaches with Velcro behind the standard vinyl back pad</p>	
A: 7" B: 11" C: 2"	AT01 020	\$77

OPTIONAL COMPONENTS

LOWER BODY SUPPORT		<p style="text-align: center;">Leg Control Pads</p> <p>Ideal for controlling windsweeping or creating abduction, each leg is guttered in a soft adjustable cup, gray PAIR.</p> <p style="text-align: center;">A: 5" - 9" B: 4" - 9" C: 3"</p>		<p style="text-align: center;">Long Leg Sitting Board</p> <p>Provides a hamstring stretch or a larger platform of support for smaller children. Support blocks can be removed or repositioned</p> <p style="text-align: center;">A: 3" B: 4" C: 4" D: 15" E: 8"</p>																												
	AT01 016	\$534	<input type="checkbox"/>	<input type="checkbox"/>																												
		<p style="text-align: center;">Size 0 Abduction Block</p> <p>Creates abducted seating position. Easily removed/attached, Size 0. Black</p> <p style="text-align: center;">A: (Front) 4" (Back) 3" B: 4" C: 6" D: 3"</p>		<p style="text-align: center;">Size 1 Abduction Block</p> <p>Creates abducted seating position- easily removed/attached, Size 1. Black</p> <p style="text-align: center;">A: 3" B: 7" C: 6" D: 3"</p>																												
	GL24	\$151	<input type="checkbox"/>	<input type="checkbox"/>																												
		<p style="text-align: center;">Booties</p> <p>Provide gentle and comfortable support for tiny feet even without shoes. Soft and padded without buckles or straps, Black. PAIR</p> <p style="text-align: center;">A: 4" B: 6" C: 3" *approx</p>		<p style="text-align: center;">Footplate Pad</p> <p>Provides a comfortable resting position for feet when no foot positioning is chosen. Gray</p> <p style="text-align: center;">A: 11" B: 9" C: 1"</p>																												
	S5	\$380	<input type="checkbox"/>	<input type="checkbox"/>																												
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Dimensions</th> <th style="width: 15%;">Sandals</th> <th style="width: 15%;">Sandal Risers</th> <th style="width: 15%;">Sandals w/Ankle Snug</th> <th style="width: 15%;">Dimensions</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">Size 1 A: 3" B: 7"</td> <td style="text-align: center;"></td> <td style="text-align: center;"></td> <td style="text-align: center;"></td> <td style="text-align: center;">Size 1 A: 5" - 9" B: 5" C: 0" - 9"</td> </tr> <tr> <td style="text-align: center;">Size 2 A: 4" B: 9"</td> <td></td> <td></td> <td></td> <td style="text-align: center;">Size 2 A: 6.5" - 11.6" B: 4.7" C: 0" - 8.9"</td> </tr> </tbody> </table>	Dimensions	Sandals	Sandal Risers	Sandals w/Ankle Snug	Dimensions	Size 1 A: 3" B: 7"				Size 1 A: 5" - 9" B: 5" C: 0" - 9"	Size 2 A: 4" B: 9"				Size 2 A: 6.5" - 11.6" B: 4.7" C: 0" - 8.9"																
Dimensions	Sandals	Sandal Risers	Sandals w/Ankle Snug	Dimensions																												
Size 1 A: 3" B: 7"				Size 1 A: 5" - 9" B: 5" C: 0" - 9"																												
Size 2 A: 4" B: 9"				Size 2 A: 6.5" - 11.6" B: 4.7" C: 0" - 8.9"																												
<p style="text-align: center;">Sandals</p> <p>Ideal for positioning feet with a firm hold. Includes toe and heel straps. Black</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Size 1</td> <td style="width: 10%; text-align: center;">\$16</td> <td style="width: 10%; text-align: center;">\$178</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Size 2</td> <td style="text-align: center;">\$2</td> <td style="text-align: center;">\$203</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Size 1	\$16	\$178	<input type="checkbox"/>	<input type="checkbox"/>	Size 2	\$2	\$203	<input type="checkbox"/>	<input type="checkbox"/>	<p style="text-align: center;">Sandal Risers</p> <p>Contains 4 blocks measuring .7" each and designed for use with 1 foot. Black</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Size 1</td> <td style="width: 10%; text-align: center;">SRB1</td> <td style="width: 10%; text-align: center;">\$133</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Size 2</td> <td style="text-align: center;">SRB2</td> <td style="text-align: center;">\$158</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Size 1	SRB1	\$133	<input type="checkbox"/>	<input type="checkbox"/>	Size 2	SRB2	\$158	<input type="checkbox"/>	<input type="checkbox"/>	<p style="text-align: center;">Sandals w/Ankle Snug</p> <p>Offers unique magnetic buckle which can be operated single handed</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Size 1</td> <td style="width: 10%; text-align: center;">\$19</td> <td style="width: 10%; text-align: center;">\$387</td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> <td style="width: 10%; text-align: center;"><input type="checkbox"/></td> </tr> <tr> <td>Size 2</td> <td style="text-align: center;">\$21</td> <td style="text-align: center;">\$438</td> <td style="text-align: center;"><input type="checkbox"/></td> <td style="text-align: center;"><input type="checkbox"/></td> </tr> </table>	Size 1	\$19	\$387	<input type="checkbox"/>	<input type="checkbox"/>	Size 2	\$21	\$438	<input type="checkbox"/>	<input type="checkbox"/>
Size 1	\$16	\$178	<input type="checkbox"/>	<input type="checkbox"/>																												
Size 2	\$2	\$203	<input type="checkbox"/>	<input type="checkbox"/>																												
Size 1	SRB1	\$133	<input type="checkbox"/>	<input type="checkbox"/>																												
Size 2	SRB2	\$158	<input type="checkbox"/>	<input type="checkbox"/>																												
Size 1	\$19	\$387	<input type="checkbox"/>	<input type="checkbox"/>																												
Size 2	\$21	\$438	<input type="checkbox"/>	<input type="checkbox"/>																												

