

BILLING/SHIPPING INFORMATION					
Date:			Name:		
PO Number:			Address:		
Branch:					
Mark for:			City:		
Dealer Contact:			State:		
Email Address:			Zip Code:		
Phone:			Phone:		
		QUOTE	<input type="checkbox"/>	ORDER	<input type="checkbox"/>

**An adjustable floor based seating system
that provides support for young children
approx. 6 months to 10 years**



KEY FEATURES:

- The sturdy "wings" and corner design help to maintain symmetry of the user's spine
- Adjustable upper support can be used to support the head or shoulders
- A simple first seat to introduce children to long-leg
 - PU material is easy to clean
 - Removable abduction block



	Size 1	Size 2	Size 3
Age Range (approx)	6 mos. - 2 yrs	6 mos. - 6 yrs	4 yrs - 10 yrs
User Weight (max)	55 lbs	77 lbs	110 lbs
Total Back Support Height	12" - 17"	17" - 24"	20" - 27"
Lower Back Support Height	7"	9"	9"
Hip Width	13"	15"	20"
Seat Depth	14"	15"	18"
Footprint	20" x 20"	21" x 22"	24" x 26"

Suggested HCPCS Code - E1399*

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.*

Corner Seat includes:

- Adjustable, removable upper support
- Adjustable abduction block
- Chest strap
- Lap Strap
- Padded, PU surfaces



Available in Green

Size 1

NR45-03	\$1,271	<input type="checkbox"/>
---------	---------	--------------------------

Size 2

NR46-03	\$1,663	<input type="checkbox"/>
---------	---------	--------------------------

Size 3

NR47-03	\$2,333	<input type="checkbox"/>
---------	---------	--------------------------

OPTIONAL COMPONENTS

NURSERY TABLE	<p>Sturdy table with a large, angle adjustable surface area. A short strap on each side can attach to the Corner Seat which would prevent tipping. Natural</p>			
	<p>Size 0 A: 15" B: 24" C: 24" D: 20" E: 8"-11" Angle: 40°-40°</p>	NR64	\$877	<input type="checkbox"/>
	<p>Size 1 A: 15" B: 24" C: 24" D: 20" E: 10"-14" Angle: 40°-40°</p>	NR67	\$881	<input type="checkbox"/>
	<p>Size 2 A: 17" B: 28" C: 28" D: 20" E: 11"-17" Angle: 40°-40°</p>	NR68	\$976	<input type="checkbox"/>

TRIANGULAR CUSHIONS	<p>For those who prefer to have a flat back with wings. The triangular cushions can also be used beside the legs to prevent external rotation. Green (Pair)</p>			
		CS2001-03	\$179	<input type="checkbox"/>

HEADALOFT

Additional External Head Support Option

HeadAloft includes:

- Adjustable, removable, padded headrest with (2) buckle attachment straps-3 different sized, flip-away head support straps -Occipital support cushion



HEADALOFT

\$490