



# Jenx Multistander - Multi-position



## Retail Price List and Order Form 2022

Billing Information:		Ship To Information:	
Date:		Name:	
PO Number:		Address:	
Branch:			
Mark for:		City:	
Dealer Contact:		State:	
Email Address:		Zip Code:	
Phone:		Phone:	
Fax:			
		Quote	<input type="radio"/>
		Order	<input type="radio"/>

Size 1



Standard column and one piece footplate



Abduction column and individual footplates

Size 2



Standard column only

**Configurations for representation only**

	Size 1 standard	Size 1 ABD	Size 2
Approximate Age Range	9 mths - 6 yrs	9 mths - 6 yrs	3 yrs - 13 yrs
Maximum User Weight	77Lbs	77Lbs	132Lbs
Maximum User Height (Approx.)	46"	46"	61"
Frame Width	22.2"	22.2"	25.6"
Frame Length	28"	28"	39"
Angle Range	10 - 90 degrees	10 - 90 degrees	0 - 80 degrees
Minimum Height (Footplate to top of upper support) Prone/Supine	19.5"/22.5"	20.5"/23.5"	27.5"/30.5"
Maximum Height (Footplate to top of upper support)	34"	34"	46"
Minimum Height (Footplate to top of head support - Supine)	26"	26"	33"
Maximum Height (Footplate to top of head support - Supine)	46"	46"	53" - 59" *
Chest Pad Width (Min/Max)	6" - 10.25"	6" - 10.25"	6.75" - 12.5"
Hip Pad Width (Min/Max)	4.75" - 10.25"	4.75" - 10.25"	6.75" - 12.5"
Footplate Angle Range	+/- 15 degrees	+/- 15 degrees	+/- 15 degrees
Knee block width (center to center) min/max	4.75" - 8"	-	5.75" - 11"
Individual Leg Abduction - 5 degree increments	-	0-30 degrees	-

\* when Multigrip™ Body Support used

**Suggested HCPCS code: E0641\*\***

\*\*It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.













### CHOOSE YOUR COLOR (CONFIGURATIONS ARE FOR REPRESENTATION ONLY)



BLUE	01			
PINK	02			
GREEN	03			
RED	04			
BLACK	06			

Please note: all codes containing \*\* will color match - other items are black

### REQUIRED COMPONENTS - Please select one option from each section

	Size 1 Standard	Size 1 ABD	Size 2
<b>Base</b>	 <p>White metal base with four locking casters and gas strut controlled angle adjustment including safety lock.</p> <p><b>MUS01</b>    <b>\$1,333</b></p>	 <p>White metal base with four locking casters and gas strut controlled angle adjustment including safety lock.</p> <p><b>MUS01</b>    <b>\$1,333</b></p>	 <p>White metal base with four locking casters and linear winding mechanism controlled angle adjustment.</p> <p><b>MUS02</b>    <b>\$1904</b></p>
<b>Column</b>	 <p>Support column for height adjustment and mounting support boards. Comes with one piece knee support receiver bracket. <i>Size 1 Base Only</i></p> <p><b>C1</b>    <b>\$662</b></p>	 <p>Support column for height adjustment and mounting support boards. Comes with individual abd legs and knee support receiver brackets. <i>Size 1 Base Only</i></p> <p><b>C1ABD</b>    <b>\$1,421</b></p>	 <p>Support column for height adjustment and mounting support boards. Comes with one piece knee support receiver bracket. <i>Size 2 Base Only</i></p> <p><b>C2</b>    <b>\$945</b></p>
<b>Footplate</b>	 <p>One piece angle adjustable footplate offering plantar and dorsa flexion with mounting slots for sandals. <i>Size 1 Base Only</i></p> <p><b>FP1</b>    <b>\$484</b></p>	 <p>Individual angle adjustable footplates offering plantar and dorsa flexion with mounting slots for sandals. PAIR <i>Size 1 Base Only</i></p> <p><b>FPABD</b>    <b>\$726</b></p>	 <p>One piece angle adjustable footplate offering plantar and dorsa flexion with mounting slots for sandals. <i>Size 2 Base Only</i></p> <p><b>FP2</b>    <b>\$692</b></p>
<b>Support Boards</b>	 <p>PU chest and hip boards in color of choice with mounting slots for trunk and pelvic support. PAIR (W) 17.5" (H) 6" (D) 1.5". <i>Size 1 Base Only</i></p> <p><b>1B-**</b>    <b>\$927</b></p>	 <p>PU chest and hip boards in color of choice with mounting slots for trunk and pelvic support. PAIR (W) 17.5" (H) 6" (D) 1.5". <i>Size 1 Base Only</i></p> <p><b>1B-**</b>    <b>\$927</b></p>	 <p>PU chest and hip boards in color of choice with mounting slots for trunk and pelvic support. PAIR (W) 20.25" (H) 7" (D) 1.5". <i>Size 2 Base Only</i></p> <p><b>2B-**</b>    <b>\$1,326</b></p>

**REQUIRED COMPONENTS - Please select one option from each section**

**Trunk Support**

	<b>Small support pads w/de-rotational strap.</b> Wedge shaped lateral pads match chosen color. (D) 3.7" (W) 3.5" (H) 3.9" <b>PAIR</b>
<b>MB06-**</b>	<b>\$410</b> <input type="checkbox"/>
	<b>Size 1 flexi support pads w/de-rotational strap.</b> Wrap around style pads with padded fabric cover. (D) 7.1" (W) 0.9" (H) 3" <b>PAIR</b>
<b>MB33-**</b>	<b>\$547</b> <input type="checkbox"/>

	<b>Medium support pads w/de-rotational strap.</b> Wedge shaped lateral pads match chosen color. (D) 5.1" (W) 3.5" (H) 3.9" <b>PAIR</b>
<b>MB07-**</b>	<b>\$455</b> <input type="checkbox"/>
	<b>Size 1 LONG flexi support pads w/de-rotational strap.</b> Wrap around style pads with padded fabric cover. (D) 8.5" (W) 0.9" (H) 3" <b>PAIR</b>
<b>MB34-**</b>	<b>\$582</b> <input type="checkbox"/>

	<b>High-rise support pads w/de-rotational strap.</b> Double height wedge shaped lateral pad. (D)3.7" (W) 3.5" (H) 7.9" <b>PAIR</b>
<b>MB02-**</b>	<b>\$649</b> <input type="checkbox"/>
	<b>Size 2 flexi support pads with de-rotational strap.</b> Wrap around style pads with padded fabric cover. (D) 8.5" (W) 0.9" (H) 4.3" <b>PAIR</b>
<b>MB35-**</b>	<b>\$641</b> <input type="checkbox"/>

**Pelvic Support**

	<b>Small support pads w/de-rotational strap.</b> Wedge shaped lateral pads match chosen color. (D) 3.7" (W) 3.5" (H) 3.9" <b>PAIR</b>
<b>MB06-**</b>	<b>\$410</b> <input type="checkbox"/>

	<b>Medium support pads w/de-rotational strap.</b> Wedge shaped lateral pads match chosen color. (D) 5.1" (W) 3.5" (H) 3.9" <b>PAIR</b>
<b>MB07-**</b>	<b>\$455</b> <input type="checkbox"/>

**Knee Support**

	<b>Size 1 knee support w/padded strap.</b> Molded PU, . (D) 2.2" (W) 3.1" (H) 3.7" <b>PAIR</b>
<b>MB09-**</b>	<b>\$364</b> <input type="checkbox"/>

	<b>Size 2 knee support w/padded strap.</b> Molded PU, color matches that chosen. (D) 3.9" (W) 3.7" (H) 4.5" <b>PAIR</b>
<b>MB04-**</b>	<b>\$422</b> <input type="checkbox"/>

**Foot Positioning**














	<b>Size 1 Sandals with toe and heel straps.</b> (L) 6.8" (W) 3.9" (D) 0.7" <b>PAIR</b>
<b>S16</b>	<b>\$155</b> <input type="checkbox"/>

	<b>Size 2 Sandals with toe and heel strap</b> (L) 8.6" (W) 4.3" (D) 0.7" <b>PAIR</b>
<b>S2</b>	<b>\$177</b> <input type="checkbox"/>

	<b>Size 3 Sandals with toe and heel strap</b> (L) 10.8" (W)4.1" (D).7" <b>PAIR</b> <b>Size 2 Base Only</b>
<b>S3</b>	<b>\$191</b> <input type="checkbox"/>

## Retail Price List and Order Form 2022

### Optional Components

Optional Components					
<b>Tray</b>		<b>Size 1 tray.</b> For use in Prone and Supine. Height, Depth and Angle adjustable. <b>Size 1 Base only</b>			
	<b>MT02</b>	<b>\$701</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Tray Elbow Blocks.</b> For use with size 1 or 2 tray. Prevent retraction of the arms on the tray support surface. Gray PU. <b>PAIR</b>			
	<b>UBEB</b>	<b>\$250</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Size 2 tray.</b> For use in Prone and Supine. Height, Depth and Angle adjustable. <b>Size 2 Base Only</b>			
	<b>MT03</b>	<b>\$965</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Head Support</b>		<b>Molded head support.</b> Comes with adjustable bracket and PU molded pad with color match to main unit			
	<b>MH07-**</b>	<b>\$514</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>1 Link Multigrip™ head support.</b> Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.			
	<b>MH06</b>	<b>\$619</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>2 Link Multigrip™ head support.</b> Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.			
	<b>MH01 01</b>	<b>\$673</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>3 Link Multigrip™ head support.</b> Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.			
	<b>MH05</b>	<b>\$696</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Thoracic Support</b>		<b>Size 2.</b> With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability.			
	<b>HB02CB</b>	<b>\$292</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Size 3.</b> With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability. <b>Size 2 base only</b>			
	<b>HB03CB</b>	<b>\$309</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Size 4.</b> With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability. <b>Size 2 base only</b>			
	<b>HB04CB</b>	<b>\$329</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Prone chest prompt.</b> Comes with fabric cover that matches main color choice. Provides sternal pressure in prone. <b>Prone use only</b>			
	<b>MB03</b>	<b>\$387</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Multigrip™ Body Support.</b> Mount to either/both the support boards to offer additional support. Adjustable w/fabric cover. <b>Size 2 Base Only</b>			
	<b>MUS2SS-XX</b>	<b>\$455</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<b>Sandal Risers</b>		<b>Size 1.</b> Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser			
	<b>SRB1</b>	<b>1 set \$116</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>2 sets \$232</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Size 2.</b> Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser			
	<b>SRB2</b>	<b>1 set \$138</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>2 sets \$276</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>Size 3.</b> Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser			
	<b>SRB3</b>	<b>1 set \$151</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
		<b>2 sets \$302</b>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>