



Jenx Multistander - Supine



Retail Price List and Order Form 2020

Billing Information:		Ship To Information:	
Date:		Name:	
PO Number:		Address:	
Branch:			
Mark for:		City:	
Dealer Contact:		State:	
Email Address:		Zip Code:	
Phone:		Phone:	
Fax:			
		Quote	<input type="radio"/>
		Order	<input type="radio"/>

Size 1




Standard column and one piece footplate



Abduction column and individual footplates

Size 2



Standard column only

Configurations for representation only

	Size 1 standard	Size 1 ABD	Size 2
Approximate Age Range	9 mths - 6 yrs	9 mths - 6 yrs	3 yrs - 13 yrs
Maximum User Weight	77Lbs	77Lbs	132Lbs
Maximum User Height (Approx.)	46"	46"	61"
Frame Width	22.2"	22.2"	25.6"
Frame Length	28"	28"	39"
Angle Range	10 - 90 degrees	10 - 90 degrees	0 - 80 degrees
Minimum Height (Footplate to top of upper support) Prone/Supine	19.5"/22.5"	20.5"/23.5"	27.5"/30.5"
Maximum Height (Footplate to top of upper support)	34"	34"	46"
Minimum Height (Footplate to top of head support - Supine)	26"	26"	33"
Maximum Height (Footplate to top of head support - Supine)	46"	46"	53" - 59" *
Chest Pad Width (Min/Max)	6" - 10.25"	6" - 10.25"	6.75" - 12.5"
Hip Pad Width (Min/Max)	4.75" - 10.25"	4.75" - 10.25"	6.75" - 12.5"
Footplate Angle Range	+/- 15 degrees	+/- 15 degrees	+/- 15 degrees
Knee block width (center to center) min/max	4.75" - 8"	-	5.75" - 11"
Individual Leg Abduction - 5 degree increments	-	0-30 degrees	-

* when Multigrip™ Body Support used

Suggested HCPCS code: E0638**

**It is the responsibility of the CRT provider to determine the HCPCS code most appropriate for each circumstance. The HCPCS codes for reference only and are not to be deemed as legal advice. Actual coding and coverage varies by funding entities and the use of HCPCS does not guarantee coverage or payment.

Retail Price List and Order Form 2020













CHOOSE YOUR COLOR (CONFIGURATIONS ARE FOR REPRESENTATION ONLY)



BLUE	01		
PINK	02		
GREEN	03		
RED	04		
BLACK	06		

Please note: all codes containing ** will color match - other items are black

REQUIRED COMPONENTS - Please select one option from each section

	Size 1 Standard	Size 1 ABD	Size 2
Base	 <p>White metal base with four locking casters and gas strut controlled angle adjustment including safety lock.</p>	 <p>White metal base with four locking casters and gas strut controlled angle adjustment including safety lock.</p>	 <p>White metal base with four locking casters and linear winding mechanism controlled angle adjustment.</p>
	<p>MUS01-SUP \$1,134</p>	<p>MUS01-SUP \$1,134</p>	<p>MUS02 \$1,629</p>
Column	 <p>Support column for height adjustment and mounting support boards. Comes with one piece knee support receiver bracket. <i>Size 1 Base Only</i></p>	 <p>Support column for height adjustment and mounting support boards. Comes with individual abd legs and knee support receiver brackets. <i>Size 1 Base Only</i></p>	 <p>Support column for height adjustment and mounting support boards. Comes with one piece knee support receiver bracket. <i>Size 2 Base Only</i></p>
	<p>C1 \$573</p>	<p>C1ABD \$1,230</p>	<p>C2 \$818</p>
Footplate	 <p>One piece angle adjustable footplate offering plantar and dorsa flexion with mounting slots for sandals. <i>Size 1 Base Only</i></p>	 <p>Individual angle adjustable footplates offering plantar and dorsa flexion with mounting slots for sandals. <i>Size 1 Base Only</i></p>	 <p>One piece angle adjustable footplate offering plantar and dorsa flexion with mounting slots for sandals. <i>Size 2 Base Only</i></p>
	<p>FP1 \$419</p>	<p>FPABD \$629</p>	<p>FP2 \$599</p>
Support Boards	 <p>Pair in color of choice with mounting slots for trunk and pelvic support. Made from wipe clean PU. (W) 17.5" (H) 6" (D) 1.5". <i>Size 1 Base Only</i></p>	 <p>Pair in color of choice with mounting slots for trunk and pelvic support. Made from wipe clean PU. (W) 17.5" (H) 6" (D) 1.5". <i>Size 1 Base Only</i></p>	 <p>Pair in color of choice with mounting slots for trunk and pelvic support. Made from wipe clean PU. (W) 20.25" (H) 7" (D) 1.5". <i>Size 2 Base Only</i></p>
	<p>1B-** \$803</p>	<p>1B-** \$803</p>	<p>2B-** \$1,148</p>

Retail Price List and Order Form 2020

REQUIRED COMPONENTS - Please select one option from each section

Trunk Support

	Small support pads w/de-rotational strap. Wedge shaped lateral pads match chosen color. (D) 3.7" (W) 3.5" (H) 3.9" PAIR
MB06-**	\$355 <input type="checkbox"/>
	Size 1 flexi support pads w/de-rotational strap. Wrap around style pads with padded fabric cover. (D) 7.1" (W) 0.9" (H) 3" PAIR
MB33-**	\$473 <input type="checkbox"/>

	Medium support pads w/de-rotational strap. Wedge shaped lateral pads match chosen color. (D) 5.1" (W) 3.5" (H) 3.9" PAIR
MB07-**	\$394 <input type="checkbox"/>
	Size 1 LONG flexi support pads w/de-rotational strap. Wrap around style pads with padded fabric cover. (D) 8.5" (W) 0.9" (H) 3" PAIR
MB34-**	\$504 <input type="checkbox"/>

	High-rise support pads w/de-rotational strap. Double height wedge shaped lateral pad. (D)3.7" (W) 3.5" (H) 7.9" PAIR
MB02-**	\$562 <input type="checkbox"/>
	Size 2 flexi support pads with de-rotational strap. Wrap around style pads with padded fabric cover. (D) 8.5" (W) 0.9" (H) 4.3" PAIR
MB35-**	\$555 <input type="checkbox"/>

Pelvic Support

	Small support pads w/de-rotational strap. Wedge shaped lateral pads match chosen color. (D) 3.7" (W) 3.5" (H) 3.9" PAIR
MB06-**	\$355 <input type="checkbox"/>

	Medium support pads w/de-rotational strap. Wedge shaped lateral pads match chosen color. (D) 5.1" (W) 3.5" (H) 3.9" PAIR
MB07-**	\$394 <input type="checkbox"/>

Knee Support

	Size 1 knee support w/padded strap. Molded PU, . (D) 2.2" (W) 3.1" (H) 3.7" PAIR
MB09-**	\$315 <input type="checkbox"/>

	Size 2 knee support w/padded strap. Molded PU, color matches that chosen. (D) 3.9" (W) 3.7" (H) 4.5" PAIR
MB04-**	\$366 <input type="checkbox"/>

Foot Positioning





	Size 1 Sandals with toe and heel straps. (L) 6.8" (W) 3.9" (D) 0.7" PAIR
S16	\$134 <input type="checkbox"/>

	Size 2 Sandals with toe and heel strap (L) 8.6" (W) 4.3" (D) 0.7" PAIR
S2	\$153 <input type="checkbox"/>

	Size 3 Sandals with toe and heel strap (L) 10.8" (W)4.1" (D).7" PAIR Size 2 Base Only
S3	\$166 <input type="checkbox"/>

Retail Price List and Order Form 2020

Optional Components

Optional Components				
Tray		Size 1 tray. For use in Prone and Supine. Height, Depth and Angle adjustable. <i>Size 1 Base only</i>		
	MT02	\$607	<input type="checkbox"/>	
		Tray Elbow Blocks. For use with size 1 or 2 tray. Prevent retraction of the arms on the tray support surface. Gray PU. PAIR		
	UBEB	\$216	<input type="checkbox"/>	
		Size 2 tray. For use in Prone and Supine. Height, Depth and Angle adjustable. <i>Size 2 Base Only</i>		
	MT03	\$835	<input type="checkbox"/>	
Head Support		Molded head support. Comes with adjustable bracket and PU molded pad with color match to main unit		
	MH07-**	\$445	<input type="checkbox"/>	
		1 Link Multigrip™ head support. Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.		
	MH06	\$536	<input type="checkbox"/>	
		2 Link Multigrip™ head support. Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.		
	MH01 01	\$583	<input type="checkbox"/>	
		3 Link Multigrip™ head support. Locking link contour system. Comes with padded fabric cover and adjustable mounting bracket.		
	MH05	\$603	<input type="checkbox"/>	
Thoracic Support		Size 2. With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability.		
	HB02CB	\$252	<input type="checkbox"/>	
		Size 3. With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability. <i>Size 2 base only</i>		
	HB03CB	\$268	<input type="checkbox"/>	
		Size 4. With zipper to allow for access to G tube feed site. Ideal for low tone children with poor trunk stability. <i>Size 2 base only</i>		
	HB04CB	\$285	<input type="checkbox"/>	
		Multigrip™ Body Support. Mount to either/both the support boards to offer additional support. Adjustable w/fabric cover. <i>Size 2 Base Only</i>		
	MUS2SS-XX	\$394	<input type="checkbox"/>	
Sandal Risers		Size 1. Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser		
	SRB1	1 set	\$99	<input type="checkbox"/>
		2 sets	\$198	<input type="checkbox"/>
		Size 2. Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser		
	SRB2	1 set	\$119	<input type="checkbox"/>
		2 sets	\$218	<input type="checkbox"/>
	Size 3. Set of 4 individual raising blocks with mounting hardware. Intended to reduce the height under one sandal. (H) 0.7" per raiser			
SRB3	1 set	\$130	<input type="checkbox"/>	
	2 sets	\$260	<input type="checkbox"/>	