



#### C-HLA-2, C-HLA-2T, C-CB-L2, C-HLA

- Innovative 6-point cradle works with an array of Hoyer slings, including 2, 4, and 6-point slings
- Extended reach for floor pick-up capabilities
- Adjustable legs (does not apply to C-CB-L2)
- 400 lb. resident capacity
- Long-wearing chrome finish maintains appearance
- Mast and base separate for easy travel or storage
- Standard 4" casters are offset for greater stability and smooth operation

## Hoyer® Classics Product Overview

Hoyer Classic Lifts	
<b>HPL700/HPL700WSC</b> 700 lb. patient capacity	Power operated base, Base clearance of 4.5", 6-point bariatric cradle design minimizes patient sway
<b>HPL402</b> 400 lb. patient capacity	Adjustable U-base, Reverse-mounted actuator for additional knee clearance, 6-point cradle works with an array of Hoyer slings
<b>P-C-HLA-2</b> 400 lb. patient capacity	Adjustable U-base with 5.5" clearance, 6-point cradle works with an array of Hoyer slings
<b>C-HLA-2, C-HLA-2T, C-HLA, C-CB-L2</b> 400 lb. patient capacity	Adjustable legs, Cradle design minimizes sway, 6-point non-sway cradle offers extra stability
<b>HML400</b> 400 lb. patient capacity	Adjustable U-base

## Hoyer® Classics Sling Overview



One Piece Commode Sling with Positioning Strap

- Durable, polyester with closed-cell foam padding for comfort
- Aperture for toileting without removing sling
- Weight Capacity: 600 lbs.



One Piece Sling with Positioning Strap

- Durable, polyester with closed-cell foam padding for comfort
- Provides easy position adjustment of patient
- Weight Capacity: 600 lbs.



U-Sling, Padded with Head Support

- Durable, polyester with closed-cell foam padding for comfort
- Nylon head support
- Weight Capacity: 600 lbs.



U-Sling, Padded

- Durable, polyester with closed-cell foam padding for comfort
- Weight Capacity: 600 lbs.

Hoyer Classic Slings	X-Large	Large	Medium	Small
<b>One Piece Commode Sling with Positioning Strap</b>	—	70056	70057	70058
<b>One Piece Sling with Positioning Strap</b>	—	70051	70052	70053
<b>U-Sling, Padded with Head Support</b>	70010	70011	70012	70013
<b>U-Sling, Padded</b>	70000	70001	70002	70003

HPL700 Slings	XX-Large	X-Large	Large	Medium
<b>Full Back</b>	NA1075	NA1074	NA1073	NA1072
<b>Quick Fit</b>	NA1012P	NA1011P	NA1010P	NA1009P



Full Back

- Weight Capacity: 850 lbs.



Quick Fit

- Weight Capacity: 850 lbs.