













2 Point Slings



Model	Item Number	Material	Chains	Straps	Head Support	Commode	HCPCS Code	Weight Capacity	C-HLA *	HML400	HPL402
9	9-D	Dacron	Yes	NA	No	Yes	E0621	400 lbs.	✓	✓	✓
	9-D-LC	Dacron	No	NA	No	Yes					
	9-P	Polydura	Yes	NA	No	Yes					
	9-P-LC	Polydura	No	NA	No	Yes					
110	110-D	Dacron	NA	Yes	Yes	No	E0621	400 lbs.	✓	✓	✓
	110-D-LS	Dacron	NA	No	Yes	No					
	110-P	Polydura	Yes	NA	Yes	No					
	110-P-LS	Polydura	No	NA	Yes	No					
	110-D-LS-510	Dacron	No	NA	Yes	No					
	110-P-LS-510	Polydura	No	NA	Yes	No					
111	111-D	Dacron	NA	Yes	Yes	No	E0621	400 lbs.	✓	✓	✓
	111-D-LS	Dacron	NA	No	Yes	No					
	111-P	Polydura	NA	Yes	Yes	No					
	111-P-LS	Polydura	NA	No	Yes	No					
	111-D-LS-510	Dacron	NA	No	No	Yes					
	111-P-LS-510	Polydura	NA	No	No	Yes					
112	112-D	Dacron	Yes	NA	No	No	E0621	400 lbs.	✓	✓	✓
	112-D-LC	Dacron	No	NA	No	No					
	112-P	Polydura	Yes	NA	No	No					
	112-P-LC	Polydura	No	NA	No	No					
113	113-D	Dacron	Yes	NA	No	Yes	E0621	400 lbs.	✓	✓	✓
	113-D-LC	Dacron	No	NA	No	Yes					
	113-P	Polydura	Yes	NA	No	Yes					
	113-P-LC	Polydura	No	NA	No	Yes					
114	114-D	Dacron	NA	Yes	No	No	E0621	400 lbs.	✓	✓	✓
	114-D-LS	Dacron	NA	No	No	No					
	114-P	Polydura	NA	Yes	No	No					
	114-P-LS	Polydura	NA	No	No	No					
115	115-D	Dacron	NA	Yes	No	Yes	E0621	400 lbs.	✓	✓	✓
	115-D-LS	Dacron	NA	No	No	Yes					
	115-P	Polydura	NA	Yes	No	Yes					
	115-P-LS	Polydura	NA	No	No	Yes					
117	117-D	Dacron	Yes	NA	Yes	Yes	E0621	400 lbs.	✓	✓	✓
	117-D-LC	Dacron	No	NA	Yes	Yes					
	117-P	Polydura	Yes	NA	Yes	Yes					
	117-P-DC	Polydura	No	NA	Yes	Yes					
118	118-D	Dacron	Yes	No	Yes	No	E0621	400 lbs.	✓	✓	✓
	118-D-LC	Dacron	No	No	Yes	No					
	118-P	Polydura	Yes	NA	Yes	No					
	118-P-LC	Polydura	No	NA	Yes	No					
127	127-P	Polydura	NA	Yes	Yes	Yes-Split	E0621	400 lbs.	✓	✓	✓
	127-P-LS	Polydura	NA	No	Yes	Yes-Split					
	127-P-LS-510	Polydura	NA	No	No	Yes-Split					
NA	020-0236	Nylon Mesh	NA	NA	Yes	No	E0621	400 lbs.	✓	✓	✓

*Includes: C-HLA-2, C-HLA-2T, P-C-HLA-2

4 Point Slings



Model	Item Number	Size	Material	Head Support	Commode	HCPCS Code	Weight Capacity	C-HLA-2	C-HLA-2T	HML400	HPL402	HPL600
Universal Padded Sling	70000	XL	Polyester	None	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	70001	L										
	70002	M										
	70003	S										
Universal Padded Sling	70010	XL	Polyester	Yes	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	70011	L										
	70012	M										
	70013	S										
Universal Mesh Sling	50020	XL	Nylon Mesh	None	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	50021	L										
	50022	M										
	50023	S										
Universal Mesh Sling	50025	L	Nylon Mesh	Yes	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	50026	M										
	50027	S										
One Piece Sling with Strap	70051	L	Polyester	No	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	70052	M										
	70053	S										
One Piece Sling with Strap and Commode	70056	L	Polyester	No	No	E0621	600 lbs.	✓	✓	✓	✓	✓
	70057	M										
	70058	S										

Bariatric Slings

Model	Item Number	Weight Capacity	Cradle	HPL700
Full Back Padded, Small (Red)	NA1072	850 lbs.	Loop	✓
Full Back Padded, Medium (Yellow)	NA1073			
Full Back Padded, Large (Green)	NA1074			
Full Back Padded, X Large (Blue)	NA1075			
Quick Fit Padded, Small (Red)	NA1009P	850 lbs	Loop	✓
Quick Fit Padded, Medium (Yellow)	NA1010P			
Quick Fit Padded, Large (Green)	NA1011P			
Quick Fit Padded, X Large (Blue)	NA1012P			

