

Ergo Flight TP (S-2512-TP)

Ultra Lightweight Ergonomic Transporter



- Removable Swing Away Footrest
- Patented S-Style Ergonomic Seat Frame
- 6061 T-6 Aircraft-grade Aluminum Only 18 lbs. (w/o footrests)
- 3-Stage handle brake: allows light to firm grip for lock
- Built in Silver Aegis Anti-bacterial Cushion
- Folding Backrest/Folding Seat for Transporting in Vehicle or Travel
- Fixed Armrests w/Concaved Armpads
- Pocket Behind Backrest & Small Carry Pouch on Each Armrest
- 14" Rear Polyurethane, High Tread, Flat Free Wheels
- 6" x 1" Polyurethane Front Casters
- 16" x 17" or 18" x 17" Seat Width
- Frame Color: Pearl Silver
- Weight Capacity: 220 lbs

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
S-2512F16S-TP	16	17	8	18	17	36	23	12	36	18	220
S-2512F18S-TP	18	17	8	18	17	36	25	12	36	18	220

**Available Options ~ Anti-Tipper, Seat belt, Memory Foam seat and back cushions, Rigidifying Push-Bar, Foldable/Adjustable Headrest*

karma KARMAN HEALTHCARE Inc.