

## ATX

## Ergonomic Class

# ULTRA LIGHTWEIGHT ACTIVE WHEELCHAIR

SPECIFICATIONS	ATX-14	ATX-16	ATX-18
Seat Width	14 in.	16 in.	18 in.
Seat Depth	15 in.	16 in.	18 in.
Armrest Height	23 in.	23 in.	23 in.
Seat Height	14 in.	14 in.	14 in.
Back Height	Varies	Varies	Varies
Overall Height	Varies	Varies	Varies
Overall Open Width	Varies	Varies	Varies
Folded Width	14 in.	14 in.	14 in.
Overall Length	Varies	Varies	Varies
Weight w/o Riggings	Varies	15.4	Varies
Weight Capacity	275 lbs.	275 lbs.	275 lbs.
HCPCS Code	N/A	N/A	N/A

(SUBJECT TO CHANGE WITHOUT NOTICE)



### THE ULTIMATE STATEMENT IS ONE WITH NO COMPROMISE...

- Foot plate height unlimited increment adj.
- 15.4 LB Starting Frame Weight
- Back Cane Height and Angle Adjustable
- One piece foldable footplate
- Fixed Footrests w/ Extra Wide Footplates
- Center of Gravity Adjustable
- S-Shape Frame (patented)
- Built in Silver Aegis Anti-bacterial Cushion
- Metallic Aegis Treated High Gloss Paint
- Foot plate height unlimited increment adj.