



S-Ergo 105

Lightweight Ergonomic Manual Wheelchair



- High strength, weighs only 27 lbs. (w footrests)
- S-Style Ergonomic Seat
- Silver ¼” removable Anti-Bacterial cushion
- Folding backrest/ folding seat for transporting vehicle or travel.
- Fixed armrest w/ wider concaved armpads
- Fixed Footrest
- “Tube-in Center” foot-plate, assures better side leg support
- Extending push to lock manual brakes
- Pouch for carrying small items attached to upholstery behind backrest.
- 24” Rear polyurethane, high profile, flat free wheels
- 7”x1” Polyurethane front caster
- 16x17 or 18x17 Seat width
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Black Breathable Mesh Bottom
- Chair comes fully assembled except for footrests, tools included
- Weight Capacity of 250 lbs

Unique Features for S-Ergonomic Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slippage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider arm pad, providing comfort/support

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
S-105F16	16	17	8	19	17	35	23	12	39	27	250
S-105F18	18	17	8	19	17	35	25	12	39	27	250

**Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Foldable Canopy, Rigidifying Push-Bar, Reflecting Spoke Guard, Foldable/Adjustable Headrest, Frog Leg Suspension*



KARMAN HEALTHCARE Inc.