

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

Select if chair quote has:  Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** \_\_\_\_\_

Effective 08/19/2020

**Liberty FT Basic** E1161 **\$2795**  
 250 lb weight capacity. Tilt range: 0° - 20°. Standard seat upholstery, Push to Lock wheel locks.

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

## Frame

**Seat Width**

*14"-16" Seat Widths ship with 300 N Assistive Gas Strut, 18"-22" Seat Widths ship with 450 N Assistive Gas Strut*

14" \_\_\_\_\_ N/C

16" \_\_\_\_\_ N/C

18" \_\_\_\_\_ N/C

20" E2201 \_\_\_\_\_ \$470

22" E2201 \_\_\_\_\_ \$470

**Seat Depth**

14" \_\_\_\_\_ N/C

16" \_\_\_\_\_ N/C

18" \_\_\_\_\_ N/C

20" E2203 \_\_\_\_\_ \$605

**Transit**

Transit System \_\_\_\_\_ \$320

Omit Transit \_\_\_\_\_ N/C

## Backrest

Back Post / Height Range / Back Height							
Back Post	Price	Height Range	Back Height				
<input type="checkbox"/> Height Adjustable 8° Bend with Push Handle	\$25	<input type="checkbox"/> Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	<input type="checkbox"/> 21	<input type="checkbox"/> 22

**Backrest Options**

Fixed (Non-Folding) \_\_\_\_\_ Std

Angle Adjustable (Folding) \_\_\_\_\_ \$95

**Back Upholstery**

*If selecting both upholstery and Axiom Back, the upholstery will be installed*

Standard \_\_\_\_\_ Std

Omit Back Upholstery \_\_\_\_\_ - \$15

## Footrest

Hanger / Front Frame Bend				
Hanger	60°	70°		Retail
Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>		Std
Pro Elevating Leg Rests (PRO ELR) E0990, Rental K0195			<input type="checkbox"/>	\$295
Omit Hangers			<input type="checkbox"/>	N/C

**Footplate**

Composite \_\_\_\_\_ Std

Aluminum Angle Adjustable k0040 \_\_\_\_\_ \$190

**Footrest continued**

**Footrest / Extension / Seat to Footrest Length**

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Seat Height recommended. Extension Tube Hanger measurements based on 70° Hanger Bend.

Footrest	Extension	Seat to Footrest Length								
Extension Tube Hangers (60°, 70°)	<input type="checkbox"/> Short	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	
	<input type="checkbox"/> Medium	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5
	<input type="checkbox"/> Long	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5								
	<input type="checkbox"/> Long	13 - 19								
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	8.5 - 17								
	<input type="checkbox"/> Long	11.5 - 17.5								

**Heel Loop**

- Heel Loops E0951 \_\_\_\_\_ \$30
- Omit Heel Loops \_\_\_\_\_ N/C

**Calf Strap**

- Velcro Adjustable K0038 \_\_\_\_\_ \$35

**Hanger Release**

- Option Required*
- 4-Way - Out \_\_\_\_\_ N/C

**Seat Height, Casters and Rear Wheels**

**Seat Height / Caster Size / Rear Wheel Size**

Select Seat Height corresponding to desired Rear Wheel and Caster Size

Rear Wheel Size	Caster Size		
	4"	5"	6"
16"	Seat Height		
	<input type="checkbox"/> 13.5 <sup>1</sup>	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 14 <sup>1</sup>	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 17
	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
20"	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 16	<input type="checkbox"/> 17
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
		<input type="checkbox"/> 18.5	
	<input type="checkbox"/> 19		

Rear Wheel Size	Caster Size		
	4"	5"	6"
22"	Seat Height		
	<input type="checkbox"/> 14	<input type="checkbox"/> 14	<input type="checkbox"/> 15
	<input type="checkbox"/> 14.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 15	<input type="checkbox"/> 15	<input type="checkbox"/> 16
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 16	<input type="checkbox"/> 17
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17	<input type="checkbox"/> 18
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18	<input type="checkbox"/> 19
24"	<input type="checkbox"/> 15	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15
	<input type="checkbox"/> 15.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5
	<input type="checkbox"/> 16	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16
		<input type="checkbox"/> 16	<input type="checkbox"/> 16.5
		<input type="checkbox"/> 16.5	<input type="checkbox"/> 17
		<input type="checkbox"/> 17	<input type="checkbox"/> 17.5
		<input type="checkbox"/> 17.5	<input type="checkbox"/> 18
		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5
		<input type="checkbox"/> 18.5	<input type="checkbox"/> 19
		<input type="checkbox"/> 19	

<sup>1</sup> Available only with Low Profile tires

**Seat Height, Casters and Rear Wheels** continued

**Caster**

		Caster Size			
		Price	4"	5"	6"
Caster Type	1" Poly	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	1.5" Poly	\$50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**Rear Wheel**

		Rear Wheel Size				
		Price	16"	20"	22"	24"
Rear Wheel Type	Mag	Std	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C				<input type="checkbox"/>
	MAXX Performance Spoke	\$180		<input type="checkbox"/>	<input type="checkbox"/>	

**Wheel Locks**

**Wheel Lock**

*Option Required*

Push to Lock \_\_\_\_\_ N/C

**Wheel Lock Options**

Extension Handles E0961 \_\_\_\_\_ \$67

**Accessories**

**Anti-Tip**

Rear Anti-Tips E0971 \_\_\_\_\_ \$150

Omit Rear Anti-Tips \_\_\_\_\_ N/C

**Other Accessories**

1.5" Auto Buckle E0978 \_\_\_\_\_ \$70

Seat Belt Mounting Hardware \_\_\_\_\_ \$30

O<sub>2</sub> Holder E2208 \_\_\_\_\_ \$295

*Mounted to the left side. Not available with Fold Down Back.*

Cane and Crutch Holder \_\_\_\_\_ \$195

**Rear Tire Type**

**Full Profile**

Street Pneumatic E2211, E2212 \_\_\_\_\_ Std

Full Poly E2218 \_\_\_\_\_ N/C

Pneumatic with Airless Insert E2211, E2213 \_\_\_\_\_ \$80

**Low Profile**

Low Poly E2218 \_\_\_\_\_ N/C

*Available only with 16" Mag Wheels*

**Handrim**

*Not Available with 16" Rear Wheels*

Aluminum Anodized \_\_\_\_\_ Std

Plastic Coated \_\_\_\_\_ \$115

**Axle**

Quick Release \_\_\_\_\_ Std

**Armrests**

**Armrest**

*Selection Required*

Height Adjustable T-Arm (8.5" - 13.5") E0973

Desk \_\_\_\_\_ \$265

Full \_\_\_\_\_ \$265

**Arm Pad**

Standard Arm Pad \_\_\_\_\_ N/C

**Frame Color**

**Gloss Colors**

Black \_\_\_\_\_ N/C

Candy Blue \_\_\_\_\_ N/C

Candy Red \_\_\_\_\_ N/C

**Matte Colors**

Charcoal Gray \_\_\_\_\_ N/C

Shadow Black \_\_\_\_\_ N/C