

LECKEY°

Distributed by **Ki** Mobility

Squiggles Stander



Promoting independence & function.

The **Squiggles Stander** is a versatile, bright and colorful 3-in-1 stander for children ages 1-5. This allows the clinician to maximize postural requirements through pelvic stability, trunk and head alignment and leg and foot positioning support for standing therapy.

Versatility




3-in-1

The **Squiggles Stander** is an extremely versatile 3-in-1 stander, offering prone, upright and supine standing in one product.



Supine

Supports the back of the body and offers lateral and anterior supports to maintain alignment. Allows for partial weight bearing and promotes control of head, trunk, pelvis, knees and feet.



Upright

Designed for users with good balance and trunk control but cannot stand independently.



Prone

Provides anterior support and stability while encouraging weight bearing and use of upper extremities.

Posture

Head & Trunk Alignment

Whether prone, supine or upright, the chest and sternum pads, lateral supports and range of head supports provide trunk and head stability.



Pelvic Stability

Height, depth, width and angle adjustable pelvic support ensures the pelvis is held in a secure and comfortable position.



Leg & Foot Positioning

As legs and feet are taking the weight of the child, proper support is essential. Fully adjustable knee cups and footplates ensure weight is distributed evenly and in the right places.



Specifications

Age	1-5 Years
User Height	29.5" – 43.3"
Weight Capacity	48.4 lbs.
Available Colors	Green, Orange, Blue, Pink

LECKEY[®]