


LECKEY[®]

Advance Bath Chair



Advance Bath Chair



Accessories



LAB/#/B Abduction Belt



LAB/1/C Head support



Advance Bath Chair

The Advance Bath Chair is available in 4 sizes for infants through to young adults. The product provides the user with appropriate levels of postural support ensuring they are safe and secure at all times during bathing.

The simple and easy to use adjustments mean that the product can be set up in a variety of positions in seconds depending on the user's postural requirements.

The height adjustable legs mean that the product can be set up at the appropriate height for the parent or carer to reduce the risk of back strain or injury.

Bath Support

The Bath Support, available in 2 sizes for children aged 1-9 years of age, complements the Advance Bath Chair, providing support for the less involved child.

The bath support easily folds flat for storage & transportation.

Shower Trolley

With the addition of the adjustable Shower Trolley, the Advance Bath Chair can be converted in seconds to a versatile shower chair.

Manufactured from medical grade stainless steel, the robust frame and high quality lockable castors ensure the trolley can be easily and safely manoeuvred over wet and slippery floors.



Product Sizing

Size	1		2		3		4	
Age	1 - 5		4 - 9		8 - 14		12 - 18	
	mm	inches	mm	inches	mm	inches	mm	inches
User Height								
Min	750	30	950	37	1150	45	1350	53
Max	1050	41	1250	49	1450	57	1650	65
Advance Bath Chair	LAB/1		LAB/2		LAB/3		LAB/4	
Head support	150	6	150	6	215	8.5	215	8.5
Back support	350	14	400	16	460	18	550	22
Seat depth	225	9	275	11	335	13	400	16
Calf support	260	10	260	10	340	13	340	13
Seat Height	50-425mm / 2-17"							
Max Weight	72kg / 159lbs							
Total Width	430mm / 17"							
Shower Trolley	LSF/1							
Width	600	24						
Length:								
Min	900	36						
Max	1150	45						
Seat Height from floor	650-700mm / 26"- 28"							
Bath Support	LBS/1		LBS/2					
Head support	150	6	150	6				
Back support	600	23.6	700	27.5				
Seat depth	300	12	400	15.7				
Seat Height	0-150mm / 0-6"							
Max Weight	62kg / 137lbs							
Total Width	400mm / 16"							

24 hour postural
care for babies,
kids & adults.
Sleeping, Sitting,
Standing, Walking,
Moving, Bathing,
Toileting.