



Easy Seat

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote PO Number: _____ City: _____ State: _____
 RTS: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____

Effective 02/04/2019

Easy Seat

Easy Seat Complete

Selection required if Easy Seat Basic is not selected. Includes Hip Belt, 1" Side Pads, Tray, Footplate, Sandal Straps, 2 Swivel and 2 Fixed Casters.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Complete - Meadow	LES/1/COMP-MW	\$1381
<input type="checkbox"/>	Size 2 Easy Seat Complete - Mondeo Blue	LES/2/COMP-MB	\$1428
<input type="checkbox"/>	Size 3 Easy Seat Complete - Mondeo Blue	LES/3/COMP-MB	\$1549
<input type="checkbox"/>	Size 4 Easy Seat Complete - Mondeo Blue	LES/4/COMP-MB	\$1654

Easy Seat Basic

Selection required if Easy Seat Complete is not selected. Includes Non-Wheeled Legs, Hip Belt and 1" Side Pads.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Basic - Meadow	LES/1/FL-MW	\$898
<input type="checkbox"/>	Size 2 Easy Seat Basic - Mondeo Blue	LES/2/FL-MB	\$940
<input type="checkbox"/>	Size 3 Easy Seat Basic - Mondeo Blue	LES/3/FL-MB	\$1045
<input type="checkbox"/>	Size 4 Easy Seat Basic - Mondeo Blue	LES/4/FL-MB	\$1132

Head Support

Selection of Flat High Back Headrest is required with Lateral Supports. Lateral selection must match Headrest size.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Flat High Back Headrest - Meadow	LES/1/L-MW	\$269
<input type="checkbox"/>	Size 2 Flat High Back Headrest - Mondeo Blue	LES/2/L-MB	\$276
<input type="checkbox"/>	Size 3 Flat High Back Headrest - Mondeo Blue	LES/3/L-MB	\$289
<input type="checkbox"/>	Size 4 Flat High Back Headrest - Mondeo Blue	LES/4/L-MB	\$292
<input type="checkbox"/>	Size 1 Lateral Supports for Flat High Back Headrest - Meadow (Pair)	LES/1/M-MW	\$107
<input type="checkbox"/>	Size 2 Lateral Supports for Flat High Back Headrest - Mondeo Blue (Pair)	LES/2/M-MB	\$114
<input type="checkbox"/>	Size 3 Lateral Supports for Flat High Back Headrest - Mondeo Blue (Pair)	LES/3/M-MB	\$126
<input type="checkbox"/>	Size 4 Lateral Supports for Flat High Back Headrest - Mondeo Blue (Pair)	LES/4/M-MB	\$131

LECKEY[®]

Trunk Support

Select only Flexible Lateral Supports or Rigid Lateral Supports

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Flexible Lateral Supports - Meadow (Pair)	LES/1/D-MW	\$177
<input type="checkbox"/>	Size 1 Flexible Lateral Supports - Mondeo Blue (Pair)	LES/1/D-MB	\$177
<input type="checkbox"/>	Size 2 Flexible Lateral Supports - Mondeo Blue (Pair)	LES/2/D-MB	\$184
<input type="checkbox"/>	Size 3 Flexible Lateral Supports - Mondeo Blue (Pair)	LES/3/D-MB	\$189
<input type="checkbox"/>	Size 4 Flexible Lateral Supports - Mondeo Blue (Pair)	LES/4/D-MB	\$193
<input type="checkbox"/>	Size 1 and 2 Rigid Lateral Supports - Meadow (Pair)	LES/2/E-MW	\$216
<input type="checkbox"/>	Size 1 and 2 Rigid Lateral Supports - Mondeo Blue (Pair)	LES/2/E-MB	\$216
<input type="checkbox"/>	Size 3 and 4 Rigid Lateral Supports - Mondeo Blue (Pair)	LES/3/E-MB	\$240

Additional Side Pads

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 1/2" Thick - Meadow (Pair)	LES/1/B1-MW	\$93
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/2/B1-MB	\$103
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/3/B1-MB	\$126
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 1/2" Thick - Mondeo Blue (Pair)	LES/4/B1-MB	\$139
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 1" Thick - Meadow (Pair)	LES/1/B2-MW	\$113
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/2/B2-MB	\$126
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/3/B2-MB	\$137
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 1" Thick - Mondeo Blue (Pair)	LES/4/B2-MB	\$139
<input type="checkbox"/>	Size 1 Easy Seat Side Pads 2" Thick - Meadow (Pair)	LES/1/B3-MW	\$113
<input type="checkbox"/>	Size 2 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/2/B3-MB	\$123
<input type="checkbox"/>	Size 3 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/3/B3-MB	\$146
<input type="checkbox"/>	Size 4 Easy Seat Side Pads 2" Thick - Mondeo Blue (Pair)	LES/4/B3-MB	\$154

Pommel

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Pommel - Meadow	LES/1/C-MW	\$113
<input type="checkbox"/>	Size 1 Easy Seat Pommel - Mondeo Blue	LES/1/C-MB	\$113
<input type="checkbox"/>	Size 2 Easy Seat Pommel - Mondeo Blue	LES/2/C-MB	\$126
<input type="checkbox"/>	Size 3 and 4 Easy Seat Pommel - Mondeo Blue	LES/3/C-MB	\$135

Potty Trainer Seat Insert

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Potty Trainer Seat Insert	LES/1/S	\$76
<input type="checkbox"/>	Size 2 Easy Potty Trainer Seat Insert	LES/2/S	\$76
<input type="checkbox"/>	Size 3 Easy Potty Trainer Seat Insert	LES/3/S	\$108
<input type="checkbox"/>	Size 4 Easy Potty Trainer Seat Insert	LES/4/S	\$108

LECKEY[®]

Padded Seat Insert

Used with Poty Trainer Seat Insert. Must match Easy Seat size.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Padded Seat Insert - Meadow	LES/1/R-MW	\$83
<input type="checkbox"/>	Size 2 Padded Seat Insert - Mondeo Blue	LES/2/R-MB	\$95
<input type="checkbox"/>	Size 3 Padded Seat Insert - Mondeo Blue	LES/3/R-MB	\$103
<input type="checkbox"/>	Size 4 Padded Seat Insert - Mondeo Blue	LES/4/R-MB	\$114

Foot Positioning

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Footplate with Sandal Straps	LES/1/N	\$130
<input type="checkbox"/>	Size 2 Footplate with Sandal Straps	LES/2/N	\$137
<input type="checkbox"/>	Size 3 Footplate with Sandal Straps	LES/3/N	\$148
<input type="checkbox"/>	Size 4 Footplate with Sandal Straps	LES/4/N	\$159
<input type="checkbox"/>	Size 1 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/1/K	\$119
<input type="checkbox"/>	Size 2 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/2/K	\$130
<input type="checkbox"/>	Size 3 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/3/K	\$142
<input type="checkbox"/>	Size 4 Easy Seat Footplate Raiser 1"	Quantity: _____ (1-3) LES/4/K	\$154

Accessories

Must match Easy Seat size

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Easy Seat Activity Tray	LES/1/A	\$277
<input type="checkbox"/>	Size 2 Easy Seat Activity Tray	LES/2/A	\$277
<input type="checkbox"/>	Size 3 Easy Seat Activity Tray	LES/3/A	\$277
<input type="checkbox"/>	Size 4 Easy Seat Activity Tray	LES/4/A	\$280
<input type="checkbox"/>	Size 1 Easy Seat Grab/Suspension Rail	LES/1/I	\$256
<input type="checkbox"/>	Size 2 Easy Seat Grab/Suspension Rail	LES/2/I	\$256
<input type="checkbox"/>	Size 3 Easy Seat Grab/Suspension Rail	LES/3/I	\$256
<input type="checkbox"/>	Size 4 Easy Seat Grab/Suspension Rail	LES/4/I	\$256
<input type="checkbox"/>	Size 1 Easy Seat Grab Posts (Pair)	LES/1/P	\$256
<input type="checkbox"/>	Size 2 Easy Seat Grab Posts (Pair)	LES/2/P	\$256
<input type="checkbox"/>	Size 3 Easy Seat Grab Posts (Pair)	LES/3/P	\$256
<input type="checkbox"/>	Size 4 Easy Seat Grab Posts (Pair)	LES/4/P	\$256
<input type="checkbox"/>	Size 1 Easy Seat Push Handle	LES/1/Q	\$150
<input type="checkbox"/>	Size 2 Easy Seat Push Handle	LES/2/Q	\$150
<input type="checkbox"/>	Size 3 Easy Seat Push Handle	LES/3/Q	\$150
<input type="checkbox"/>	Size 4 Easy Seat Push Handle	LES/4/Q	\$150
<input type="checkbox"/>	Size 1 Caster Legs (2 Fixed and 2 Swivel)	LES/1/H	\$300
<input type="checkbox"/>	Size 2 Caster Legs (2 Fixed and 2 Swivel)	LES/2/H	\$300
<input type="checkbox"/>	Size 3 Caster Legs (2 Fixed and 2 Swivel)	LES/3/H	\$300
<input type="checkbox"/>	Size 4 Caster Legs (2 Fixed and 2 Swivel)	LES/4/H	\$300

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3	Size 4
User Height (Min / Max)	30" / 40"	37" / 46"	40" / 52"	46" / 58"
Maximum User Weight	66 lbs	66 lbs	88 lbs	88 lbs
Seat Width (Min / Max)	5" / 9"	6" / 10"	9" / 13"	11" / 15"
Seat Depth (Min / Max)	7" / 10"	9" / 12"	10" / 14"	12" / 16"
Backrest Height (Min / Max)	9" / 11"	10" / 14"	13" / 17"	16" / 20"
Seat to Footrest (Min / Max)	6" / 9"	9" / 12"	10" / 14"	12" / 17"
Seat to Armrest (Min / Max)	4" / 6"	5" / 7"	6" / 8"	7" / 10"
Tray Dimensions (Length x Width)	18" L x 16.5" W		23" L x 20.5" W	

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3	Size 4
User Height (Min / Max)	760 mm / 1000 mm	930 mm / 1150 mm	1000 mm / 1300 mm	1150 mm / 1450 mm
Maximum User Weight	30 Kgs	30 Kgs	40 Kgs	40 Kgs
Seat Width (Min / Max)	125 mm / 225 mm	150 mm / 250 mm	225 mm / 325 mm	275 mm / 375 mm
Seat Depth (Min / Max)	175 mm / 250 mm	225 mm / 330 mm	250 mm / 350 mm	300 mm / 400 mm
Backrest Height (Min / Max)	225 mm / 285 mm	265 mm / 340 mm	325 mm / 425 mm	400 mm / 500 mm
Seat to Footrest (Min / Max)	150 mm / 230 mm	220 mm / 300 mm	250 mm / 350 mm	300 mm / 420 mm
Seat to Armrest (Min / Max)	110 mm / 150 mm	100 mm / 180 mm	160 mm / 210 mm	180 mm / 240 mm
Tray Dimensions (Length x Width)	450 mm L x 420 mm W		590 mm L x 520 mm W	