

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

 Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Effective 03/09/2020

Freestander

HCPC: E0638

Freestander

Selection required. Includes Chest and Hip Pad, Flexible Lateral, Lap Belt, PU Knee Supports, Footplate, Sandals, Tray Cover and Activity Bowl.

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Freestander - Meadow	LFS/1-MW-EXP	\$2491
<input type="checkbox"/>	Size 2 Freestander - Mondeo Blue	LFS/2-MB-EXP	\$2699
<input type="checkbox"/>	Size 3 Freestander - Mondeo Blue	LFS/3-MB-EXP	\$3053

Additional Positioning

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Lap Belt - Meadow	LFS/1/B-MW	\$181
<input type="checkbox"/>	Size 1 Lap Belt - Mondeo Blue	LFS/1/B-MB	\$167
<input type="checkbox"/>	Size 2 Lap Belt - Meadow	LFS/2/B-MW	\$230
<input type="checkbox"/>	Size 2 Lap Belt - Mondeo Blue	LFS/2/B-MB	\$188
<input type="checkbox"/>	Size 3 Lap Belt - Mondeo Blue	LFS/3/B-MB	\$204
<input type="checkbox"/>	Size 1 High Back Strap - Meadow	LFS/1/D-MW	\$113
<input type="checkbox"/>	Size 2 High Back Strap - Mondeo Blue	LFS/2/D-MB	\$121
<input type="checkbox"/>	Size 3 High Back Strap - Mondeo Blue	LFS/3/D-MB	\$134
<input type="checkbox"/>	Size 1 Solid Pelvic Plate - Meadow <i>Pelvic Band selection required</i>	LFS/1/L-MW	\$293
<input type="checkbox"/>	Size 2 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>	LFS/2/L-MB	\$331
<input type="checkbox"/>	Size 3 Solid Pelvic Plate - Mondeo Blue <i>Pelvic Band selection required</i>	LFS/3/L-MB	\$370
<input type="checkbox"/>	Size 1 Pelvic Band - Meadow	LPS/1/B-MW	\$172
<input type="checkbox"/>	Size 2 Pelvic Band - Mondeo Blue	LPS/2/B-MB	\$188
<input type="checkbox"/>	Size 3 Pelvic Band - Mondeo Blue	LPS/3/B-MB	\$204

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

LECKEY[®]

Additional Knee Support

Description	Item Code	Retail (Each)
<input type="checkbox"/> Size 1 Upholstered Knee Cups - Meadow	LPS/1/K-MW	\$117
<input type="checkbox"/> Size 2 Upholstered Knee Cups - Mondeo Blue	LPS/2/K-MB	\$117
<input type="checkbox"/> Size 3 Upholstered Knee Cups - Mondeo Blue	LPS/3/K-MB	\$117
<input type="checkbox"/> Size 4 Upholstered Knee Cups - Mondeo Blue	LPS/4/K-MB	\$117
Size 1 Knee Spacer (Each) <input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/1/F	\$60
Size 2 Knee Spacer (Each) <input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/2/F	\$60
Size 3 Knee Spacer (Each) <input type="checkbox"/> Quantity 1 <input type="checkbox"/> Quantity 2	LPS/3/F	\$60

Additional Foot Positioning

Sandal selection required with Sandal Raiser

Description	Item Code	Retail (Each)
<input type="checkbox"/> Small Sandals with Straps (Pair)	152-1600	\$241
<input type="checkbox"/> Medium Sandals with Straps (Pair)	152-2600	\$263
<input type="checkbox"/> Large Sandals with Straps (Pair)	152-3600	\$284
<input type="checkbox"/> Size 1 Sandal Raiser (Each) Quantity: _____ (1-6)	LFS/1/J	\$69
<input type="checkbox"/> Size 2 Sandal Raiser (Each) Quantity: _____ (1-6)	LFS/2/J	\$69
<input type="checkbox"/> Size 3 Sandal Raiser (Each) Quantity: _____ (1-6)	LFS/3/J	\$69

Headrest

Description	Item Code	Retail
<input type="checkbox"/> Headrest Interface Frame <i>Required with Contoured Headrest</i>	LIF/1	\$504
Headrest Bracket - Left and Right	LIF/7	\$123
<input type="checkbox"/> Size 1 Contoured Headrest - Meadow	LIF/2/A-MW	\$342
<input type="checkbox"/> Size 1 Contoured Headrest - Mondeo Blue	LIF/2/A-MB	\$342
<input type="checkbox"/> Size 1 Contoured Headrest - Vinyl	LIF/2/A	\$342
<input type="checkbox"/> Size 2 Contoured Headrest - Mondeo Blue	LIF/2/B-MB	\$342
<input type="checkbox"/> Size 2 Contoured Headrest - Vinyl	LIF/2/B	\$342
<input type="checkbox"/> Size 3 Contoured Headrest - Mondeo Blue	LIF/2/C-MB	\$342
<input type="checkbox"/> Size 3 Contoured Headrest - Vinyl	LIF/2/C	\$342

Accessories

Description	Item Code	Retail
<input type="checkbox"/> Size 1 Grab/Suspension Rail	LPS/1/C	\$192
<input type="checkbox"/> Size 2 Grab/Suspension Rail	LPS/2/C	\$192
<input type="checkbox"/> Size 3 Grab/Suspension Rail	LPS/3/C	\$192

LECKEY[®]

Technical Specifications

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55.1" / 72.8"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Pad Height (Min / Max)	20.7" / 33.5"	26.8" / 43.3"	38.4" / 56.1"
Width Between Chest Supports (Min / Max)	6" / 10"	8" / 12"	10" / 14"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Specification	Size 1	Size 2	Size 3
Age (approximate)	1 1/2 - 5 years	4 - 10 years	9 - 18 years
User Body Height (Min / Max)	800 mm / 1151 mm	1100 mm / 1448 mm	1400 mm / 1849 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Pad Height (Min / Max)	470 mm / 680 mm	680 mm / 975 mm	975 mm / 1350 mm
Width Between Chest Supports (Min / Max)	152 mm / 254 mm	203 mm / 305 mm	254 mm / 356 mm
Activity Tray	400 mm x 505 mm	450 mm x 575 mm	500 mm x 625 mm