

Company: \_\_\_\_\_

Ship To: \_\_\_\_\_

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Address: \_\_\_\_\_

 Quote  PO Number: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_

RTS: \_\_\_\_\_

Zip: \_\_\_\_\_ Phone: \_\_\_\_\_

Marked For: \_\_\_\_\_

Email/Fax: \_\_\_\_\_

Effective 03/31/2020

## Mygo Stander

**HCPC: E0641**
*Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs*

### Packages

*Package includes Mygo Stander, Stander Covers, Shoulder Pad and Cover, Contoured Headrest Hardware and Cushion - Black, Fixed Chest Laterals, Chest Wraparound Harness, Supine Derotation Hip Belt, Knee Cup Height Extension Brackets, Basic Knee Straps, Tray, Sandals with Straps and Product Assembly. Substitutions are not available for packages.*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Package - Blue <i>Includes Medium Sandals</i>	PMSTB1	E0641	\$7010
<input type="checkbox"/>	Size 1 Mygo Stander Package - Pink <i>Includes Medium Sandals</i>	PMSTP1	E0641	\$7010
<input type="checkbox"/>	Size 2 Mygo Stander Package - Blue <i>Includes Large Sandals</i>	PMSTB2	E0641	\$7113
<input type="checkbox"/>	Size 2 Mygo Stander Package - Pink <i>Includes Large Sandals</i>	PMSTP2	E0641	\$7113

### Mygo Stander

*Selection required. Includes manual angle adjustment, polyurethane chest pad, hip pad and knee supports, two piece footplate and chassis*

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	<b>Size 1 Mygo Stander</b>	143-1600		
	Mygo Stander Frame	143-1600 - 1	E0641	\$1511
	Pivot Base	143-1600 - 2		\$1200
	Chest Pad Assembly	143-1600 - 3		\$150
	Hip Assembly	143-1600 - 4		\$150
	Knee Supports	143-1600 - 5		\$200
	Positioning Pad Covers	143-1600 - 6		\$355
	Swivel Casters	143-1600 - 7		\$395
	Two Piece Footplate	143-1600 - 8		\$350
<input type="checkbox"/>	<b>Size 2 Mygo Stander</b>	143-2600		
	Mygo Stander Frame	143-2600 - 1	E0641	\$1511
	Pivot Base	143-2600 - 2		\$1200
	Chest Pad	143-2600 - 3		\$150
	Hip Pad	143-2600 - 4		\$150
	Knee Supports	143-2600 - 5		\$200
	Positioning Pad Covers	143-2600 - 6		\$355
	Swivel Casters	143-2600 - 7		\$395
	Two Piece Footplate	143-2600 - 8		\$350

*HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions.*

*Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

# LECKEY<sup>o</sup>

## Cover

**Selection required. Must match Stander size. Includes chest and hip pad covers and knee support covers.**

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Blue	143-1770-08	N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Orange	143-1770-07	N/C
<input type="checkbox"/>	Size 1 Mygo Stander Covers - Pink	143-1770-09	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Blue	143-2770-08	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Orange	143-2770-07	N/C
<input type="checkbox"/>	Size 2 Mygo Stander Covers - Pink	143-2770-09	N/C

## Shoulder Support

**Cover selection required with Shoulder Pad. Shoulder Support requires Head Support.**

	Description	Item Code	Retail
<input type="checkbox"/>	Mygo Stander Shoulder Pad	143-1638	\$345
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Blue	143-1768-08	N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Orange	143-1768-07	N/C
<input type="checkbox"/>	Mygo Stander Shoulder Pad Cover - Pink	143-1768-09	N/C
<input type="checkbox"/>	No Cover		N/C

## Head Support

**Cushion selection required with hardware. Head Support requires Shoulder Support.**

	Description	Item Code	Retail
<input type="checkbox"/>	Contoured Headrest Support Hardware	120-859	\$340
<input type="checkbox"/>	Contoured Headrest Cushion Infection Control - Black	143-1787	\$160
<input type="checkbox"/>	Contoured Headrest Cushion - Black	137-800-05	\$156
<input type="checkbox"/>	Flat Headrest Support Hardware	117-697	\$209
<input type="checkbox"/>	Flat Headrest Cushion - Black	117-850-05	\$106
<input type="checkbox"/>	Flat Headrest Laterals with Covers - Black (Pair) <i>Requires Flat Headrest. Not available with Contoured Headrest.</i>	120-893-05	\$298
<input type="checkbox"/>	Thoracic Headrest Adapter Kit <i>Required when selecting a headrest without Shoulder Support</i>	111214	\$335

# LECKEY<sup>o</sup>

## Wraparound Harness

*Harnesses may be used as chest harness or hip harness. If used as chest harness, chest laterals are required. If used as hip harness, hip laterals are required.*

Description	Quantity 1	Quantity 2	Item Code	Retail (Each)	Retail (Pair)
Size 1 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-08	\$248	\$496
Size 1 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-07	\$248	\$496
Size 1 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-1780-09	\$248	\$496
Size 1 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-1775-14	\$248	\$496
Size 2 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-08	\$255	\$510
Size 2 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-07	\$255	\$510
Size 2 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-2780-09	\$255	\$510
Size 2 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-2775-14	\$255	\$510
Size 3 Chest/Hip Wraparound Harness - Blue	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-08	\$262	\$524
Size 3 Chest/Hip Wraparound Harness - Orange	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-07	\$262	\$524
Size 3 Chest/Hip Wraparound Harness - Pink	<input type="checkbox"/>	<input type="checkbox"/>	143-3780-09	\$262	\$524
Size 3 Chest/Hip Wraparound Harness Infection Control - Black	<input type="checkbox"/>	<input type="checkbox"/>	143-3775-14	\$262	\$524

## Lateral Support

*Chest Laterals required if ordering Wraparound Harness to be used as chest harness. Hip Laterals required if ordering Wraparound Harness to be used as hip harness.*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Fixed Chest Laterals (Pair)	143-1619	\$149
<input type="checkbox"/>	Size 2 Fixed Chest Laterals (Pair)	143-2619	\$149
<input type="checkbox"/>	Size 1 and 2 Flip-Away Chest Laterals (Pair)	143-1675	\$480
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Right <i>Not available if ordering Wraparound Harness to be used as chest harness</i>	143-1674	\$281
<input type="checkbox"/>	Size 1 and 2 Chest Complex Lateral Vertical with Cover Left <i>Not available if ordering Wraparound Harness to be used as chest harness</i>	143-1673	\$281
<input type="checkbox"/>	Hip Laterals (Pair) <i>Not available with Derotation Belt</i>	143-1663	\$149

# LECKEY<sup>o</sup>

## Hip Belt

*Not available with Hip Laterals or Hip Wraparound Harness*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Blue	143-1798-08	\$281
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Orange	143-1798-07	\$281
<input type="checkbox"/>	Size 1 Prone Derotation Hip Belt with Hardware - Pink	143-1798-09	\$281
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Blue	143-2798-08	\$318
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Orange	143-2798-07	\$318
<input type="checkbox"/>	Size 2 Prone Derotation Hip Belt with Hardware - Pink	143-2798-09	\$318
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Blue	143-3798-08	\$351
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Orange	143-3798-07	\$351
<input type="checkbox"/>	Size 3 Prone Derotation Hip Belt with Hardware - Pink	143-3798-09	\$351
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Blue	143-1783-08	\$281
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Orange	143-1783-07	\$281
<input type="checkbox"/>	Size 1 Supine Derotation Hip Belt with Hardware - Pink	143-1783-09	\$281
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Blue	143-2783-08	\$318
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Orange	143-2783-07	\$318
<input type="checkbox"/>	Size 2 Supine Derotation Hip Belt with Hardware - Pink	143-2783-09	\$318
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Blue	143-3783-08	\$351
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Orange	143-3783-07	\$351
<input type="checkbox"/>	Size 3 Supine Derotation Hip Belt with Hardware - Pink	143-3783-09	\$351

## Hip Infill

	Description	Item Code	Retail
<input type="checkbox"/>	Hip Infill Pad	143-1666	\$308
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Blue	143-1799-08	\$129
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Orange	143-1799-07	\$129
<input type="checkbox"/>	Size 1 Hip Infill Pad Cover - Pink	143-1799-09	\$129
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Blue	143-2799-08	\$148
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Orange	143-2799-07	\$148
<input type="checkbox"/>	Size 2 Hip Infill Pad Cover - Pink	143-2799-09	\$148

**LECKEY**°

**Posterior Support**

*Selection of Posterior Support and Pommel requires selection of both Posterior Support Cover AND Pommel Support Cover. Not available with Hip Laterals or Hip Wraparound Harness.*

	Description	Item Code	Retail
<input type="checkbox"/>	Size 1 Posterior Support and Pommel	143-1624	\$730
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08	\$136
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07	\$136
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09	\$136
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14	\$136
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08	\$114
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07	\$114
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09	\$114
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14	\$114
<input type="checkbox"/>	Size 2 Posterior Support and Pommel	143-2624	\$768
<input type="checkbox"/>	Posterior Support Cover - Blue	143-1790-08	\$136
<input type="checkbox"/>	Posterior Support Cover - Orange	143-1790-07	\$136
<input type="checkbox"/>	Posterior Support Cover - Pink	143-1790-09	\$136
<input type="checkbox"/>	Posterior Support Infection Control Cover - Black	143-1790-14	\$136
<input type="checkbox"/>	Pommel Support Cover - Blue	143-1792-08	\$114
<input type="checkbox"/>	Pommel Support Cover - Orange	143-1792-07	\$114
<input type="checkbox"/>	Pommel Support Cover - Pink	143-1792-09	\$114
<input type="checkbox"/>	Pommel Support Infection Control Cover - Black	143-1792-14	\$114

**Knee**

	Description	Item Code	Retail
<input type="checkbox"/>	Knee Cup Height Extension Brackets (Pair) <i>Increases distance from footplate to top of knee by 4"</i>	143-1690	\$217
<input type="checkbox"/>	Size 1 Basic Knee Straps - Blue (Pair)	143-1776-08	\$169
<input type="checkbox"/>	Size 1 Basic Knee Straps - Orange (Pair)	143-1776-07	\$169
<input type="checkbox"/>	Size 1 Basic Knee Straps - Pink (Pair)	143-1776-09	\$169
<input type="checkbox"/>	Size 1 Basic Knee Straps Infection Control - Black (Pair)	143-1776-14	\$169
<input type="checkbox"/>	Size 2 Basic Knee Straps - Blue (Pair)	143-2776-08	\$207
<input type="checkbox"/>	Size 2 Basic Knee Straps - Orange (Pair)	143-2776-07	\$207
<input type="checkbox"/>	Size 2 Basic Knee Straps - Pink (Pair)	143-2776-09	\$207
<input type="checkbox"/>	Size 2 Basic Knee Straps Infection Control - Black (Pair)	143-2776-14	\$207
<input type="checkbox"/>	Split Knee Straps - Blue (Pair)	143-1784-08	\$312
<input type="checkbox"/>	Split Knee Straps - Orange (Pair)	143-1784-07	\$312
<input type="checkbox"/>	Split Knee Straps - Pink (Pair)	143-1784-09	\$312

# LECKEY<sup>o</sup>

## Accessories

	Description	Item Code	Retail	
<input type="checkbox"/>	Size 1 Mygo Stander Tray	143-1609	\$407	
<input type="checkbox"/>	Size 2 Mygo Stander Tray	143-2609	\$407	
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Bowl Insert	143-1633	\$421	
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Bowl Insert	143-2633	\$421	
<input type="checkbox"/>	Size 1 Mygo Stander Tray with Clear Insert	143-1610	\$463	
<input type="checkbox"/>	Size 2 Mygo Stander Tray with Clear Insert	143-2610	\$463	
<input type="checkbox"/>	Grab Rail/Toy Bar for Tray	143-1649	\$262	
<input type="checkbox"/>	Small Sandals with Straps (Pair)	152-1600	\$241	
<input type="checkbox"/>	Medium Sandals with Straps (Pair)	152-2600	\$263	
<input type="checkbox"/>	Large Sandals with Straps (Pair)	152-3600	\$284	
<input type="checkbox"/>	Sandal Raiser 1" (Each) <i>Requires Sandal selection</i>	Quantity: _____ (1-6)	143-1615	\$119

## Assembly

<input type="checkbox"/>	Product Assembly Fee		\$120
--------------------------	----------------------	--	-------

# LECKEY<sup>o</sup>

## Technical Specifications

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	37" / 57.1"	53.1" / 66.9"
Maximum User Weight	110 lbs	132 lbs
Chest Support Height	29.9" - 42.5"	38.6" - 51.2"
Chest Width	6.3" - 11"	7.9" - 12.6"
Hip Support Height	18.1" - 28.7"	24.4" - 34.4"
Knee Cup Height	11" - 15"	12.25" - 16.25"
Knee Cup Height with Extension Bracket	15" - 19"	16.25" - 20.25"
Hip Width	8.7" - 13"	11" - 15.4"
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0" - 6"	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	48 lbs	55 lbs
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	38.5" x 24.8"	44.5" x 25.6"

## Technical Specifications - Metric

Specification	Size 1	Size 2
Age (approximate)	4 - 10 years	8 - 14 years
User Height (Min / Max)	1080 mm / 1450 mm	1350 mm / 1700 mm
Maximum User Weight	50 kgs	60 kgs
Chest Support Height	760 mm - 1080 mm	980 mm - 1300 mm
Chest Width	160 mm - 280 mm	200 mm - 320 mm
Hip Support Height	460 mm - 730 mm	620 mm - 875 mm
Knee Cup Height	280 mm - 380 mm	310mm - 410 mm
Knee Cup Height with Extension Bracket	381 mm - 481 mm	411 mm - 511 mm
Hip Width	220 mm - 330 mm	280 mm - 390 mm
Footplate Angle: Plantarflexion / Dorsiflexion	10° / 10°	
Maximum Knee Contracture	25°	
Tray Height Range	0 mm - 150 mm	
Tray Angle Adjustment: Prone / Supine	45° / 45°	
Product Weight	22 kgs	25 kgs
Product Angle Range: Prone / Supine	0° - 80°	
Stander Footprint (length x width)	978 mm - 630 mm	1130 mm x 650 mm