



ORDER FORM

Horizon Stander

Robust and durable 3-in-1 standing frame for both young children and adult users.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Horizon Stander

Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 220 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Frame			\$
Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals. (Knee supports and Footplates separate)			
<input type="checkbox"/> LHZ/2/198-GRB-EUR	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E1399	\$1,910
<input type="checkbox"/> LHZ/3/198-GRB-EUR	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E1399	\$1,990
Base			\$
Selection Required. MUST match Frame size. Includes worldwide adapters.			
<input type="checkbox"/> LHZ/2/177-EXP	Size 2 Horizon Base/Chassis	E0641	\$3,528
<input type="checkbox"/> LHZ/3/177-EXP	Size 3 Horizon Base/Chassis	E0641	\$4,403
Headrest			\$
Head Support Assembly selection required if Head Lateral Support/Extension Bracket is selected.			
<input type="checkbox"/> LHZ/2/169	Size 2 and 3 Head Support Assembly - Black	E1399	\$287
<input type="checkbox"/> LHZ/2/117	Size 2 and 3 Head Lateral Supports - Black (Pair)	E1399	\$248
<input type="checkbox"/> 101-1959	Headrest 4" Horizontal Extension Bracket	E1399	\$138
<input type="checkbox"/> LHZ/1/180	Size 1-3 Contoured Head Support	E1399	\$628
Hip Laterals			\$
Selection required if Hip Belt is not selected.			
<input type="checkbox"/> LHZ/2/164	Size 2 Horizon Hip Laterals - Black (Pair)		\$368
<input type="checkbox"/> LHZ/3/164	Size 3 Horizon Hip Laterals - Black (Pair)		\$383
Hip Belts			\$
Selection required if Hip Laterals are not selected. Multiple Selections can be made.			
<input type="checkbox"/> LHZ/2/181	Size 2 Supine Belt - Black		\$258
<input type="checkbox"/> LHZ/3/181	Size 3 Supine Belt - Black		\$284
Knee Supports			\$
Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" adjustment. Long Knee Bracket has 7" - 9" adjustment.			
<input type="checkbox"/> LHZ/2/167	Size 2 and 3 Standard Knee Assembly - Black (Pair)		\$516
<input type="checkbox"/> 101-2224	Medium Knee Adjustment Bracket (Pair)		\$119
<input type="checkbox"/> 101-2225	Long Knee Adjustment Bracket (Pair)		\$122
Additional Foot Positioning			\$
<input type="checkbox"/> 152-1600	Small Sandals with Straps (Pair)		\$273
<input type="checkbox"/> 152-2600	Medium Sandals with Straps (Pair)		\$298
<input type="checkbox"/> 152-3600	Large Sandals with Straps (Pair)		\$321
<input type="checkbox"/> LHZ/2/161	Size 2 and 3 Flip Up Footplate Assembly (Pair)		\$815
<input type="checkbox"/> LHZ/3/160	Size 1-3 3D Footplate Assembly (Pair) - One size fits all		\$1,332
Trays			\$
Must match Horizon Stander Size			
<input type="checkbox"/> LHZ/2/143	Size 2 and 3 Supine Tray		\$432
<input type="checkbox"/> LHZ/2/158	Size 2 and 3 Prone Tray		\$572

Technical Specifications - Imperial

Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

Technical Specifications - Metric

Frame Sizing Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm

Sunrise Medical (US) LLC

2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel.: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com

