



# Horizon Stander

January 2022

**Mark For:**

Date: \_\_\_\_\_

Dealer Acct #: \_\_\_\_\_

Dealer: \_\_\_\_\_

Dealer Contact: \_\_\_\_\_

Dealer Address: \_\_\_\_\_

Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Dealer Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Confirmation Email: \_\_\_\_\_

Confirm Via:  Fax  Email

**Submitting for:****Quote****Order**

PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_

Attention: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_

Ship To Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

**Horizon Stander****HCPCS: E0641****Frame**

*Selection Required. Includes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals, Standard Knees and Basic Footplate.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/178BBEXP	Size 2 Horizon Support Frame - Blue	E1399	\$2,670
<input type="checkbox"/>	LHZ/2/178PBEXP	Size 2 Horizon Support Frame - Pink	E1399	\$2,670
<input type="checkbox"/>	LHZ/3/178BBEXP	Size 3 Horizon Support Frame - Blue	E1399	\$3,302
<input type="checkbox"/>	LHZ/3/178PBEXP	Size 3 Horizon Support Frame - Pink	E1399	\$3,302

**Base**

*Selection Required. Must match Frame size. Includes Worldwide Adapters.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/177-EXP	Size 2 Horizon Base/Chassis	E0641	\$3,283
<input type="checkbox"/>	LHZ/3/177-EXP	Size 3 Horizon Base/Chassis	E0641	\$4,097

**Headrest**

*Head Support Assembly or Extended Head Support Assembly selection required if Head Lateral Supports is selected. Head Support Assembly Size must match Head Lateral Supports size.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/169-BB	Size 2 and 3 Head Support Assembly - Blue	E1399	\$267
<input type="checkbox"/>	LHZ/2/169-PB	Size 2 and 3 Head Support Assembly - Pink	E1399	\$267
<input type="checkbox"/>	LHZ/2/117-BB	Size 2 and 3 Head Lateral Supports - Blue (Pair)	E1399	\$230
<input type="checkbox"/>	LHZ/2/117-PB	Size 2 and 3 Head Lateral Supports - Pink (Pair)	E1399	\$230
<input type="checkbox"/>	LHZ/2/187-BB	Size 2 and 3 Extended 4" Additional Depth Head Support Assembly - Blue	E1399	\$358
<input type="checkbox"/>	LHZ/2/187-PB	Size 2 and 3 Extended 4" Additional Depth Head Support Assembly - Pink	E1399	\$358
<input type="checkbox"/>	LHZ/1/180	Size 1-3 Contoured Head Support	E1399	\$585

**Additional Chest Supports**

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/168-BB	Size 2 Flexible Chest Laterals - Blue (Pair)		\$480
<input type="checkbox"/>	LHZ/2/168-PB	Size 2 Flexible Chest Laterals - Pink (Pair)		\$480
<input type="checkbox"/>	LHZ/3/168-BB	Size 3 Flexible Chest Laterals - Blue (Pair)		\$536
<input type="checkbox"/>	LHZ/3/168-PB	Size 3 Flexible Chest Laterals - Pink (Pair)		\$536
<input type="checkbox"/>	LHZ/3/154-BB	Size 3 Flip Away Chest Laterals - Blue (Pair)		\$589
<input type="checkbox"/>	LHZ/3/154-PB	Size 3 Flip Away Chest Laterals - Pink (Pair)		\$589

**Hip Laterals**

*Selection required if Hip Belt is not selected*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/164-BB	Size 2 Horizon Hip Laterals - Blue (Pair)		\$343
<input type="checkbox"/>	LHZ/2/164-PB	Size 2 Horizon Hip Laterals - Pink (Pair)		\$343
<input type="checkbox"/>	LHZ/3/164-BB	Size 3 Horizon Hip Laterals - Blue (Pair)		\$356
<input type="checkbox"/>	LHZ/3/164-PB	Size 3 Horizon Hip Laterals - Pink (Pair)		\$356
<input type="checkbox"/>	LHZ/2/162-BB	Size 3 Horizon Hip Laterals for Size 2 Frame - Blue (Pair)		\$356
<input type="checkbox"/>	LHZ/2/162-PB	Size 3 Horizon Hip Laterals for Size 2 Frame - Pink (Pair)		\$356
<input type="checkbox"/>	LHZ/3/153-BB	Size 3 Flip Away Hip Laterals - Blue (Pair)		\$477
<input type="checkbox"/>	LHZ/3/153-PB	Size 3 Flip Away Hip Laterals - Pink (Pair)		\$477

**Hip Belts**

*Selection required if Hip Laterals are not selected . Multiple selections can be made.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/181-BB	Size 2 Supine Belt - Blue		\$205
<input type="checkbox"/>	LHZ/2/181-PB	Size 2 Supine Belt - Pink		\$205
<input type="checkbox"/>	LHZ/3/181-BB	Size 3 Supine Belt - Blue		\$226
<input type="checkbox"/>	LHZ/3/181-PB	Size 3 Supine Belt - Pink		\$226
<input type="checkbox"/>	LHZ/2/182-BB	Size 2 Prone Belt - Blue		\$205
<input type="checkbox"/>	LHZ/2/182-PB	Size 2 Prone Belt - Pink		\$205
<input type="checkbox"/>	LHZ/3/182-BB	Size 3 Prone Belt - Blue		\$226
<input type="checkbox"/>	LHZ/3/182-PB	Size 3 Prone Belt - Pink		\$226

**Additional Knee Supports**

*Standard Knee Assembly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Extended Knee Assembly has 5" - 7" adjustment and Long Extended Knee Assembly has 7" - 9" adjustment.*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/167-BB	Size 2 and 3 Standard Knee Assembly - Blue (Pair)		\$480
<input type="checkbox"/>	LHZ/2/167-PB	Size 2 and 3 Standard Knee Assembly - Pink (Pair)		\$480
<input type="checkbox"/>	LHZ/2/183-BB	Size 2 and 3 Medium Extended Knee Assembly - Blue (Pair)		\$632
<input type="checkbox"/>	LHZ/2/183-PB	Size 2 and 3 Medium Extended Knee Assembly - Pink (Pair)		\$632
<input type="checkbox"/>	LHZ/2/184-BB	Size 2 and 3 Long Extended Knee Assembly - Blue (Pair)		\$669
<input type="checkbox"/>	LHZ/2/184-PB	Size 2 and 3 Long Extended Knee Assembly - Pink (Pair)		\$669
<input type="checkbox"/>	LHZ/2/185	Size 1-3 Medium Extended Knee Brackets Only (Pair)		\$532

**Lower Leg Positioning**

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	152-1600	Small Sandals with Straps (Pair)		\$254
<input type="checkbox"/>	152-2600	Medium Sandals with Straps (Pair)		\$277
<input type="checkbox"/>	152-3600	Large Sandals with Straps (Pair)		\$299
<input type="checkbox"/>	LHZ/2/161	Size 2 and 3 Flip Up Footplate Assembly (Pair)		\$758
<input type="checkbox"/>	LHZ/3/160	Size 1-3 3D Footplate Assembly (Pair) <i>One size fits all</i>		\$1,239

**Accessories***Must match Horizon Stander Size*

	Item Code	Description	HCPCS	Retail
<input type="checkbox"/>	LHZ/2/143	Size 2 and 3 Supine Tray		\$402
<input type="checkbox"/>	LHZ/2/158	Size 2 and 3 Prone Tray		\$532
<input type="checkbox"/>	LHZ/2/192	Size 2 and 3 Horizon Grab Rail		\$264
<input type="checkbox"/>	LHZ/2/193	Size 2 and 3 Horizon Grab Posts (Pair)		\$264

**Technical Specifications - Imperial**

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	50.4" x 25.2" x 31.7"	50.4" x 25.2" x 31.7"
Weight of Product	110.6 lbs	130 lbs
Maximum User Weight	176 lbs	220 lbs
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	7"	7"

**Technical Specifications - Metric**

Specification	Size 2	Size 3
Age Range	9 -- 18	14 -- Adult
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm
Weight of Product	50 kgs	95 kgs
Maximum User Weight	80 kgs	100 kgs
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm
Width Between Chest Laterals (Min / Max)	240 mm / 380 mm	340 mm / 460 mm
Footplate Angle Adjustable + / -	10°	10°
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm

