



REACH SE™ U.S. Order Form



June 2020

Date: _____ QUOTE ORDER - PO# _____

Company Name: _____

Ship To: _____

Contact Name: _____

Attention: _____

ATP Name: _____

Address: _____

Billing Address: _____

City: _____ ST: _____ ZIP: _____

City: _____ ST: _____ ZIP: _____

Phone: _____ Fax: _____

Email: _____

Mark for: _____

Model shown with accessories.



REACH SE™ INCLUDES:

- Fixed Seat Tilt 10°, 20°, 30°
- Hip Angle Adjustment
- QR 12" Poly Rear Wheels
- 7" Poly Front Casters
- Height Adjustable Push Handle
- Height Adjustable Foot Plate
- Foot Activated Hub Lock
- Under Seat Storage Bag
- Sunshade

Specifications & Adjustability

	REACH™ 12	REACH™ 14	REACH™ 16
Max Occupant Weight	80 lbs/36 kg	110 lbs/50 kg	150 lbs/68 kg
Seat Width	8-12 in/20-30 cm	10-14 in/25-35 cm	12-16 in/30-40 cm
Seat Depth	8-13 in/20-33 cm	11-16 in/28-40 cm	13-18 in/33-46 cm
Back Height	11-18 in/28-46 cm	13-21 in/33-53 cm	16-24 in/40-61 cm
Seat to Footplate	8-13 in/20-33 cm	11-16 in/28-40 cm	13-18 in/33-46 cm
Seat to Floor	19.5 in/50 cm	19.5 in/50 cm	19.5 in/50 cm
Fixed Tilt Positions	10, 20, or 30°		
Seat to Back Angle (Range depends on seat angle)			
10° seat angle	85, 90, 95, 100, 105, 110, 115°		
20° seat angle	85, 90, 95, 100, 105°		
30° seat angle	85, 90, 95°		
Weight of Chair	27 lbs/12 kg	29 lbs/13 kg	31 lbs/14 kg
Chair Width	21.5 in/55cm	23.5 in/60 cm	25.5 in/65 cm
Folded Dimensions	31 x 21.5 x 18 in 79 x 55 x 46 cm	31 x 23.5 x 18 in 79 x 60 x 46 cm	33 x 25.5 x 18 in 84 x 65 x 46 cm

MEASUREMENT WORKSHEET

A. Seat Depth _____

B. Trunk Depth _____

C. Trunk Width _____

D. Hip Width _____

E. Seat to Top of Shoulder _____

F. Lower Leg Length _____

G. Seat to Axilla _____



REACH™ Chair

REACH SE™

<input type="checkbox"/>	3004-1004	E1238	REACH SE™ 12 - 12 in./30.5 cm width	\$1,790.00
<input type="checkbox"/>	3004-1005	E1238	REACH SE™ 14 - 14 in./35.5 cm width	\$1,790.00
<input type="checkbox"/>	3004-1006	E1238	REACH SE™ 16 - 16 in./40.5 cm width	\$1,790.00

Upholstery Color

Upholstery color applies to back and seat cushions, sunshade.

<input type="checkbox"/>	3004-0101	BIG SKY BLUE - Blue/Black	N/C
<input type="checkbox"/>	3004-0102	ROSSO RED - Red/Black	N/C
<input type="checkbox"/>	3004-0103	ICE SILVER - Silver/Black	N/C
<input type="checkbox"/>	3004-0104	SNOWBERRY PINK - Pink/Black	N/C
<input type="checkbox"/>	3004-0105	HUCKLEBERRY PURPLE - Purple/Black	N/C

Transit Option & Transit Harness Selection

Transit Option

Requires belt/harness, back rest and a head rest.

<input type="checkbox"/>	3004-7001	Transit Option (RESNA WC-19)	\$250.00
--------------------------	-----------	------------------------------	----------

This option is not retrofittable.

Transit Belt or Harness

Transit 5-Point Harness NOT COMPATIBLE with Reach 16

<input type="checkbox"/>	3004-7002	K0108	Transit 5-Point Harness for 12/14	\$165.00
		<input type="checkbox"/>	Add pad set for 5-point harness (3002-7022)	\$90.00
<input type="checkbox"/>	3001-7008	K0108	Transit Sure-Lok Lap Belt for 12/14	\$90.00
<input type="checkbox"/>	3004-7004	K0108	Transit Sure-Lok Lap Belt for 16	\$90.00

Select 5-PT Harness for occupants less than 50 lbs./23kg

Select Sure-Lok Lap Belt for occupants over 50 lbs./23kg

The transit option includes 4 installed transportation anchors & has been tested in compliance with WC-19 standards. For more information on RESNA WC-19 please visit www.wc-transportation-safety.com.

Seat Pan and Front Rigging

A. Standard Front Rigging and Solid Seat Pan

Select both items to order the standard version

<input type="checkbox"/>	3004-3001	E2231	Solid Seat Pan	\$165.00
<input type="checkbox"/>	3004-7032	E0990	Angle Adjustable Legrest	\$255.00

B. Swing-Away Leg Hangers and Solid Seat Pan

Seat Pan 3004-3002 Required for use with Swing-Away Leg Hangers

<input type="checkbox"/>	3004-3002	E2231	Solid Seat Pan for Swing-Aways	\$285.00
--------------------------	-----------	-------	--------------------------------	----------

Select rigging style, then select length from table to the right.

<input type="checkbox"/>	3004-5007		70° Swing Away Leg Hangers	n/c
<input type="checkbox"/>	3004-5003		80° Swing Away Leg Hangers	n/c
<input type="checkbox"/>	3004-5009	E0990	Swing Away ELR	\$255.00

Option B is Not Compatible with Reach 12

Front Rigging Length, Foot Plate, Calf Pad (only if B is selected)

Append the rigging length combination's two digit part number extension to the rigging part number.

70°	Short Extension	Standard Extension
Short Upper*	10-14"/25.4-35.56 cm -11 <input type="checkbox"/>	11-16"/27.94-40.64 -12 <input type="checkbox"/>
Medium Upper	11.5-16"/29.21-40 cm -21 <input type="checkbox"/>	11.5-18"/19.21-45.72 -22 <input type="checkbox"/>
Long Upper	n/a	13-19"/33-48.26 cm -32 <input type="checkbox"/>

80°	Standard Extension
Short Upper*	10-14"/25-35.56 cm -12 <input type="checkbox"/>
Medium Upper	11.5-14"/29.21-35.56 -22 <input type="checkbox"/>
Long Upper	11-17"/28-43 cm -32 <input type="checkbox"/>

Certain measurements may cause caster interference

Foot Plate and Front Rigging Accessories

<input type="checkbox"/>	3004-5199		Omit Foot Plate	n/c
<input type="checkbox"/>	3004-5010	K0040	Flip Up, Height/Angle Adjustable, Alum.	\$185.00
<input type="checkbox"/>	3004-5011	E0951	Heel Loop for Foot Plate	\$40.00
<input type="checkbox"/>	3004-5013	E0995	Calf Pad, Left	\$40.00
<input type="checkbox"/>	3004-5014	E0995	Calf Pad, Right	\$40.00

**Calf Pad Not Compatible with Short Upper Extension*

Ankle Supports

<input type="checkbox"/>	3001-6007	E0951	LEFT, Extra Small (Pad 2x8"), REACH 12/14	\$72.00
<input type="checkbox"/>	3001-6008	E0951	RIGHT, Extra Small (Pad 2x8"), REACH 12/14	\$72.00
<input type="checkbox"/>	3001-6009	E0951	LEFT, Small (Pad 3x10"), REACH 14/16	\$72.00
<input type="checkbox"/>	3001-6010	E0951	RIGHT, Small (Pad 3x10"), REACH 14/16	\$72.00

Small ankle supports are compatible with REACH 14 and 16 ONLY

XS ankle supports are NOT compatible with swing away leg hangers

Seat Cushion

LT Seat Cushion				
<input type="checkbox"/>	3002-7099		Omit Seat Cushion	N/C
<input type="checkbox"/>	3005-3001	E2292	REACH™ 12 Planar Seat Cushion	\$295.00
<input type="checkbox"/>	3005-3002	E2292	REACH™ 14 Planar Seat Cushion	\$295.00
<input type="checkbox"/>	3005-3003	E2292	REACH™ 16 Planar Seat Cushion	\$295.00

EVA Positioning Wedges for Seat and/or Back

<input type="checkbox"/>	3001-7011	K0108	Wedges QTY 4 and Lumbar QTY 1	\$50.00
--------------------------	-----------	-------	-------------------------------	---------

Lateral Pelvic/Thigh Supports

<input type="checkbox"/>	3005-6010	E0956	Reach Adj. Thigh Support, LEFT	\$120.00
<input type="checkbox"/>	3005-6011	E0956	Reach Adj. Thigh Support, RIGHT	\$120.00

Armrests & Positioning Tray

<input type="checkbox"/>	3004-6009	E0973	Reach Height/Angle Adj. Armrests	\$195.00
--------------------------	-----------	-------	----------------------------------	----------

<input type="checkbox"/>	3004-6113	E0950	REACH™ 12 Clear Positioning Tray	\$150.00
<input type="checkbox"/>	3004-6114	E0950	REACH™ 14 Clear Positioning Tray	\$150.00
<input type="checkbox"/>	3004-6115	E0950	REACH™ 16 Clear Positioning Tray	\$150.00
<input type="checkbox"/>	3004-6008	E1028	Tray Hardware, Adj./Removable	\$209.00

3004-6009 Armrests & 3004-6008 Hardware are required when ordering Tray Platform

Anterior Posture Supports & Pelvic Positioning

SP™ Pelvic Positioning Belts

<input type="checkbox"/>	IFPB5141	E0978	4-Point, Push Button 1", XS	\$130.00
<input type="checkbox"/>	IFPB5242	E0978	4-Point, Push Button 1.5", SM	\$130.00
<input type="checkbox"/>	IFPB5243	E0978	4-Point, Push Button 1.5", M	\$130.00
<input type="checkbox"/>	IFPB5111	E0978	2-Point, Push Button 1", XS	\$86.00
<input type="checkbox"/>	IFPB5112	E0978	2-Point, Push Button 1", SM	\$86.00

Abductor/Medial Knee Block

Must select hardware & one size of abductor pad. Wedge sizes are width x depth x height.

<input type="checkbox"/>	3005-7001	E1028	Leggero Flip Down Medial Hdw.	\$240.00
<input type="checkbox"/>	MTMD-12	E0957	Small Wedge Pad (3x5x3")	\$150.00
<input type="checkbox"/>	MTMD-13	E0957	Medium Wedge Pad (5x5x3")	\$150.00

Harness/Belt Accessories

<input type="checkbox"/>	3001-6011	K0108	Strap Risers, Pair for 1.25" Webbing	\$45.00
<input type="checkbox"/>	IFMDBC10		SP™ Buckle Cover Kit 1" Buckle	\$25.00
<input type="checkbox"/>	IFMDBC15		SP™ Buckle Cover Kit 1.5" Buckle	\$25.00

Anterior Posture Supports

<input type="checkbox"/>	IFPS1101-02	E0960	SP™ Classic, Dynamic, XS	\$142.00
<input type="checkbox"/>	IFPS1102-02	E0960	SP™ Classic, Dynamic, S	\$142.00
<input type="checkbox"/>	IFPS1103-02	E0960	SP™ Classic Dynamic, M	\$142.00
<input type="checkbox"/>	IFPS1201-02	E0960	SP™ Classic, Structured, XS	\$142.00
<input type="checkbox"/>	IFPS1202-02	E0960	SP™ Classic, Structured, S	\$142.00
<input type="checkbox"/>	IFPS1203-02	E0960	SP™ Classic, Structured, M	\$142.00
<input type="checkbox"/>	IFPS2101-02	E0960	SP™ Contoured, Dynamic, XS	\$142.00
<input type="checkbox"/>	IFPS2102-02	E0960	SP™ Contoured, Dynamic, S	\$142.00
<input type="checkbox"/>	IFPS2103-02	E0960	SP™ Contoured, Dynamic, M	\$142.00
<input type="checkbox"/>	IFPS2201-02	E0960	SP™ Contoured, Structured, XS	\$142.00
<input type="checkbox"/>	IFPS2202-02	E0960	SP™ Contoured, Structured, S	\$142.00
<input type="checkbox"/>	IFPS2203-02	E0960	SP™ Contoured, Structured, M	\$142.00
<input type="checkbox"/>	3005-6016	E0960	Leggero 'H' Harness, S	\$165.00
<input type="checkbox"/>	3005-6017	E0960	Leggero 'H' Harness, M	\$165.00
<input type="checkbox"/>	3005-6018	E0960	Leggero 'H' Harness, L	\$165.00
<input type="checkbox"/>	3005-6014-12	E0960	5-Point Positioning Harness	\$165.00
<input type="checkbox"/>	3005-6014-14	E0960	5-Point Positioning Harness	\$165.00
<input type="checkbox"/>	3005-6014-16	E0960	5-Point Positioning Harness	\$165.00

<input type="checkbox"/>	Add pad set for 5-point harness (3002-7022)	\$90.00
--------------------------	---	---------

Back Rest & Accessories

Select a system from options A, B, or C. Then continue your selections to the right to complete the package.

A. Reach Planar Back

<input type="checkbox"/>	3005-4004	E2291	REACH™ 12 Planar Back	\$265.00
<input type="checkbox"/>	3005-4005	E2291	REACH™ 14 Planar Back	\$265.00
<input type="checkbox"/>	3005-4006	E2291	REACH™ 16 Planar Back	\$265.00
<input type="checkbox"/>	3003-4099		Omit Back	N/C

When omitting back, crossbars between canes will remain in place.

B. LT Contoured Back

<input type="checkbox"/>	3005-4001	E2293	REACH™ 12 Contoured Adj. Back	\$525.00
<input type="checkbox"/>	3005-4002	E2293	REACH™ 14 Contoured Adj. Back	\$525.00
<input type="checkbox"/>	3005-4003	E2293	REACH™ 16 Contoured Adj. Back	\$525.00

Soft Lateral w/ LT 'H' Harness Options (Only available with B)

<input type="checkbox"/>	3005-6005	E0960	Soft Laterals w/ SM 'H' Harness	\$165.00
<input type="checkbox"/>	3005-6006	E0960	Soft Laterals w/ MED 'H' Harness	\$165.00
<input type="checkbox"/>	3005-6007	E0960	Soft Laterals w/ LRG 'H' Harness	\$165.00

Soft Laterals' are pads that integrate with back and harness

LT Lumbar Support (Only available with B)

<input type="checkbox"/>	3005-6002	K0108	REACH™ 12 Adj. Lumbar Support	\$68.00
<input type="checkbox"/>	3005-6003	K0108	REACH™ 14 Adj. Lumbar Support	\$68.00
<input type="checkbox"/>	3005-6004	K0108	REACH™ 16 Adj. Lumbar Support	\$68.00

C. REACH RCB™ Rigid Contoured Back

Swing Away Lateral Hardware & Pads (Only available with C)

Must select hardware & pad. Leggero Swing In/Out Lateral includes rail and track slider.

REACH RCB™ 12 Contoured Back Options (includes hardware, plate, and cushion)

<input type="checkbox"/>	3004-4212	E2293	12" W x 12" H (Adj. 12-14")	\$525.00
<input type="checkbox"/>	3004-4214	E2293	12" W x 14" H (Adj. 14-16")	\$525.00
<input type="checkbox"/>	3004-4216	E2293	12" W x 16" H (Adj. 16-18")	\$525.00
<input type="checkbox"/>	3004-4218	E2293	12" W x 18" H (Fixed Height)	\$525.00

REACH RCB™ 14 Contoured Back Options (includes hardware, plate, and cushion)

<input type="checkbox"/>	3004-4414	E2293	14" W x 14" H (Adj. 14-16")	\$525.00
<input type="checkbox"/>	3004-4416	E2293	14" W x 16" H (Adj. 16-18")	\$525.00
<input type="checkbox"/>	3004-4418	E2293	14" W x 18" H (Adj. 18-20")	\$525.00
<input type="checkbox"/>	3004-4420	E2293	14" W x 20" H (Adj. 20-21")	\$525.00

REACH RCB™ 16 Contoured Back Options (includes hardware, plate, and cushion)

<input type="checkbox"/>	3004-4616	E2293	16" W x 16" H (Adj. 16-18")	\$525.00
<input type="checkbox"/>	3004-4618	E2293	16" W x 18" H (Adj. 18-20")	\$525.00
<input type="checkbox"/>	3004-4620	E2293	16" W x 20" H (Adj. 20-22")	\$525.00
<input type="checkbox"/>	3004-4622	E2293	16" W x 22" H (Adj. 22-24")	\$525.00
<input type="checkbox"/>	3004-4624	E2293	16" W x 24" H (Fixed Height)	\$525.00

<input type="checkbox"/>	3004-6111	E1028	LEFT, Leggero Swing In/Out Lateral	\$225.00	
<input type="checkbox"/>	3004-6112	E1028	RIGHT, Leggero Swing In/Out Lateral	\$225.00	
L	R	3004-6035	E0956	Lateral Pad 3x5", Each	\$75.00
L	R	3004-6045	E0956	Lateral Pad 4x5", Each	\$75.00
L	R	3004-6046	E0956	Lateral Pad 4x6", Each	\$75.00

Leggero and Stealth Products Head Positioning

Leggero Head Positioning

<input type="checkbox"/>	3005-6008	E0966	LT Planar Rectangle Headrest	\$200.00
<input type="checkbox"/>	3001-6005	E0955	8" Oval Headrest	\$240.00
<input type="checkbox"/>	3002-6003	E0955	RESPOND™ Cradle Headrest and Hardware	\$240.00
<input type="checkbox"/>	3005-6099		Omit Head Positioning	N/C

Hardware is included with LT, 8" Multi-Adjustable & RESPOND™ Headrests

A. Link Headrest Hardware 1" Ball

<input type="checkbox"/>	TWB480-1-01	E1028	Link Hardware 1" Ball	\$231.00
--------------------------	-------------	-------	-----------------------	----------

B. Link Headrest Hardware 17mm Ball/Flip Down

<input type="checkbox"/>	TWB480-17-01	E1028	Link Hdw. 17mm Ball/Flip Down	\$320.00
--------------------------	--------------	-------	-------------------------------	----------

Comfort Plus Headrest (Only if A is selected)

<input type="checkbox"/>	CP250-SR	E0955	Comfort Plus 10"	\$152.00
<input type="checkbox"/>	CP450-SR	E0955	Comfort Plus 8"	\$152.00

i2i Head Positioning System (Only if B is selected)

<input type="checkbox"/>	i2i250	K0108	i2i Low Pro Nino XS, Pad Only	\$348.00
<input type="checkbox"/>	i2i350	K0108	i2i Low Pro Nino SM, Pad Only	\$348.00
<input type="checkbox"/>	i2i450	K0108	i2i Low Pro Adult SM, Pad Only	\$348.00

Chin Prompts (Only if i2i Head Positioning is selected)

<input type="checkbox"/>	i2i460	K0108	i2i SM Chin Prompt	\$80.00
<input type="checkbox"/>	i2i560	K0108	i2i MED Chin Prompt	\$80.00

Combo Series Headrest

<input type="checkbox"/>	CB315	E0955	Combo Series, Small	\$295.00
<input type="checkbox"/>	CB305	E0955	Combo Series, Medium	\$295.00
<input type="checkbox"/>	CB325	E0955	Combo Series, Large	\$295.00

Accessories

<input type="checkbox"/>	3004-7020		Rain Cover, REACH™ 12	\$195.00
<input type="checkbox"/>	3004-7021		Rain Cover, REACH™ 14	\$195.00
<input type="checkbox"/>	3004-7022		Rain Cover, REACH™ 16	\$195.00

<input type="checkbox"/>	3004-7011	K0108	Calf Support, REACH™ 12	\$75.00
<input type="checkbox"/>	3004-7027	K0108	Calf Support, REACH™ 14	\$75.00
<input type="checkbox"/>	3004-7028	K0108	Calf Support, REACH™ 16	\$75.00

May interfere with folding

<input type="checkbox"/>	3004-5002		Custom set lower leg length to (F) dimension	N/C
--------------------------	-----------	--	--	-----

When 3004-5002 is selected, most leg lengths can be accommodated.

<input type="checkbox"/>	3004-7012		Grab Bar, REACH™ 12	\$135.00
<input type="checkbox"/>	3004-7014		Grab Bar, REACH™ 14	\$135.00
<input type="checkbox"/>	3004-7016		Grab Bar, REACH™ 16	\$135.00

<input checked="" type="checkbox"/>	3004-7009	E0971	Rear Anti-Tips (required)	\$140.00
<input type="checkbox"/>	3004-7025	K0105	IV Pole	\$150.00
<input type="checkbox"/>	3004-7026	E2208	Oxygen Tank Holder (D size)	\$175.00

The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions and are not to be deemed as legal advice. The use of HCPCS CODES does not ensure coverage or payment. Appropriate billing codes are the responsibility of DME providers. For coverage information, verify the appropriate payer's coverage policy.