



TRAK SE™ U.S. Order Form



August 2019

Date: _____ QUOTE ORDER - PO# _____

Account: _____ Ship To: _____

Contact Name: _____ Attention: _____

ATP Name: _____ Address: _____

Billing Address: _____

City: _____ ST: _____ ZIP: _____

Phone: _____ Fax: _____

City: _____ ST: _____ ZIP: _____

Email: _____

Mark for ATP: _____

TRAK SE™ INCLUDES:



Model shown with some accessories


- Integrated Suspension
- 45° Tilt-in-Space Frame
- 30° Hip Angle Adjustment
- QR 16" Poly Rear Wheels
- 7" Poly Front Casters
- Height Adjustable Push Handle
- Height Adjustable Foot Plate
- Foot Activated Wheel Lock
- Caster Swivel Locks
- Under Seat Storage Bag
- Sunshade

Specifications

Maximum Chair Capacity	(TRAK 12 & TRAK 14) 80 lbs / 36 kg	
Push Handle Height	40-50 in / 101-127 cm	
Weight	29 lbs / 13 kg	
Rear Wheel Diameter	16 in / 41 cm	
Front Caster Diameter	7 in / 25 cm	
Reversing Width	41 in / 104 cm	
Maximum Safe Slope	5°	
	TRAK™ 12	TRAK™ 14
Overall Width	22.5 in / 57 cm	24.5 in / 62 cm
Folded without rear wheels	35x17x22.5 in 89x43x57 cm	35x17x24.5 in 89x43x62 cm
Folded with wheels	36x21x22.5 in 92x51x57 cm	36x21x24.5 in 92x51x62 cm

MEASUREMENT WORKSHEET

- A. Seat Depth _____
- B. Trunk Depth _____
- C. Trunk Width _____
- D. Hip Width _____
- E. Seat to Top of Shoulder _____
- F. Lower Leg Length _____
- G. Seat to Axilla _____



Adjustability

	TRAK™ 12	TRAK™ 14
Seat Depth	6.5-15 in / 17-38 cm	6.5-15 in / 17-38 cm
Seat Width	6-12 in / 15-30 cm	8-14 in / 20-36 cm
Lower Leg Length	4-17 in / 10-43 cm	4-17 in / 10-43 cm
Back Height	12-19 in / 30-48 cm	13-19 in / 33-48 cm
Hip Angles	90°, 105°, 120°	90°, 105°, 120°

TRAK™ Chair

TRAK SE™

<input type="checkbox"/>	3002-1004	E1234	TRAK SE™ 12	\$2,395.00
<input type="checkbox"/>	3002-1002	E1234	TRAK SE™ 14	\$2,395.00

Upholstery Color

Upholstery Color applies to RESPOND™ back and seat cushions, sunshade.

<input type="checkbox"/>	3002-0101	BIG SKY BLUE - Blue/Black	N/C
<input type="checkbox"/>	3002-0102	ROSSO RED - Red/Black	N/C
<input type="checkbox"/>	3002-0103	ICE SILVER - Silver/Black	N/C
<input type="checkbox"/>	3002-0104	SNOWBERRY PINK - Pink/Black	N/C
<input type="checkbox"/>	3002-0105	HUCKLEBERRY PURPLE - Purple/Black	N/C

Transit Option & Transit Harness Selection

Transit Option

Select belt or harness to right, 3003-4001 back rest and a head rest.

<input type="checkbox"/>	3002-7001	Transit Option (RESNA WC-19)	\$250.00
--------------------------	-----------	------------------------------	----------

This option is not retrofittable.

Transit Belt or Harness

<input type="checkbox"/>	3002-7002	K0108	Transit - 5-Point Harness	\$165.00
		<input type="checkbox"/>	Add pad set for 5-point harness (3002-7022)	\$90.00
<input type="checkbox"/>	3001-7008	K0108	Transit - Sure-Lok Lap Belt	\$90.00

Select 5-PT Harness for occupants less than 50 lbs./23kg

Select Sure-Lok Lap Belt for occupants over 50 lbs./23kg

The transit option includes 4 installed transportation anchors & has been tested in compliance with WC-19 standards. For more information on RESNA WC-19 please visit www.wc-transportation-safety.com.

Leggero 5-Point Positioning Harness

Not compatible with transit option.

<input type="checkbox"/>	3005-6015-12	<input type="checkbox"/>	3005-6015-14	5-Point Positioning Harness (-## indicates size)	\$165.00
		<input type="checkbox"/>		Add pad set for 5-point harness (3002-7022)	\$90.00

Front Rigging & Foot Plate

Front Rigging

<input type="checkbox"/>	3002-7015		90° Fixed Leg Rest	N/C
<input type="checkbox"/>	3002-7006	E0990	Angle Adjustable Leg Rest	\$225.00

Foot Plate

<input type="checkbox"/>	3002-5001		Standard Height Adjustable	N/C
<input type="checkbox"/>	3002-5002	K0040	Height & Angle Adjustable	\$165.00

RESPOND™ Seat Cushions

<input type="checkbox"/>	3002-3006	E2231	Solid Seat Pan	\$165.00	<input type="checkbox"/>	3002-7099		Omit Seat Cushion	N/C
--------------------------	-----------	-------	----------------	----------	--------------------------	-----------	--	-------------------	-----

Select from A or B. If B is chosen, continue selection to right.

A. RESPOND™ Planar Seat Cushion

<input type="checkbox"/>	3002-3004	E2292	TRAK 12 Cascading Planar Cushion	\$295.00
<input type="checkbox"/>	3002-7024	E0981	TRAK 12 Extra Planar Cushion Cover	\$70.00
<input type="checkbox"/>	3002-3001	E2292	TRAK 14 Cascading Planar Cushion	\$295.00
<input type="checkbox"/>	3002-7010	E0981	TRAK 14 Extra Planar Cushion Cover	\$70.00

B. RESPOND™ Contoured Seat Cushions

<input type="checkbox"/>	3002-3005	E2292	TRAK 12 Cascading Contoured Cushion	\$365.00
<input type="checkbox"/>	3002-7025	E0981	TRAK 12 Extra Cushion Cover	\$70.00
<input type="checkbox"/>	3002-3002	E2292	TRAK 14 Cascading Contoured Cushion	\$365.00
<input type="checkbox"/>	3002-7011	E0981	TRAK 14 Extra Cushion Cover	\$70.00
<input type="checkbox"/>	3001-7011	K0108	Lumbar EVA & 4 EVA Wedges	\$50.00

Ischial Insert Options (Only if B is selected)

<input type="checkbox"/>	3003-3010		TRAK 12 PU Foam Ischial Insert	N/C
<input type="checkbox"/>	3003-3006		TRAK 14 PU Foam Ischial Insert	N/C
<input type="checkbox"/>	3003-3007	K0108	TRAK 14 Soft Mem. Foam Ischial Insert	\$50.00
<input type="checkbox"/>	3003-3008	K0108	TRAK 14 Med. Mem. Foam Ischial Insert	\$50.00

Please continue on next page.

RESPOND™ Back

<input type="checkbox"/>	3003-4002	E2293	TRAK 12 Contoured Config. Back	\$525.00
<input type="checkbox"/>	3003-3013	E0982	TRAK 12 Extra Back Cushion Cover	\$70.00
<input type="checkbox"/>	3003-4001	E2293	TRAK 14 Contoured Config. Back	\$525.00
<input type="checkbox"/>	3003-3009	E0982	TRAK 14 Extra Back Cushion Cover	\$70.00
<input type="checkbox"/>	3001-7011	K0108	Lumbar EVA & 4 EVA Wedges	\$50.00

<input type="checkbox"/>	3003-4099		Omit Back	N/C
--------------------------	-----------	--	-----------	-----

When omitting back, crossbars between canes will remain in place.

Lateral/Pelvic Thigh Supports

<input type="checkbox"/>	3003-6006	E0956	LEFT, Width/Depth/Angle Adj.	\$120.00
<input type="checkbox"/>	3003-6007	E0956	RIGHT, Width/Depth/Angle Adj.	\$120.00

Armrests & Positioning Tray

3002-6017 Armrests & 3002-6016 Hardware are required when ordering Tray Platform

<input type="checkbox"/>	3002-6017	E0973	Height/Angle Adj. Armrests	\$195.00
--------------------------	-----------	-------	----------------------------	----------

<input type="checkbox"/>	3002-6021	E0950	TRAK 12 Clear Positioning Tray Platform	\$150.00
<input type="checkbox"/>	3002-6015	E0950	TRAK 14 Clear Positioning Tray Platform	\$150.00
<input type="checkbox"/>	3002-6016	E1028	Tray Hdw., Height & Angle Adjustable	\$209.00

Ankle Supports

<input type="checkbox"/>	3001-6007	E0951	LEFT, Extra Small (Pad 2x8")	\$72.00
<input type="checkbox"/>	3001-6008	E0951	RIGHT, Extra Small (Pad 2x8")	\$72.00
<input type="checkbox"/>	3001-6009	E0951	LEFT, Small (Pad 3x10")	\$72.00
<input type="checkbox"/>	3001-6010	E0951	RIGHT, Small (Pad 3x10")	\$72.00

Trunk & Pelvic Positioning

Stealth i-FIT Pelvic Positioning Belts

<input type="checkbox"/>	IFPB5141	E0978	Push Button 1", XS Pads	\$130.00
<input type="checkbox"/>	IFPB5242	E0978	Push Button 1.5", SM Pads	\$130.00
<input type="checkbox"/>	IFPB5111	E0978	2-Point, Push Button 1", 4.5" Pads	\$86.00
<input type="checkbox"/>	IFPB5112	E0978	2-Point, Push Button 1", 5.5" Pads	\$86.00

Buckle Covers (compatible with Stealth positioning belts)

<input type="checkbox"/>	IFMDBC10		Buckle Cover Kit 1" Buckle	\$25.00
<input type="checkbox"/>	IFMDBC15		Buckle Cover Kit 1.5" Buckle	\$25.00

Abductor/Medial Knee Block

Must select hardware & one size of abductor pad. Wedge sizes are width x depth x height.

<input type="checkbox"/>	3005-7001	E1028	Leggero Flip Down Medial Hdw.	\$240.00
<input type="checkbox"/>	MTMD-12	E0957	Small Wedge Pad (3x5x3")	\$150.00
<input type="checkbox"/>	MTMD-13	E0957	Medium Wedge Pad (5x5x3")	\$150.00

Anterior Posture Supports

<input type="checkbox"/>	IFPS1101-02	E0960	i-Fit Classic, Dynamic, XS	\$142.00
<input type="checkbox"/>	IFPS1102-02	E0960	i-Fit Classic, Dynamic, S	\$142.00
<input type="checkbox"/>	IFPS1103-02	E0960	i-Fit Classic Dynamic, M	\$142.00
<input type="checkbox"/>	IFPS1201-02	E0960	i-Fit Classic, Structured, XS	\$142.00
<input type="checkbox"/>	IFPS1202-02	E0960	i-Fit Classic, Structured, S	\$142.00
<input type="checkbox"/>	IFPS1203-02	E0960	i-Fit Classic, Structured, M	\$142.00
<input type="checkbox"/>	IFPS2101-02	E0960	i-Fit Contoured, Dynamic, XS	\$142.00
<input type="checkbox"/>	IFPS2102-02	E0960	i-Fit Contoured, Dynamic, S	\$142.00
<input type="checkbox"/>	IFPS2103-02	E0960	i-Fit Contoured, Dynamic, M	\$142.00
<input type="checkbox"/>	IFPS2201-02	E0960	i-Fit Contoured, Structured, XS	\$142.00
<input type="checkbox"/>	IFPS2202-02	E0960	i-Fit Contoured, Structured, S	\$142.00
<input type="checkbox"/>	IFPS2203-02	E0960	i-Fit Contoured, Structured, M	\$142.00
<input type="checkbox"/>	3005-6016	E0960	Leggero 'H' Harness, S	\$165.00
<input type="checkbox"/>	3005-6017	E0960	Leggero 'H' Harness, M	\$165.00
<input type="checkbox"/>	3005-6018	E0960	Leggero 'H' Harness, L	\$165.00

Swing Away Lateral Hardware & Pads

Requires RESPOND™ Back 3003-4001. Must select hardware & pad.

<input type="checkbox"/>	3002-6111	E1028	LEFT, Leggero Swing In/Out Lateral	\$225.00	
<input type="checkbox"/>	3002-6112	E1028	RIGHT, Leggero Swing In/Out Lateral	\$225.00	
<input type="checkbox"/>	TWBTL-L-01	E1028	LEFT, Stealth Thin Lateral	\$236.00	
<input type="checkbox"/>	TWBTL-R-01	E1028	RIGHT, Stealth Thin Lateral	\$236.00	
L	R	3004-6035	E0956	Lateral Pad 3x5", Each	\$75.00
L	R	3004-6045	E0956	Lateral Pad 4x5", Each	\$75.00
L	R	3004-6046	E0956	Lateral Pad 4x6", Each	\$75.00

Fixed Trunk Supports

Requires RESPOND™ Back 3003-4001

<input type="checkbox"/>	3003-6001	E0956	LEFT, 3x4" (7.5x10cm)	\$110.00
<input type="checkbox"/>	3003-6002	E0956	RIGHT, 3x4" (7.5x10cm)	\$110.00

Posture Support Strap Risers

<input type="checkbox"/>	3001-6011	K0108	Strap Risers, Pair	\$45.00
--------------------------	-----------	-------	--------------------	---------

Strap Risers are not compatible with LTV Bracket Accessories

Head Positioning

Select a system from options **A, B, or C**. Then continue your selections to the right to complete the package.

A. Leggero Head Positioning Options

<input type="checkbox"/>	3001-6005	E0955	8" Oval Headrest	\$240.00
<input type="checkbox"/>	3002-6003	E0955	RESPOND™ Cradle Headrest	\$240.00
<input type="checkbox"/>	3005-6099		Omit Head Positioning	N/C

Hardware is included with 8" Multi-Adjustable & RESPOND™ Headrests

B. Stealth Link Hardware and Comfort Plus Headrest

<input type="checkbox"/>	TWB480-1-02	E1028	Link Hdw. 1" Ball for CP Series	\$231.00
--------------------------	-------------	-------	---------------------------------	----------

Comfort Plus Headrest (Only if B is selected)

<input type="checkbox"/>	CP250-SR	E0955	Comfort Plus 10"	\$152.00
<input type="checkbox"/>	CP450-SR	E0955	Comfort Plus 8"	\$152.00

C. Stealth Link Hardware and i2i Head Positioning System

<input type="checkbox"/>	TWB480-17-02	E1028	Link Hdw. 17mm Ball/Flip Down	\$320.00
--------------------------	--------------	-------	-------------------------------	----------

Combo Series and i2i Head Positioning System (Only if C is selected)

<input type="checkbox"/>	i2i250	K0108	i2i Low Pro Nino XS, Pad Only	\$348.00
<input type="checkbox"/>	i2i350	K0108	i2i Low Pro Nino SM, Pad Only	\$348.00
<input type="checkbox"/>	i2i450	K0108	i2i Low Pro Adult SM, Pad Only	\$348.00
<input type="checkbox"/>	CB315	E0955	Combo Series, Small	\$295.00
<input type="checkbox"/>	CB305	E0955	Combo Series, Medium	\$295.00
<input type="checkbox"/>	CB325	E0955	Combo Series, Large	\$295.00

Chin Prompts (Only if i2i Head Positioning is selected)

<input type="checkbox"/>	i2i460	K0108	i2i SM Chin Prompt	\$80.00
<input type="checkbox"/>	i2i560	K0108	i2i MED Chin Prompt	\$80.00

Medical Mounting Accessories

Must be factory installed.

<input type="checkbox"/>	3002-7016	K0105	IV Pole	\$150.00
<input type="checkbox"/>	3002-7017	E2208	Oxygen Tank Holder (D size)	\$175.00
<input type="checkbox"/>	3002-7019	E1029	LTV Bracket for LTV950 ¹	\$250.00

<input type="checkbox"/>	3002-7027	E1029	TRAK 12 Equipment Tray	\$350.00
<input type="checkbox"/>	3002-7028	E1029	TRAK 12 LTV Bracket ¹ and Equip. Tray	\$495.00
<input type="checkbox"/>	3002-7018	E1029	TRAK 14 Equipment Tray	\$350.00
<input type="checkbox"/>	3002-7021	E1029	TRAK 14 LTV Bracket ¹ and Equip. Tray	\$495.00

Combined weight of accessories & occupant must not exceed chair capacity.

Ventilator, oxygen tank, and any other accessories must be removed and secured off chair during transit in a motor vehicle.

¹ LTV Bracket options are not compatible with 3001-6011 Strap Risers.

Anti-Tips and Accessories

<input checked="" type="checkbox"/>	3002-7007	E0971	Rear Anti-Tips	\$140.00
<input checked="" type="checkbox"/>	3002-7004		Under Seat Storage	N/C
<input checked="" type="checkbox"/>	3002-7012		Deluxe Sunshade	N/C
<input type="checkbox"/>	3002-6113		TRAK 12 Grab Bar	\$135.00
<input type="checkbox"/>	3002-7026		TRAK 12 Rain Cover	\$195.00
<input type="checkbox"/>	3002-6019		TRAK 14 Grab Bar	\$135.00
<input type="checkbox"/>	3002-7014		TRAK 14 Rain Cover	\$195.00

The **HCPCS CODES** provided are based on code verification through PDAC or our interpretation of the code definitions and are not to be deemed as legal advice. The use of **HCPCS CODES** does not ensure coverage or payment. Appropriate billing codes are the responsibility of DME providers. For coverage information, verify the appropriate payer's coverage policy.