



## Lopital *Tango XL* Shower Chair

More and more care facilities are turning to showering for resident's personal hygiene. The **Tango XL** is an affordable bariatric shower chair providing resident care in a responsible manner. The Tango XL provides optimum comfort for the caregiver as well as the resident and can be used as both a shower chair or wheeled over the toilet.

The Tango XL is a true "intermediate" shower chair, purpose designed to fill the gap between the Tango and the Tango XXL for care recipients too big for the Tango but do not require the size or capacity of the Tango XXL.



Hygienic split seat

PN: 5100 5750



Excellent over toilet clearance



### Features

- ◆ Wide wheelbase for stability.
- ◆ Large diameter swivel castors with brakes.
- ◆ Equipped with handles that improve maneuverability and drivability.
- ◆ Pivot away padded armrests.
- ◆ Rail for plastic commode pan.
- ◆ Swing away and removable footrests.
- ◆ Padded backrest and padded split seat seat.
- ◆ Suitable for use in the shower and can be moved into position over toilet.
- ◆ Corrosion resistant frame and polymer parts.
- ◆ Light and durable frame.

### Technical specifications

- ◆ Weight capacity: 550 lbs.
- ◆ Overall dimensions: 39" L x 26" W x 39. H
- ◆ Seat height: 21.6"
- ◆ Seat dimensions: 22" W / 19.6" Depth
- ◆ Castor diameter: 5"

