

SmartDrive

Featuring
SwitchControl

permobil





SwitchControl

SwitchControl gives users another programmable option for controlling their SmartDrive. With the push of a button SmartDrive can be engaged for a momentary burst of power or activated in latched mode for consistent power over

extended distances. This empowers a more diverse group of people to experience SmartDrive with increased confidence and when paired with the PushTracker E2 offers added flexibility, convenience and control options.

SmartDrive Specifications

Dimensions	Length: 15.3" (389 mm) Width: 5.7" (141 mm) Height: 9.5" (242 mm)
Product weight	12.5 lbs (5.7 kg) Drive Unit
User weight	30 to 331 lbs (14 to 150 kg)
Driving range*	Up to 12.3 mi (19.8 km)
Motor power	250 W Brushless DC
Maximum speed	5.5 mph (8.8 km/h)
Operating temp	-13°F to 122°F (-25°C to 50°C)
Chair type compatibility	Rigid, Folding, Tilt-In Space, One-Arm Drive and others
Drive wheel diameter	22", 24", 25", 26" (501, 540, 559, 590 mm)
Warranty	2-Year Limited Warranty

Bluetooth® LE	FCC ID QOQBT113 IC ID 5123A-BGTBLE113
HCPCS	E0986

**Performance measurements will vary based on user and chair characteristics, driving and battery conditions.*

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

NOTE: A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.

Please refer to your owner's manual for more detailed operation and care instructions.

Since U.S. Medicare coding and coverage criteria is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. ©2019 Permobil. Permobil has a policy of continual product improvement and reserves the right to amend this document.