

## Cast Padding

McKesson Cast Padding is designed to help protect skin under casting materials. Offered in four sizes, with minimal application time. Our superior quality cast padding provides excellent comfortability and cushioning. Not made with natural rubber latex.

### Features and Benefits

- Non-woven polyester padding provides protection and patient comfort while helping reduce water retention for a drier cast
- Strong and durable padding is easy to tear
- Conforms easily to body contours
- Porous to allow the skin to breathe
- Suitable with all casting materials, plaster and synthetic
- Color-coded packaging for easy size identification



### ORDERING INFORMATION

DESCRIPTION	MFR #	QTY
Cast Padding, 2", Synthetic, Non-Sterile	16-CP2	20ea/bg, 4bg/cs
Cast Padding, 3", Synthetic, Non-Sterile	16-CP3	20ea/bg, 4bg/cs
Cast Padding, 4", Synthetic, Non-Sterile	16-CP4	20ea/bg, 4bg/cs
Cast Padding, 6", Synthetic, Non-Sterile	16-CP6	20ea/bg, 4bg/cs

#### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

#### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.