

## McKesson casting tapes provide your patients the mobility they need.

Help your patients return to their normal daily activities with McKesson Medi-Pak™ Performance casting tape. Our casting tape combines the benefits of fiberglass with the handling ease of plaster. High and early strength allows for weight bearing in 20 minutes and provides a durable cast with less chance of breakdown and the need for costly repairs or replacements.

### Easy To Use

Our casting products have the same performance as national brands so you won't have to learn new application techniques. Water-activated resin lets you control working time according to application needs.

### All You Need

McKesson offers a variety of colors to meet patient preferences and a comprehensive line of accessories to support your casting needs.

**LATEX  
FREE!**

The product information contained in this flier, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers, not McKesson Medical-Surgical or its affiliates ("McKesson") and have not been independently verified by McKesson. McKesson is not responsible for errors or omissions in the product information. All trademarks and registered trademarks are the property of their respective owners.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the flier, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print fliers by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this flier is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

# Casting Products



Medi-Pak™ Performance  
Non-Sterile Tubular Stockinettes

- Knitted from high-grade yarn and inspected to guarantee consistent knit construction, cleanliness, and proper roll dimensions
- Used under plaster casts for skin protection
- Provides pressure to hold dressing in place
- Use with orthopedic appliances
- Packaged in an easy-to-use dispenser box on a continuous roll, allowing user to determine length



Synthetic

Synthetic stockinettes (polyester) designed for use with synthetic casting systems because of wicking action and fast-drying features.

Size	Color	Quantity	Catalog No.	Item No.
2" x 25 yds.	White	1 roll/bx, 12 bx/cs	16-4T-225	457643
3" x 25 yds.	White	1 roll/bx, 12 bx/cs	16-4T-325	457646
4" x 25 yds.	White	1 roll/bx, 10 bx/cs	16-4T-425	457644
6" x 25 yds.	White	1 roll/bx, 6 bx/cs	16-4T-825	457645

Cotton

Size	Color	Quantity	Catalog No.	Item No.
2" x 25 yds.	Beige	1 roll/bx, 12 bx/cs	16-4T-225	457629
3" x 25 yds.	Beige	1 roll/bx, 12 bx/cs	16-4T-325	457630
4" x 25 yds.	Beige	1 roll/bx, 10 bx/cs	16-4T-425	457631
6" x 25 yds.	Beige	1 roll/bx, 6 bx/cs	16-4T-625	457632
8" x 25 yds.	Beige	1 roll/bx, 4 bx/cs	16-4T-825	457633

Medi-Pak™ Performance Cast Tape

- Excellent lamination
- Lightweight for easy mobility
- Stretches to conform well to contours for a more comfortable, better-fitting cast
- Porous weave allows the skin to breathe
- Low tack resin—easy to unwind during application and molding using standard examination gloves

Color	Size	Quantity	Catalog No.	Item No.
White	2" x 4 yards	10 rolls/box	115-2	733683
	3" x 4 yards	10 rolls/box	115-3	733690
	4" x 4 yards	10 rolls/box	115-4	733697
	5" x 4 yards	10 rolls/box	115-5	733700
Black	2" x 4 yards	10 rolls/box	115-2A	733684
	3" x 4 yards	10 rolls/box	115-3A	733691
	4" x 4 yards	10 rolls/box	115-4A	733698
Blue	2" x 4 yards	10 rolls/box	115-2B	733685
	3" x 4 yards	10 rolls/box	115-3B	733692
	4" x 4 yards	10 rolls/box	115-4B	733699
Green	2" x 4 yards	10 rolls/box	115-2G	733686
	3" x 4 yards	10 rolls/box	115-3G	733693
Red	2" x 4 yards	10 rolls/box	115-2R	733687
	3" x 4 yards	10 rolls/box	115-3R	733694
Purple	2" x 4 yards	10 rolls/box	115-2U	733688
	3" x 4 yards	10 rolls/box	115-3U	733695
Bright Pink	2" x 4 yards	10 rolls/box	115-2X	733689
	3" x 4 yards	10 rolls/box	115-3X	733696

Medi-Pak™ Performance Synthetic Cast Padding

- Non-woven polyester padding provides protection and patient comfort while helping reduce water retention for a drier cast
- Strong and durable padding is easy to tear
- Conforms easily to body contours
- Porous to allow the skin to breathe
- Suitable with all casting materials, plaster and synthetic



Size	Color	Quantity	Catalog No.	Item No.
2"	White	20 rolls/bg, 4 bg/cs	115-CP2	734224
3"	White	20 rolls/bg, 4 bg/cs	115-CP3	734225
4"	White	20 rolls/bg, 4 bg/cs	115-CP4	734226
6"	White	20 rolls/bg, 4 bg/cs	115-CP6	734227

Casting tape available in a variety of colors:\*



\* Colors may vary slightly from those pictured.