

## Bath Transfer Systems Measurement Guide

Bath/ commode chairs and slider systems are known for their ability to fit in cramped spaces others simply can't. While we have yet to find a bathroom we can't accommodate, a few simple measurements will insure the best fit of your chair. Follow the descriptions below, and enter the measurements. Round up or down to 1/2-inch increments, and use the decimal system (for example, 10.5 for 10 and-half inches).

**A** Overall length of bathtub:

\_\_\_\_\_

**B** Length of bathtub floor:

\_\_\_\_\_

**C** Width of bathtub floor (flat portion only - do not include any rounded area)

\_\_\_\_\_

**D** Overall width of bathtub:

\_\_\_\_\_

**E** Distance between edge of toilet bowl and bathtub:

\_\_\_\_\_

**F** Widest portion of toilet bowl (including "wedge" shape up to face tank)

\_\_\_\_\_

**G** Distance between edge of toilet bowl to closest objects/wall:

\_\_\_\_\_

**H** Distance from front edge of toilet bowl to closest objects/wall:

\_\_\_\_\_

**K** Distance from center of toilet bowl to center of bathtub drain/spout:

\_\_\_\_\_

**L** Distance from top edge of bathtub to floor of bathtub:


\_\_\_\_\_

**M** Distance from top edge of bathtub to bathroom floor:

\_\_\_\_\_

**N** Distance from top of toilet bowl to bathroom floor (with seat and lid in raised position)

\_\_\_\_\_



Note: Please verify the measurements of the product, to ensure proper fitting for the client and adequate space availability in the user environment. medicaleshop will not be responsible for incorrect measurements/fitting and inadequate space availability in the user environment.

Need help with measurements, Please call us...

**Medicaleshop Inc.**

87 Danbury Road, Unit 1, New Milford, CT 06776, Phone : (866) 563-6812 | Fax : (860) 838-4671

Email : [csteam@medicaleshop.com](mailto:csteam@medicaleshop.com)