

## Guardian recliners

- Infinite position reclining up to 180°
- Durable hammertone frame

## Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Elevating	KR168N21E	Black	300 lbs. (136.1 kg)	16" x 18"	K0006+E0990+E1226	N	1 ea
		KR188N21E	Black	300 lbs. (136.1 kg)	18" x 18"	K0006+E0990+E1226	N	1 ea
		KR208N21E	Black	300 lbs. (136.1 kg)	20" x 18"	K0006+E0990+E1226	N	1 ea
Full-length removable	Elevating	KR188N11E	Black	300 lbs. (136.1 kg)	18" x 18"	K0006+E0990+E1226	N	1 ea

## Reclining wheelchairs

- Infinite position reclining from 90°–140°
- Durable carbon steel frame
- Hemi-height adjustable

## Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Elevating	MDS808350	Chrome	300 lbs. (136.1 kg)	16" x 17"	K0003+E0966+E1225+K0195+E0971	N	1 ea
		MDS808450	Chrome	300 lbs. (136.1 kg)	18" x 17"	K0003+E0966+E1225+K0195+E0972	N	1 ea
		MDS808550	Chrome	300 lbs. (136.1 kg)	20" x 17"	K0003+E0966+E1225+K0195+E0973+E2201	N	1 ea
		MDS808650	Chrome	300 lbs. (136.1 kg)	22" x 17"	K0007 + E1226 + E2201 + E0971 + E0966 + K0195 or E0990	N	1 ea

## Hybrid wheelchairs

- Combination wheelchair and transport chair
- Quick release large wheel

## Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions	HCPCS	Retail Pkg.	Pkg.
Desk-length, removable	Swing-away	MDS806250NH2	Blue	300 lbs. (136.1 kg)	16" x 16"	K0004 or E1038	N	1 ea
		MDS806250H2	Blue	300 lbs. (136.1 kg)	18" x 16"	K0004 or E1038	N	1 ea
	Elevating	MDS806300H2	Blue	300 lbs. (136.1 kg)	18" x 16"	K0004 or E1038+K0195 or E0990	N	1 ea

KR168N21E



MDS808550



MDS806250H2

